



FDA Home<sup>3</sup> Medical Devices<sup>4</sup> Databases<sup>5</sup>

### Class 2 Device Recall OERPro Endoscope Reprocessor



510(k)<sup>7</sup>|DeNovo<sup>8</sup>|Registration & Listing<sup>6</sup>|Adverse Events<sup>10</sup>|Recalls<sup>11</sup>|PMA<sup>12</sup>|HDE<sup>13</sup>|Classification<sup>14</sup>|Standards<sup>15</sup>  
CFR Title 21<sup>16</sup>|Radiation-Emitting Products<sup>17</sup>|X-Ray Assembler<sup>18</sup>|Medsun Reports<sup>19</sup>|CLIA<sup>20</sup>|TPLC<sup>21</sup>

New Search

#### Class 2 Device Recall OERPro Endoscope Reprocessor

<b>Date Initiated by Firm</b>	June 23, 2016
<b>Create Date</b>	October 19, 2016
<b>Recall Status<sup>1</sup></b>	Open <sup>3</sup> , Classified
<b>Recall Number</b>	Z-0138-2017
<b>Recall Event ID</b>	74542 <sup>23</sup>
<b>510(K)Number</b>	K103264 <sup>24</sup>
<b>Product Classification</b>	Accessories, cleaning, for endoscope <sup>25</sup> - Product Code FEB <sup>26</sup>
<b>Product</b>	OER-Pro Endoscope Reprocessor, Endoscope washer/disinfectant
<b>Code Information</b>	Model: OER-Pro, All serial numbers affected
<b>Recalling Firm/ Manufacturer</b>	Olympus Corporation of the Americas 3500 Corporate Pkwy PO Box 610 Center Valley PA 18034-0610
<b>For Additional Information Contact</b>	Laura Storms 484-896-5688
<b>Manufacturer Reason for Recall</b>	New warning requiring Users to conduct duodenoscope precleaning and manual cleaning even when using an Automated Endoscope Reprocessor (AEI) may indicate a user could forego certain steps in precleaning and manual cleaning of the endoscopes.
<b>FDA Determined Cause<sup>2</sup></b>	Under Investigation by firm
<b>Action</b>	Olympus America Inc. (OAI) mailed a letter to customers informing them that they are issuing an updated Operation Manual for the OER-Pro Automated following actions were asked to be taken: -Olympus has discontinued previously distributed copies of the OER-Pro Operation Manual. - Inspect your inventory and discard any existing inventory of OER-Pro Operation Manuals. - Implement use of the enclosed OER-Pro Operation Manual, which contains the new 80. - Ensure all reprocessing personnel are completely knowledgeable and thoroughly trained on the requirement to perform precleaning and manual clean duodenoscopes prior to placing the duodenoscopes in the OER-Pro for high level disinfection. - Additional copies of the new OER-Pro Operation Manual contacting our Technical Assistance Center at 1-800-848-9024, option 1, or by indicating on the enclosed questionnaire. Additional Operation Manuals will Please indicate on the enclosed questionnaire that you have received this notification. Fax the completed form to (484) 896-7128.
<b>Quantity in Commerce</b>	2, 686 units
<b>Distribution</b>	Distributed to: AK,AL,AR,AZ,CA,CO,CT,DC,DE,FL,GA,HI,IA,ID,IL,IN,KS,KY,LA,MA,MD,ME,MI,MN,MO,MS,MT,NC,ND,NE,NH,NJ,NM,NV,NY,OH,OK,OR,PA,RI,SC,TX
<b>Total Product Life Cycle</b>	<u>TPLC Device Report<sup>27</sup></u>

<sup>1</sup> A record in this database is created when a firm initiates a correction or removal action. The record is updated if the FDA identifies a violation and classifies the action as a recall, and when the recall is terminated. Learn more about [medical device recalls<sup>28</sup>](#).

<sup>2</sup> Per FDA policy, recall cause determinations are subject to modification up to the point of termination of the recall.

<sup>3</sup> The manufacturer has initiated the recall and not all products have been corrected or removed. This record will be updated as the status changes.

**510(K) Database** 510(K)s with Product Code = FEB and Original Applicant = OLYMPUS MEDICAL SYSTEMS CORPORATION<sup>29</sup>

#### Links on this page:

1. <http://www.addthis.com/bookmark.php?u508=true&v=152&username=fdamain>
2. <http://www.addthis.com/bookmark.php>
3. <http://www.fda.gov/default.htm>
4. <http://www.fda.gov/MedicalDevices/default.htm>
5. <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Databases/default.htm>
6. </scripts/cdrh/devicesatfda/index.cfm>
7. </scripts/cdrh/cfdocs/cfPMN/pmn.cfm>
8. </scripts/cdrh/cfdocs/cfpmn/denovo.cfm>
9. </scripts/cdrh/cfdocs/cfRL/rl.cfm>
10. </scripts/cdrh/cfdocs/cfMAUDE/TextSearch.cfm>
11. </scripts/cdrh/cfdocs/cfRES/res.cfm>
12. </scripts/cdrh/cfdocs/cfPMA/pma.cfm>
13. </scripts/cdrh/cfdocs/cfHDE/hde.cfm>
14. </scripts/cdrh/cfdocs/cfPCD/classification.cfm>
15. </scripts/cdrh/cfdocs/cfStandards/search.cfm>
16. </scripts/cdrh/cfdocs/cfCFR/CFRSearch.cfm>