

# URGENT MEDICAL DEVICE PRODUCT FIELD ACTION NOTIFICATION LFIT<sup>TM</sup> Anatomic CoCr V40<sup>TM</sup> Femoral Heads

August 29, 2016

Product Field Action Number: RA2016-028

Description: LFIT™ Anatomic CoCr V40™ Femoral Heads

Catalog Number(s): 6260-9-236, 6260-9-240, 6260-9-244, 6260-9-340, 6260-9-344, 6260-9-440, 6260-

9-444

Lot Code(s): See attached

Dear XXX,

Stryker has initiated a voluntary medical device product field action for the following Femoral Heads.

The intent of this letter is to describe all potential hazards associated with the below noted issue, and any risk mitigation factors associated with the use of the product.

Our records indicate that you have received the above referenced product. It is Stryker's responsibility as the manufacturer to ensure that customers who may have received these affected products also receive this important communication.

### Reason for the Voluntary Product Field action:

Stryker has received higher than expected complaints of taper lock failure for specific lots of the following certain sizes of LFIT™ Anatomic CoCr V40TM Femoral Heads manufactured prior to 2011.

Catalog Number	<b>Head Diameter</b>	Offset
6260-9-236	36mm	+5
6260-9-240	40mm	+4
6260-9-244	44mm	+4
6260-9-340	40mm	+8
6260-9-440	40mm	+12
6260-9-344	44mm	+8
6260-9-444	44mm	+12

### Potential Hazards may include:

- Disassociation of femoral head from hip stem
- Fractured hip stem trunnion
- Excessive metallic debris
- Insufficient ROM
- Insufficient soft tissue tension
- Noise
- Loss of implant: bone fixation strength
- Excessive wear debris (polymeric)
- Implant construct with a shortened neck length

The aforementioned potential hazards may result in one or more of the following potential patient harms:

- User annoyance
- Loss of mobility
- Pain requiring revision
- Inflammatory response
- Adverse local tissue reaction
- Dislocation
- Joint instability
- Revision to alleviate hazardous situation
- Pain associated with implant loosening
- Periprosthetic fracture
- Leg length discrepancy

#### Follow up:

Implanted patients with LFIT<sup>™</sup> Anatomic CoCr V40<sup>™</sup> Femoral Heads as described above should continue to be followed per the normal protocol established by his/her surgeon.

#### Required actions:

Our records indicate that you have received at least one of the subject devices and you are therefore affected by this action.

We request that you read this notice carefully and complete the following actions:

- 1. Please be aware that all affected products are either expired or already implanted. Check your internal inventory and in case you still have any product, quarantine all subject devices pending return to Stryker.
- 2. Circulate this Field Safety Notice internally to all interested/affected parties.
- 3. Maintain awareness of this notice internally until all required actions have been completed within your facility.
- 4. Inform Stryker if any of the subject devices have been distributed to other organisations.
  - a) Please provide contact details so that Stryker can inform the recipients appropriately.
  - b) If you are a Distributor, note that you are responsible for notifying your affected customers.
- 5. Please inform Stryker of any adverse events concerning the use of the subject devices.
  - a) Please comply with any local laws or regulations concerning the notification of adverse events to your National Competent Authority.
- 6. Complete the attached customer response form. Please complete even if you no longer have any of the subject devices in your physical inventory.
- 7. Return the completed form to your nominated Stryker Representative (indicated below) for this PFA

Your designated contact pe matter please do not hesitat	_	below.	Should you have any queries concerning this		
Name:	Position:	email:			
We request that you respond to this notice within XXX calendar days from the date of receipt. The target date for completion of this action is XXX and your timely response will enable us to ensure that we meet this target.					
- Local PM email and d		e find bel	low the contact for your country:		
	Please be advised that the ser		Il center in the US that works for calls from ars will be according to the US time zone and		
	_		ance document Ref 2.12-1, we can confirm that ent Authority for your country.		
date and regret any inconve	nience that may be caused. \	We woul	oport in completing this action within the target ld like to reassure you that Stryker is committed al quality standards, remain on the market.		
Yours					
Sincerely,					

# URGENT MEDICAL DEVICE RECALL NOTIFICATION LFIT™ Anatomic CoCr V40™ Femoral Heads

August 29, 2016

Product Field Action Number: Description: Catalog Number(s):	LFIT™ Anatomic CoCr V40™	<sup>n</sup> Femoral Heads 260-9-244, 6260-9-340, 6260-9-344, 6260-9-440, 6260
Lot Code(s):	See attached	
I have received the product re a voluntary, lot-specific produ		August 29, 2016 stating that the company has initiated ced product.
Stryker Branch / Agent / Hosp (Signature)	oital Representative	Date
Stryker Branch / Agent / Hosp (Print)	 pital Representative	Stryker Branch / Agency/Hospital Name

## PLEASE COMPLETE THIS FORM WITHIN 5 BUSINESS DAYS AND RETURN IT BY USING THE EMAIL OR FAX LISTED BELOW:

email: <u>strykerortho8402@stericycle.com</u>

fax: **888-912-8457**