SHANGHAI CHANNELMED IMPORT & EXPORT CO.,LTD

Tel:+8621-58354522

Fax:+8621-58354018 E-mail: info@channelmed.com

Http://www.channelmed.com

PAGE 1 OF EF015 FIELD SAFETY NOTICE European Form

URGENT FIELD SAFETY NOTICE

STERILE GAUZE SWAB WITH X-RAY 4"X4"

Type of Action: Removal

Date: Jan.24,2017 Ref: 201701-1

1. Manufacturer: Shanghai Channelmed Import and Export Co., Ltd.

2. Device Type: Single-use Devices

3. Description: Sterile Gauze Swabs with X-Ray 4" x 4" 16 Ply, 10pcs/pkt + Photo

4. Reason of Field Safety Correction Action :

The Saudi Food And Drugs Authority reflect the problem that the Sterile Gauze Swab with X-Ray products found those problem as below:

- 1. Impurities, dust particles, plankton, insect inside
- 2. Quantity inside each packet not fixed some times less or more than 10 pieces
- 3. Some of the piece of X-Ray gauze did not have X-Ray line detector inside.
- 5. Affected Code Information: Lot No.: 0020114-CJ14128 & 0020114-CJ14148 14558 & 14364 & 3760811&15216

6. Potential Hazard:

1. If forget the blue line detector not exist, the operating surgeon may use the Gauze swab for the operations and he can notidiscover it by X-Ray.

2. If the quantity is more offices, the users can not prepare the correct quantities of Gauze swab which them have se.

 If plankton and dusta how he countries the impurity substance may attached the wound.

SHANGHAI CHANNELMED IMPORT & EXPORT CO.,LTD

Tel:+8621-58354522

Fax:+8621-58354018 E-mail:info@channelmed.com

Http://www.channelmed.com

EF015 European	FIELD SAFETY NOTICE	PAGE 2	
Form		OF 2	

7. The Current Action:

· Please discard the affected Sterile Gauze Swab with X-Ray products and give compensation to customers free.

Please ensure that the information in this notice is made available to all affected technical aid centers, service companies and dealers within your organization, and that the advice stated in theses Manuals is adhered to.

8. Future Action:

- 1. We are trying to use new looms, and use the new type of gauze folding machine and it will avoid the debris falling.
- 2. We will add automatic weighing procedure by measuring machine to confirm the quantities of pieces after artificial calculation.
- 3. We will check it carefully during packaging process, and discard those pieces which do not have X-Ray line detector.

9. Confirmation of National Competent Authorities:

Please kindly note that this recall is made with the cooperation of Saudi Food and Drugs Authority because there are some remarks on the product,.

