

-		
	Catalog Number	Lot Number
	4400322X	82144115
	4400602S	82144141
	4400515X	82144947
	4400308S	82148810
	4400508S	82148811
	4400604S*	82144604*
	4400604S*	82144617*
	4400804S*	82144499*

<u>Update</u>: Urgent Field Safety Notice (Removal) Cordis[®] POWERFLEX[®] PRO PTA Dilatation Catheter

September 06, 2018, Updated August 23, 2019

Dear Valued Customer,

Previously, Cordis had notified you of a Field Safety Notice (Removal) in September 2018 regarding five (5) lots of Cordis[®] POWERFLEX[®] PRO Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter product (reference number Cordis20180906-OUS). The removal of these five lots from the impacted markets has been completed.

Cordis was recently informed that three (3) additional lots (see lot numbers marked with an asterisk "*" in tables above and below) are also affected by this action.

	related to shaft burst or leakage.
	Cordis has not received any complaints related to POWERFLEX® PRO that are
	There is no safety concern for patients that are treated successfully using product from these lots.
	A shaft leakage/burst that occurs during balloon inflation would likely create an inability to inflate or maintain pressure of the balloon component. The user may experience inflation difficulty/deflation difficulty of the balloon component. The most likely occurrence would result in a procedural delay, however damage to healthy intima, vessel spasm, ischemia, or additional intervention could result, but only occasionally and/or under unusual circumstances.
Recall Overview:	Cordis has determined that eight total lots of POWERFLEX® PRO PTA Dilatation Catheters have not met an internal manufacturing specification for shaft burst strength, though it meets the label claim (18 ATM).

Details on	Product involved					
Affected Devices,	Eight (8) lots are affected:					
to assist in	Catalog Number	Lot Number	Balloon Diameter	Balloon Length		
identification of	4400322X	82144115	3mm	22cm		
the product	4400602S	82144141	6mm	2cm		
involved:	4400515X	82144947	5mm	15cm		
	4400308S	82148810	3mm	8cm		
	4400508S	82148811	5mm	8cm		
	4400604S*	82144604*	6mm	4cm		

	4400604S*	82144617*	6mm	4cm
	4400804S*	82144499*	8mm	4cm
	femoral, ilio-femoral treatment of obstruc fistulae. The device self-expanding stent <u>Identification</u>	PRO PTA Catheter , popliteal, infra poplit tive lesions of native is also indicated for p is in the peripheral va g below is provided t	eal and renal arterio or synthetic arterio ostdilation of balloo sculature.	es and for the venous dialysis on-expandable and
Details on	Identification (Con			
Affected Devices, to assist in	Example lab	peling:		
identification of the product involved (Continued):	Ø8 mm X 4 cm 80 c 4 cm Ø8 mm Ø8 mm Ø	Control With Control of the Control	Angelopatie (PTA) / Conterne distance per PTA Properties (Contentionalises of the PTA / Benefacteer for PTA / Distance angelopaties (PTA / No (No 600 Contentionalises (PTA / Distance angelopaties (PTA / Distance ange	Af Andrew Tayle Control to PFA / International Control to PFA / International Control to PFA / International Control to PFA / A T SK 59 29 / / / / / / / / / / / / / / / / / /
	Cordis Corpsretion St201 N.W 60th Ave North Lakes, FL 30014, US4	REP Cerdis Casher Cehir Road, Cashel Ce Tippenry, Ireland 0086 1252 10.4e	Res. A_0001 of 0272	Assembled in Costs Rica

being contacted:	purchased the POWERFLEX [®] PRO lot numbers indicated in this letter.		
Actions requested	1) Read this Field Safety Notice (Removal) letter.		
on your part:			
	2) Immediately check your inventory to confirm whether you have any units from the affected lots in your possession. Identify and set aside any units from the affected lots in a manner that ensures the affected product will not be used. Check all storage and usage locations.		
	3) Review, complete, sign and return the enclosed Acknowledgement Form in accordance with the directions on the form.		
	 Return all affected product to the Cardinal Health distribution center. Please contact your local sales representative to facilitate return of the affected product, if necessary. Your sales representative will inform you of the product replacement or credit options. 		

5)	Share this letter with others in your facility who need to be made aware of this recall.
6)	Please contact any other facility who may have received the affected units of POWERFLEX [®] PRO product from your facility. If any units of the affected lots are found to be at the other facility, please arrange the return of the units.
7)	Maintain awareness of this notice until all affected product has been returned to Cordis.
8)	Keep a copy of this notice with the affected product.

Description of the problem:	What is the issue? Cordis became aware from the manufacturer that three additional lots of product may not meet the shaft subassembly burst strength specification. The lots were previously not identified by the manufacturer and were discovered by follow-up investigation.
	Why are we recalling this product? A shaft leakage/burst that occurs during balloon inflation would likely create an inability to inflate or maintain pressure of the balloon component. Additionally, the user may experience inflation difficulty/deflation difficulty of the balloon component. The most likely occurrence would result in a procedural delay, however damage to healthy intima, vessel spasm, ischemia, or additional intervention could result, but only occasionally and/or under unusual circumstances.
	There is no safety concern for patients that are treated successfully using product from these lots.
	What other actions is Cordis taking? Cordis has performed a root cause investigation and taken immediate corrective action. Cordis has not identified any other lots that may be affected. In keeping with our commitment to provide customers with quality products, Cordis has voluntarily decided to recall these eight (8) lots.
Available	If you have any questions regarding this recall please contact your local sales

Available Assistance:	If you have any questions regarding this recall, please contact your local sales representative or local sales office.

Additional Information:	Regulatory Notification
	The applicable regulatory agencies and notified body are being notified that Cordis is voluntarily taking this action.

We apologize for any inconvenience this communication may cause. We know that you place high value in our products, and we appreciate your cooperation in this matter. Cordis is committed to maintaining your confidence in the safety and quality of the products that Cordis supplies.

Respectfully yours,

Cordis Corporation



Urgent Field Safety Notice (Removal) Cordis[®] POWERFLEX[®] PRO PTA Dilatation Catheter

Catalog Number	Lot Number
4400322X	82144115
4400602S	82144141
4400515X	82144947
4400308S	82148810
4400508S	82148811

September 06, 2018

Dear Valued Customer,

The purpose of this communication is to inform you Cordis is recalling (removing) five (5) lots of Cordis[®] POWERFLEX[®] PRO Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter product.

Recall Overview:	Cordis has determined that five lots of POWERFLEX® PRO PTA Dilatation Catheters have not met an internal manufacturing specification for shaft burst strength, though it meets the label claim (18 ATM).
	A shaft leakage/burst that occurs during balloon inflation would likely create an inability to inflate or maintain pressure of the balloon component. The user may experience inflation difficulty/deflation difficulty of the balloon component. The most likely occurrence would result in a procedural delay, however damage to healthy intima, vessel spasm, ischemia, or additional intervention could result, but only occasionally and/or under unusual circumstances.
	There is no safety concern for patients that are treated successfully using product from these lots.
	Cordis has not received any complaints related to POWERFLEX® PRO that are related to shaft burst or leakage.

Details on	Product involved			
Affected Devices,	• Five (5) lots are	affected:		
to assist in	Catalog Number	Lot Number	Balloon Diameter	Balloon Length
identification of	4400322X	82144115	3mm	22cm
the product	4400602S	82144141	6mm	2cm
involved:	4400515X	82144947	5mm	15cm
	4400308S	82148810	3mm	8cm
	4400508S	82148811	5mm	8cm
	femoral, ilio-femoral, treatment of obstruct fistulae. The device i self-expanding stent	popliteal, infra poplit tive lesions of native is also indicated for p s in the peripheral va	or synthetic arteriove ostdilation of balloon	s and for the nous dialysis -expandable and

Details on Affected Devices, to assist in identification of the product involved (Continued):	Identification (Continued) Example labeling:
	Ø3 mm X 8 cm Ø3 mm X 8 cm Ø3 mm X 8 cm Ø3 mm X 8 cm
	80 cm 5F 15 mm 15 mm 03 mm 11 3.16 1 3 3.16 15 3.16 1 3 3.16 15 3.16 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 15 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26 1 3 3.16 16 3.26

Why you are	You are receiving this letter because our records indicate that you have
being contacted:	purchased the POWERFLEX [®] PRO lot numbers indicated in this letter.
Actions requested	1) Read this Field Safety Notice (Removal) letter.
on your part:	
	2) Immediately check your inventory to confirm whether you have any units from the affected lots in your possession. Identify and set aside any units from the affected lots in a manner that ensures the affected product will not be used. Check all storage and usage locations.
	3) Review, complete, sign and return the enclosed Acknowledgement Form in accordance with the directions on the form.
	4) Return all affected product to the Cardinal Health distribution center. Please contact your local sales representative to facilitate return of the affected product, if necessary. Your sales representative will inform you of the product replacement or credit options.
	5) Share this letter with others in your facility who need to be made aware of this recall.
	6) Please contact any other facility who may have received the affected units of POWERFLEX [®] PRO product from your facility. If any units of the affected lots are found to be at the other facility, please arrange the return of the units.
	7) Maintain awareness of this notice until all affected product has been returned to Cordis.
	8) Keep a copy of this notice with the affected product.


Description of the	What is the issue?
problem:	Cordis became aware that the product may not meet the shaft
	subassembly burst strength specification.
	subassembly buist strength specification.
	Why are we recalling this product?
	A shaft leakage/burst that occurs during balloon inflation would likely
	create an inability to inflate or maintain pressure of the balloon
	component. Additionally, the user may experience inflation
	difficulty/deflation difficulty of the balloon component. The most likely
	occurrence would result in a procedural delay, however damage to
	healthy intima, vessel spasm, ischemia, or additional intervention could
	result, but only occasionally and/or under unusual circumstances.
	······································
	There is no safety concern for patients that are treated successfully
	using product from these lots.
	What other actions is Cordis taking?
	Cordis has performed a root cause investigation and taken immediate
	corrective action. Cordis has not identified any other lots that may be
	affected. In keeping with our commitment to provide customers with
	quality products, Cordis has voluntarily decided to recall these five (5)
	lots.
	1013.
	If you have any guartiana garanding this goodly places contact your local cales

Available	If you have any questions regarding this recall, please contact your local sales
Assistance:	representative or local sales office.

Additional Information:	Regulatory Notification
	The applicable regulatory agencies and notified body are being notified that Cordis is voluntarily taking this action.

We apologize for any inconvenience this communication may cause. We know that you place high value in our products and we appreciate your cooperation in this matter. Cordis is committed to maintaining your confidence in the safety and quality of the products that Cordis supplies.

Respectfully yours,



Cordis Corporation