



Urgent Safety Communication

Intra-Aortic Balloon Pumps-Battery Power May Deplete Earlier than Expected

Device/ Product Name:	Intra-Aortic Balloon Pumps (IABPs): CARDIOSAVE(Hybrid, Rescue), CS300 and CS100/CS100i.	
Lot numbers/Serials:	Not specified	
Manufacturer:	Maquet Datascope Corp. / Getinge Group	
Problem:	 SFDA would like to bring your attention to the hazard associated with the battery power of the Intra-Aortic Balloon Pumps model: CARDIOSAVE Hybrid, CARDIOSAVE Rescue. The battery operating time may not last the expected time. Either 	
	 If the CARDIOSAVE runs out of battery power and line power or a replacement battery is not readily available, the patient may be seriously injured. 	
	This Safety Communication is being issued to ensure that all hospitals are aware of the issue and that adequate action is taken to mitigate potential risk to the patients.	
Recommendation/Actions:	 Health care providers shall follow each device's Operating Instructions Manual for recommendations on usage, charging, maintenance and storage of the system batteries, since battery run times and discharge cycles vary between IABP models. Ensure the IABP system is plugged into an AC outlet whenever possible during patient use to prevent the battery from depleting. Ensure the IABP system is plugged into an AC outlet when the system is not in use. The batteries should be kept at a full charge even when the IABP is not in use. When transporting patients within or between facilities, please refer to the system's Operating Instructions Manual for recommendations for portable/battery operation. For example: Prior to portable operation, the battery should be fully charged 	
	 Additional charged batteries should be on hand during transport, if applicable for the system 	





	 Ensure the batteries are properly seated in the battery compartment/charger 		
	 Periodically check battery run time and replace batteries as required, as recommended in each system's Operating Instructions Manual. A reduction in run time can occur over a battery's life for reasons such as age, storage temperature and discharge depth. Batteries should be replaced: After reaching the maximum number of charge-discharge cycles 		
	 When the battery provides less than the minimum expected run time 		
	- When the labeled lifetime of the battery is reached		
	- If the battery is broken, cracked, leaking or damaged		
	 Report to SFDA any adverse events suspected to be associated with affected Intra-Aortic Balloon Pumps model: CARDIOSAVE Hybrid, CARDIOSAVE Rescue through <u>National Center for Medical Devices</u> <u>Reporting (NCMDR)</u> Or <u>Saudi Vigilance</u> Prompt reporting of adverse events can help the SFDA identify and better understand the risks associated with medical devices. 		
Devices/Products photo:	<image/> <image/> <image/>		
Authorized Representative	Company name:	Al-Jeel Medical & Trading Co. LTD	
Details	Contact Person:	Meshal Aleshaiwy	
	Phone:	+966112168222 Ext:253, 234	
	Email:	aljeel.sfda@aljeel.com	





You **should** be aware of the mentioned risks in the notice and **contact** the Authorized Representative for corrective action.

Healthcare Professionals should **report** any adverse events suspected to be associated with affected devices above (or other Medical Devices) to:

National Center for Medical Devices Reporting. Medical Devices Sector Saudi Food and Drug Authority Postal Address: Saudi Arabia - Saudi Food and Drug Authority (3292) North Ring Road - Al Nafal Unit (1) Riyadh 13312 - 6288 Tel: +966 (11) 2038222 Ext: 2406, 2412 Fax: +966 (11) 2757245

SFDA committed to keep the identities of the reporter and the patient confidential.

For latest published Recalls/Alerts, please visit (NCMDR Website)

Sincerely, NCMDR Team

