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Therapy tissue ablation, cardiac ablation system.  
Thermocool SmartTouch catheter.

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- Therapy tissue ablation, cardiac ablation system. Thermocool SmartTouch catheter.
- Trade Mark: Biosense Webster Inc.
- Local Representative:

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Medicine and Health Care Products Regulatory Agency (UK) MHRA

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**URGENT FIELD SAFETY NOTICE  
MEDICAL DEVICE – VOLUNTARY FIELD REMOVAL**

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September 3, 2013

Dear Valued Customer,

The purpose of this communication is to inform you that Biosense Webster, Inc., a division of Johnson & Johnson Medical NV/SA (“Biosense Webster”) is initiating a voluntary field removal of the THERMOCOOL® SMARTTOUCH™ Catheter (All catalog numbers). This letter provides important information about the affected products and instructions on how you can return the product to Biosense Webster.

**Overview:**

After conducting a root cause analysis resulting from increased reports of irrigation fluid flow interruption, Biosense Webster has identified an issue in the production process of the THERMOCOOL® SMARTTOUCH™ Catheter, which can lead to occlusion of the irrigation fluid lumen. As a result, Biosense Webster is initiating a voluntary field removal of all affected products. To date, there have been no patient injuries or adverse events reported as a result of this defect. There is no concern for patients who have already been successfully treated with the device.

**Details on Affected Devices:**

**Indications for Use:**

The THERMOCOOL® SMARTTOUCH™ Catheter and related accessories are indicated for catheter-based cardiac electrophysiological mapping (stimulating and recording) and, when used in conjunction with a radiofrequency generator, for cardiac ablation.

**Actions Requested on Your Part:**

- Read the “Description of the Problem” section below carefully.
- Immediately identify and set aside all affected products in a manner that ensures the product will not be used. Maintain a copy of this letter with the affected THERMOCOOL® SMARTTOUCH™ Catheter until all units are returned to Biosense Webster.
- Sign and return the attached Voluntary Field Removal Certification Form in accordance with the instructions listed on the form.
- Arrange for return of all units of the THERMOCOOL® SMARTTOUCH™ Catheter that you may have in your inventory per the instructions on the Voluntary Field Removal Certification Form.
- Pass on this notice to anyone in your facility that needs to be informed.
- Maintain awareness of this notice until all affected products have been returned to Biosense Webster.
- If any of the affected THERMOCOOL® SMARTTOUCH™ Catheters have been forwarded to another facility, contact that facility and arrange for the return.

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Biosense Webster, Inc. THERMOCOOL® SMARTTOUCH™ Catheter