


14 October 2013

## Urgent Field Safety Notice – RECALL

### Incorrect Fluid Filter in Sight Chamber of Lifeshield, Latex-Free, Non-DEHP, Primary Plumset, 15 Micron Filter In Sight Chamber, Prepierced Y-Site, 272 CM



**Lifeshield, Latex-Free, Non-DEHP, Primary Plumset, 15 Micron Filter In Sight Chamber, Prepierced Y-Site, 272 CM**

**Lot Numbers: 14000-28 & 14001-38**

**Lot Numbers: 27112-5H & 30079-5H (14000-92-28),  
28187-5H (14001-92-38)**


**Q.FA.EMEA.2013.028**

**14<sup>th</sup> October 2013**

**Dear Healthcare Professional and Hospira Customer,**

Hospira Inc. (Hospira) has become aware that some 15 Micron Filter PlumSet IV Administration Sets have been manufactured with an incorrect fluid filter.

**Issue:** Hospira is issuing this urgent product recall for the PlumSets listed above as three (3) lots of the device identified above have been manufactured with the incorrect fluid filter. The correct fluid filter is identified by a white clip ring that holds the filter in place. The incorrect filter has a transparent clip ring. To date, there have been no reported complaints related to this issue.



**Risk to Health:**

Use of PlumSets manufactured with the incorrect fluid filter may generate allergic reactions, particularly when administering parenteral lipid emulsions. Severe allergic reactions may result in serious injury or death.

