



XX July 2021

**URGENT: FIELD SAFETY NOTICE – MDS-21-4144**

**BD SmartSite™ Needle Free Valve**

**REF & Lot Numbers:** See Appendix 1

**Type of Action:** Advisory

**Attention: Clinical Personnel, Risk Managers, Biomedical Personnel**

This letter contains important information which requires your **immediate** attention.

Dear Customer,

BD is issuing an advisory Field Safety Notice for specific lots of BD SmartSite™ Needle Free Valve products, as listed in Appendix 1, and according to our distribution records your organisation may have received the impacted product, which was shipped from 23<sup>rd</sup> January 2020.

**Description of the Problem**

Based on customer feedback, BD has identified potential issues with the valve of the Needle-free Connector for the products listed in Appendix 1 resulting in difficulty to flush, flow issues, and partial or total occlusions on the Extension Sets with SmartSite™ and Standalone SmartSite™ connectors. Only the lot numbers listed in Appendix 1 are affected by this issue.

BD has an online tool to support the identification of impacted lot numbers <https://www.bd.com/MDS-21-4144>.

**Clinical Impact**

The health consequences associated with the use of an occluded connector is a delay in therapy that can potentially lead to an injury that requires medical intervention.

BD has not identified any reports of serious adverse events to date that could be associated with this Field Safety Corrective Action. No specific patient follow-up activities are required if the product has already been used or already safely connected to the female luer of connecting device.

**Advice for Clinical Users**

As a result of this feedback and to mitigate the occurrence of the connector occluding, BD is recommending the following steps be followed prior to clinical use of the needle-free connector

- Depress the syringe plunger whilst the syringe is still attached to the needle-free connector.
- If unable to depress the plunger, pull back on the syringe plunger and depress the syringe plunger again. Repeating this step may assist in opening the needle-free connector valve.



- If the valve still is not open; without disconnecting the syringe male luer tip from the BD SmartSite™ valve surface, carefully unthread and rethread the syringe to partially disengage and reengage the syringe.
- Repeating this step up to three times may assist in opening the needle-free valve.  
*Note: To maintain sterility, do not allow the syringe and needle-free valve to completely disconnect*
- If previous steps do not open the valve, then replace with another BD SmartSite needle-free connector.
- After the steps above, proceed with the existing instructions for use for product.

**Actions by BD:**

1. BD has identified the root cause and initiated corrective actions to prevent recurrence of this issue.

**Actions for Customers to take:**

1. Circulate this Field Safety Notice to all those within your organisation that may use the BD SmartSite™ Needle Free Valve devices listed in Appendix 1.
2. Where a potential delay in therapy cannot be clinically tolerated (e.g., in an emergency situation) consideration should be given to administering the critical drug directly into the catheter hub or through a different route.
3. If you have further distributed the product, please identify those users and notify them at once of this Field Safety Notice.
4. If you continue to experience issues with the BD SmartSite™ Needle Free Valve report as a complaint to BD per your normal process.
5. Complete the customer response form on page 3 and return it to <<insert contact details here>> as soon as possible or no later than <<xx August 2021>>. If you no longer use the product, it is still important that you return the Customer Response Form for our reconciliation purposes.

**There is no requirement for customers to return any BD SmartSite™ Needle Free Valve devices to BD. These products can continue to be used in accordance with the guidance in this safety notice.**

**Contact Reference Person**

If you have any questions about this, please contact your local BD representative or the local BD office on <<insert telephone details here>> or e-mail <<insert contact email address here>>.

We confirm that the appropriate regulatory agencies have been informed of these actions.

BD is committed to advancing the world of health. Our primary objectives are patient safety and user safety and providing you with quality products. We apologise for the inconvenience this situation may



BD Switzerland Sàrl  
Terre Bonne Park – A4  
Route de Crassier 17  
1262 Eysins – Switzerland  
Tél: +41 21 556 30 Fax:  
+41 21 556 30 99  
www.BD.com

cause you and thank you in advance for helping BD to resolve this matter as quickly and effectively as possible.

Sincerely,

A handwritten signature in black ink, appearing to read 'L. Darrock'.

Lorna Darrock

Senior Manager Post Market Quality, EMEA



---

## Customer Response Form – MDS-21-4144

### BD SmartSite™ Needle Free Valve

REF & Lot Numbers: See Appendix 1

---

Please read in conjunction with Field Safety Notice MDS-21-4144 and return the completed and signed form as soon as possible or **no later than <<xx August 2021>>** to <<insert fax/email address here>>.

**By signing below, you confirm this notice has been read, understood and that all recommended actions have been implemented as required.**

<b>Account/Organisation Name:</b>	
<b>Department (if applicable):</b>	
<b>Address:</b>	
<b>Postcode:</b>	<b>City:</b>
<b>Contact Name:</b>	
<b>Job Title:</b>	
<b>Contact Telephone Number:</b>	<b>Contact E-mail Address:</b>
<b>Signature:</b>	<b>Date:</b>

*This form must be returned to BD before this action can be considered closed for your account.*



## Appendix 1 – Impacted Product Codes / Lot numbers

Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
10010511	20035964	10885403233470	13-Mar-20	13-Mar-23
10010912	20015157	10885403233531	07-Jan-20	07-Jan-23
10010912	20025207	10885403233531	05-Feb-20	05-Feb-23
10010912	20025243	10885403233531	05-Feb-20	05-Feb-23
10010912	20025244	10885403233531	05-Feb-20	05-Feb-23
10010912	20025245	10885403233531	06-Feb-20	06-Feb-23
10010912	20025246	10885403233531	06-Feb-20	06-Feb-23
10010912	20106552	10885403233531	20-Oct-20	20-Oct-23
10010912	20115468	10885403233531	02-Nov-20	02-Nov-23
10010912	20115469	10885403233531	02-Nov-20	02-Nov-23
10010912	21039510	10885403233531	03-Mar-21	03-Mar-24
10010983	20015758	10885403233579	10-Jan-20	10-Jan-23
10010983	20055744	10885403233579	28-May-20	28-May-23
10010983	20055745	10885403233579	29-May-20	29-May-23
10010983	20075637	10885403233579	07-Jul-20	07-Jul-23
10010983	20075638	10885403233579	07-Jul-20	07-Jul-23
10010983	20075639	10885403233579	07-Jul-20	07-Jul-23
10010983	20075673	10885403233579	07-Jul-20	07-Jul-23
10010983	20096017	10885403233579	09-Sep-20	09-Sep-23
10010983	20096018	10885403233579	09-Sep-20	09-Sep-23
10010983	21016323	10885403233579	22-Jan-21	22-Jan-24
10011087	20125462	10885403233593	10-Dec-20	10-Dec-23
10011253	20016553	10885403233609	19-Jan-20	19-Jan-23
10011253	20056551	10885403233609	23-May-20	23-May-23
10012564	20026597	10885403233746	20-Feb-20	20-Feb-23
10012564	20026598	10885403233746	20-Feb-20	20-Feb-23
10012564	20096477	10885403233746	14-Sep-20	14-Sep-23
10012564	20106218	10885403233746	23-Oct-20	23-Oct-23
10012564	20115237	10885403233746	02-Nov-20	02-Nov-23
10012564	20126022	10885403233746	11-Dec-20	11-Dec-23
10012564	20126023	10885403233746	11-Dec-20	11-Dec-23
10013902	20076494	10885403233814	23-Jul-20	23-Jul-23
10014692	20075342	10885403233869	03-Jul-20	03-Jul-23
10014692	21026223	10885403233869	26-Feb-21	26-Feb-24
10015012	20055717	10885403233951	28-May-20	28-May-23
10015012	20066810	10885403233951	20-Jun-20	20-Jun-23
10015012	20066818	10885403233951	21-Jun-20	21-Jun-23
10015012	20066819	10885403233951	22-Jun-20	22-Jun-23
10015012	20077318	10885403233951	28-Jul-20	28-Jul-23
10015012	20077321	10885403233951	28-Jul-20	28-Jul-23
10015012	20077322	10885403233951	29-Jul-20	29-Jul-23
10015012	20077522	10885403233951	29-Jul-20	29-Jul-23
10015012	20077523	10885403233951	30-Jul-20	30-Jul-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
10015012	20077541	10885403233951	30-Jul-20	30-Jul-23
10015012	20086584	10885403233951	14-Aug-20	14-Aug-23
10015012	20106957	10885403233951	26-Oct-20	26-Oct-23
10015012	20106958	10885403233951	26-Oct-20	26-Oct-23
10015012	20106959	10885403233951	27-Oct-20	27-Oct-23
10015488	20075341	10885403234026	03-Jul-20	03-Jul-23
10015488	21016571	10885403234026	28-Jan-21	28-Jan-24
10037031	20065590	07613203010399	05-Jun-20	05-Jun-23
10037031	20065591	07613203010399	06-Jun-20	06-Jun-23
10037031	20075066	07613203010399	01-Jul-20	01-Jul-23
10037031	20105945	07613203010399	07-Oct-20	07-Oct-23
10037031	20105946	07613203010399	07-Oct-20	07-Oct-23
10037031	20105947	07613203010399	08-Oct-20	08-Oct-23
10037031	20106408	07613203010399	13-Oct-20	13-Oct-23
10037031	21015533	07613203010399	08-Jan-21	08-Jan-24
10037031	21015552	07613203010399	08-Jan-21	08-Jan-24
10037032	20106612	10885403244056	16-Oct-20	16-Oct-23
10037032	20115540	10885403244056	18-Nov-20	18-Nov-23
10037032	21046676	10885403244056	29-Apr-21	29-Apr-24
10798696	20026497	10885403234354	25-Feb-20	25-Feb-23
10798703	20025590	10885403234378	10-Feb-20	10-Feb-23
10798703	20076421	10885403234378	13-Jul-20	13-Jul-23
10798703	21015406	10885403234378	07-Jan-21	07-Jan-24
10798703	21015407	10885403234378	07-Jan-21	07-Jan-24
10817920	20025114	10885403234521	01-Feb-20	01-Feb-23
10821753	20026226	07613203019187	15-Feb-20	15-Feb-23
10821753	20026227	07613203019187	17-Feb-20	17-Feb-23
10821753	20026322	07613203019187	17-Feb-20	17-Feb-23
10821753	20026382	07613203019187	18-Feb-20	18-Feb-23
11607787	20017423	10885403234781	31-Jan-20	31-Jan-23
11613191	20096285	10885403276026	11-Sep-20	11-Sep-23
2000E	20015128	07613203011938	07-Jan-20	07-Jan-23
2000E	20015129	07613203011938	07-Jan-20	07-Jan-23
2000E	20015130	07613203011938	08-Jan-20	08-Jan-23
2000E	20015131	07613203011938	08-Jan-20	08-Jan-23
2000E	20015132	07613203011938	08-Jan-20	08-Jan-23
2000E	20015133	07613203011938	08-Jan-20	08-Jan-23
2000E	20015134	07613203011938	08-Jan-20	08-Jan-23
2000E	20015257	07613203011938	10-Jan-20	10-Jan-23
2000E	20015258	07613203011938	10-Jan-20	10-Jan-23
2000E	20015259	07613203011938	10-Jan-20	10-Jan-23
2000E	20015260	07613203011938	10-Jan-20	10-Jan-23
2000E	20015261	07613203011938	10-Jan-20	10-Jan-23
2000E	20015908	07613203011938	10-Jan-20	10-Jan-23
2000E	20015909	07613203011938	10-Jan-20	10-Jan-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E	20015910	07613203011938	11-Jan-20	11-Jan-23
2000E	20015911	07613203011938	11-Jan-20	11-Jan-23
2000E	20016199	07613203011938	15-Jan-20	15-Jan-23
2000E	20016200	07613203011938	22-Jan-20	22-Jan-23
2000E	20016201	07613203011938	24-Jan-20	24-Jan-23
2000E	20016202	07613203011938	24-Jan-20	24-Jan-23
2000E	20016545	07613203011938	24-Jan-20	24-Jan-23
2000E	20016546	07613203011938	24-Jan-20	24-Jan-23
2000E	20016547	07613203011938	24-Jan-20	24-Jan-23
2000E	20016548	07613203011938	24-Jan-20	24-Jan-23
2000E	20016549	07613203011938	25-Jan-20	25-Jan-23
2000E	20016872	07613203011938	25-Jan-20	25-Jan-23
2000E	20016873	07613203011938	25-Jan-20	25-Jan-23
2000E	20016874	07613203011938	25-Jan-20	25-Jan-23
2000E	20016876	07613203011938	25-Jan-20	25-Jan-23
2000E	20017016	07613203011938	28-Jan-20	28-Jan-23
2000E	20017017	07613203011938	28-Jan-20	28-Jan-23
2000E	20017018	07613203011938	28-Jan-20	28-Jan-23
2000E	20017019	07613203011938	28-Jan-20	28-Jan-23
2000E	20017077	07613203011938	28-Jan-20	28-Jan-23
2000E	20017078	07613203011938	29-Jan-20	29-Jan-23
2000E	20017079	07613203011938	29-Jan-20	29-Jan-23
2000E	20017080	07613203011938	29-Jan-20	29-Jan-23
2000E	20017081	07613203011938	29-Jan-20	29-Jan-23
2000E	20017082	07613203011938	29-Jan-20	29-Jan-23
2000E	20017431	07613203011938	31-Jan-20	31-Jan-23
2000E	20017432	07613203011938	31-Jan-20	31-Jan-23
2000E	20017445	07613203011938	31-Jan-20	31-Jan-23
2000E	20017446	07613203011938	31-Jan-20	31-Jan-23
2000E	20017447	07613203011938	31-Jan-20	31-Jan-23
2000E	20025394	07613203011938	06-Feb-20	06-Feb-23
2000E	20025395	07613203011938	06-Feb-20	06-Feb-23
2000E	20025396	07613203011938	06-Feb-20	06-Feb-23
2000E	20025397	07613203011938	06-Feb-20	06-Feb-23
2000E	20025398	07613203011938	06-Feb-20	06-Feb-23
2000E	20025399	07613203011938	06-Feb-20	06-Feb-23
2000E	20025400	07613203011938	06-Feb-20	06-Feb-23
2000E	20025401	07613203011938	06-Feb-20	06-Feb-23
2000E	20025402	07613203011938	09-Feb-20	09-Feb-23
2000E	20025403	07613203011938	09-Feb-20	09-Feb-23
2000E	20025537	07613203011938	10-Feb-20	10-Feb-23
2000E	20025904	07613203011938	11-Feb-20	11-Feb-23
2000E	20025905	07613203011938	11-Feb-20	11-Feb-23
2000E	20025906	07613203011938	12-Feb-20	12-Feb-23
2000E	20025907	07613203011938	12-Feb-20	12-Feb-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E	20025908	07613203011938	13-Feb-20	13-Feb-23
2000E	20026416	07613203011938	17-Feb-20	17-Feb-23
2000E	20026417	07613203011938	17-Feb-20	17-Feb-23
2000E	20026422	07613203011938	17-Feb-20	17-Feb-23
2000E	20026612	07613203011938	19-Feb-20	19-Feb-23
2000E	20026613	07613203011938	19-Feb-20	19-Feb-23
2000E	20026614	07613203011938	19-Feb-20	19-Feb-23
2000E	20026615	07613203011938	19-Feb-20	19-Feb-23
2000E	20026655	07613203011938	25-Feb-20	25-Feb-23
2000E	20026656	07613203011938	25-Feb-20	25-Feb-23
2000E	20026657	07613203011938	22-Feb-20	22-Feb-23
2000E	20026658	07613203011938	22-Feb-20	22-Feb-23
2000E	20035497	07613203011938	06-Mar-20	06-Mar-23
2000E	20035498	07613203011938	06-Mar-20	06-Mar-23
2000E	20035499	07613203011938	06-Mar-20	06-Mar-23
2000E	20035584	07613203011938	06-Mar-20	06-Mar-23
2000E	20035585	07613203011938	07-Mar-20	07-Mar-23
2000E	20035693	07613203011938	07-Mar-20	07-Mar-23
2000E	20035694	07613203011938	07-Mar-20	07-Mar-23
2000E	20035695	07613203011938	07-Mar-20	07-Mar-23
2000E	20035696	07613203011938	07-Mar-20	07-Mar-23
2000E	20035831	07613203011938	09-Mar-20	09-Mar-23
2000E	20035832	07613203011938	09-Mar-20	09-Mar-23
2000E	20035833	07613203011938	09-Mar-20	09-Mar-23
2000E	20035954	07613203011938	10-Mar-20	10-Mar-23
2000E	20035974	07613203011938	10-Mar-20	10-Mar-23
2000E	20035976	07613203011938	10-Mar-20	10-Mar-23
2000E	20035977	07613203011938	10-Mar-20	10-Mar-23
2000E	20035978	07613203011938	11-Mar-20	11-Mar-23
2000E	20035979	07613203011938	11-Mar-20	11-Mar-23
2000E	20036241	07613203011938	13-Mar-20	13-Mar-23
2000E	20045312	07613203011938	03-Apr-20	03-Apr-23
2000E	20045446	07613203011938	03-Apr-20	03-Apr-23
2000E	20045447	07613203011938	03-Apr-20	03-Apr-23
2000E	20045448	07613203011938	04-Apr-20	04-Apr-23
2000E	20045449	07613203011938	04-Apr-20	04-Apr-23
2000E	20045714	07613203011938	08-Apr-20	08-Apr-23
2000E	20045715	07613203011938	08-Apr-20	08-Apr-23
2000E	20045791	07613203011938	09-Apr-20	09-Apr-23
2000E	20045792	07613203011938	09-Apr-20	09-Apr-23
2000E	20045793	07613203011938	09-Apr-20	09-Apr-23
2000E	20045794	07613203011938	09-Apr-20	09-Apr-23
2000E	20045795	07613203011938	09-Apr-20	09-Apr-23
2000E	20045796	07613203011938	11-Apr-20	11-Apr-23
2000E	20045797	07613203011938	11-Apr-20	11-Apr-23





Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E	20045798	07613203011938	11-Apr-20	11-Apr-23
2000E	20046047	07613203011938	13-Apr-20	13-Apr-23
2000E	20046246	07613203011938	15-Apr-20	15-Apr-23
2000E	20046320	07613203011938	16-Apr-20	16-Apr-23
2000E	20046321	07613203011938	16-Apr-20	16-Apr-23
2000E	20046322	07613203011938	16-Apr-20	16-Apr-23
2000E	20046323	07613203011938	16-Apr-20	16-Apr-23
2000E	20046324	07613203011938	16-Apr-20	16-Apr-23
2000E	20046325	07613203011938	17-Apr-20	17-Apr-23
2000E	20046371	07613203011938	17-Apr-20	17-Apr-23
2000E	20046372	07613203011938	17-Apr-20	17-Apr-23
2000E	20046373	07613203011938	17-Apr-20	17-Apr-23
2000E	20046374	07613203011938	17-Apr-20	17-Apr-23
2000E	20046539	07613203011938	20-Apr-20	20-Apr-23
2000E	20046540	07613203011938	20-Apr-20	20-Apr-23
2000E	20046618	07613203011938	21-Apr-20	21-Apr-23
2000E	20046619	07613203011938	21-Apr-20	21-Apr-23
2000E	20046620	07613203011938	21-Apr-20	21-Apr-23
2000E	20046621	07613203011938	21-Apr-20	21-Apr-23
2000E	20046645	07613203011938	21-Apr-20	21-Apr-23
2000E	20046646	07613203011938	21-Apr-20	21-Apr-23
2000E	20046647	07613203011938	22-Apr-20	22-Apr-23
2000E	20046648	07613203011938	22-Apr-20	22-Apr-23
2000E	20046649	07613203011938	22-Apr-20	22-Apr-23
2000E	20085924	07613203011938	11-Aug-20	11-Aug-23
2000E	20085925	07613203011938	11-Aug-20	11-Aug-23
2000E	20086110	07613203011938	12-Aug-20	12-Aug-23
2000E	20086111	07613203011938	12-Aug-20	12-Aug-23
2000E	20086573	07613203011938	11-Aug-20	11-Aug-23
2000E	20087250	07613203011938	26-Aug-20	26-Aug-23
2000E	20095322	07613203011938	07-Sep-20	07-Sep-23
2000E	20095323	07613203011938	07-Sep-20	07-Sep-23
2000E	20095557	07613203011938	08-Sep-20	08-Sep-23
2000E	20095558	07613203011938	08-Sep-20	08-Sep-23
2000E	20105992	07613203011938	15-Oct-20	15-Oct-23
2000E	20106349	07613203011938	15-Oct-20	15-Oct-23
2000E	20106350	07613203011938	16-Oct-20	16-Oct-23
2000E	20106351	07613203011938	16-Oct-20	16-Oct-23
2000E	20106352	07613203011938	16-Oct-20	16-Oct-23
2000E	20106353	07613203011938	16-Oct-20	16-Oct-23
2000E	20106354	07613203011938	16-Oct-20	16-Oct-23
2000E	20106355	07613203011938	16-Oct-20	16-Oct-23
2000E	20115584	07613203011938	03-Nov-20	03-Nov-23
2000E	20115585	07613203011938	03-Nov-20	03-Nov-23
2000E	20115586	07613203011938	03-Nov-20	03-Nov-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E	20115587	07613203011938	03-Nov-20	03-Nov-23
2000E	20115589	07613203011938	04-Nov-20	04-Nov-23
2000E	20115590	07613203011938	04-Nov-20	04-Nov-23
2000E	20115591	07613203011938	04-Nov-20	04-Nov-23
2000E	20115592	07613203011938	04-Nov-20	04-Nov-23
2000E	20115596	07613203011938	04-Nov-20	04-Nov-23
2000E	20115597	07613203011938	05-Nov-20	05-Nov-23
2000E	20115599	07613203011938	06-Nov-20	06-Nov-23
2000E	20115602	07613203011938	07-Nov-20	07-Nov-23
2000E	20115603	07613203011938	09-Nov-20	09-Nov-23
2000E	20115604	07613203011938	09-Nov-20	09-Nov-23
2000E	20115605	07613203011938	09-Nov-20	09-Nov-23
2000E	20115710	07613203011938	09-Nov-20	09-Nov-23
2000E	20116994	07613203011938	29-Nov-20	29-Nov-23
2000E	20116995	07613203011938	29-Nov-20	29-Nov-23
2000E	20116996	07613203011938	29-Nov-20	29-Nov-23
2000E	20116997	07613203011938	30-Nov-20	30-Nov-23
2000E	20116998	07613203011938	30-Nov-20	30-Nov-23
2000E	20125044	07613203011938	02-Dec-20	02-Dec-23
2000E	20125045	07613203011938	02-Dec-20	02-Dec-23
2000E	20125046	07613203011938	02-Dec-20	02-Dec-23
2000E	20125251	07613203011938	03-Dec-20	03-Dec-23
2000E	20125252	07613203011938	03-Dec-20	03-Dec-23
2000E	20125253	07613203011938	03-Dec-20	03-Dec-23
2000E	20125254	07613203011938	03-Dec-20	03-Dec-23
2000E	20125639	07613203011938	04-Dec-20	04-Dec-23
2000E	20125640	07613203011938	04-Dec-20	04-Dec-23
2000E	20125641	07613203011938	04-Dec-20	04-Dec-23
2000E	20125642	07613203011938	04-Dec-20	04-Dec-23
2000E	20125683	07613203011938	04-Dec-20	04-Dec-23
2000E	20125684	07613203011938	04-Dec-20	04-Dec-23
2000E	20125685	07613203011938	04-Dec-20	04-Dec-23
2000E	21016238	07613203011938	20-Jan-21	20-Jan-24
2000E	21016239	07613203011938	20-Jan-21	20-Jan-24
2000E	21016240	07613203011938	20-Jan-21	20-Jan-24
2000E	21025046	07613203011938	20-Feb-21	20-Feb-24
2000E	21025047	07613203011938	16-Feb-21	16-Feb-24
2000E	21025048	07613203011938	16-Feb-21	16-Feb-24
2000E	21025049	07613203011938	20-Feb-21	20-Feb-24
2000E	21025158	07613203011938	20-Feb-21	20-Feb-24
2000E	21025284	07613203011938	20-Feb-21	20-Feb-24
2000E	21025286	07613203011938	22-Feb-21	22-Feb-24
2000E	21025287	07613203011938	23-Feb-21	23-Feb-24
2000E	21025288	07613203011938	23-Feb-21	23-Feb-24
2000E	21025387	07613203011938	23-Feb-21	23-Feb-24



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E	21025388	07613203011938	24-Feb-21	24-Feb-24
2000E	21025389	07613203011938	24-Feb-21	24-Feb-24
2000E	21035545	07613203011938	06-Mar-21	06-Mar-24
2000E	21035546	07613203011938	03-Mar-21	03-Mar-24
2000E	21035547	07613203011938	06-Mar-21	06-Mar-24
2000E	21035549	07613203011938	06-Mar-21	06-Mar-24
2000E	21035550	07613203011938	06-Mar-21	06-Mar-24
2000E	21035551	07613203011938	07-Mar-21	07-Mar-24
2000E	21035552	07613203011938	07-Mar-21	07-Mar-24
2000E	21035553	07613203011938	07-Mar-21	07-Mar-24
2000E	21035554	07613203011938	08-Mar-21	08-Mar-24
2000E	21035555	07613203011938	08-Mar-21	08-Mar-24
2000E	21035556	07613203011938	08-Mar-21	08-Mar-24
2000E	21035557	07613203011938	08-Mar-21	08-Mar-24
2000E	21035558	07613203011938	09-Mar-21	09-Mar-24
2000E	21035559	07613203011938	09-Mar-21	09-Mar-24
2000E	21035768	07613203011938	09-Mar-21	09-Mar-24
2000E	21035769	07613203011938	11-Mar-21	11-Mar-24
2000E	21035770	07613203011938	11-Mar-21	11-Mar-24
2000E	21035771	07613203011938	11-Mar-21	11-Mar-24
2000E	21035772	07613203011938	12-Mar-21	12-Mar-24
2000E	21035773	07613203011938	12-Mar-21	12-Mar-24
2000E	21035774	07613203011938	12-Mar-21	12-Mar-24
2000E	21036160	07613203011938	17-Mar-21	17-Mar-24
2000E	21036161	07613203011938	17-Mar-21	17-Mar-24
2000E	21036162	07613203011938	18-Mar-21	18-Mar-24
2000E	21036163	07613203011938	18-Mar-21	18-Mar-24
2000E	21036164	07613203011938	19-Mar-21	19-Mar-24
2000E	21036165	07613203011938	19-Mar-21	19-Mar-24
2000E	21036265	07613203011938	22-Mar-21	22-Mar-24
2000E	21036266	07613203011938	24-Mar-21	24-Mar-24
2000E	21036267	07613203011938	24-Mar-21	24-Mar-24
2000E	21036268	07613203011938	24-Mar-21	24-Mar-24
2000E	21036269	07613203011938	26-Mar-21	26-Mar-24
2000E	21036270	07613203011938	26-Mar-21	26-Mar-24
2000E	21036685	07613203011938	26-Mar-21	26-Mar-24
2000E	21036708	07613203011938	26-Mar-21	26-Mar-24
2000E	21036709	07613203011938	27-Mar-21	27-Mar-24
2000E	21036793	07613203011938	29-Mar-21	29-Mar-24
2000E CHINA	20025240	10885403238345	05-Feb-20	05-Feb-23
2000E CHINA	20025241	10885403238345	05-Feb-20	05-Feb-23
2000E CHINA	20025242	10885403238345	05-Feb-20	05-Feb-23
2000E CHINA	20056225	10885403238345	14-May-20	14-May-23
2000E CHINA	20056229	10885403238345	27-May-20	27-May-23
2000E CHINA	20056230	10885403238345	28-May-20	28-May-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E CHINA	20056231	10885403238345	14-May-20	14-May-23
2000E CHINA	21025820	10885403238345	08-Feb-21	08-Feb-24
2000E CHINA	21025823	10885403238345	09-Feb-21	09-Feb-24
2000E CHINA	21025824	10885403238345	09-Feb-21	09-Feb-24
2000E CHINA	21025825	10885403238345	09-Feb-21	09-Feb-24
2000E CHINA	21025829	10885403238345	10-Feb-21	10-Feb-24
2000E-04	20036690	10885403480928	20-Mar-20	20-Mar-23
2000E-04	20036691	10885403480928	20-Mar-20	20-Mar-23
2000E-04	20036692	10885403480928	20-Mar-20	20-Mar-23
2000E-04	20036694	10885403480928	20-Mar-20	20-Mar-23
2000E-04	20036983	10885403480928	25-Mar-20	25-Mar-23
2000E-04	20046920	10885403480928	25-Apr-20	25-Apr-23
2000E-04	20046921	10885403480928	25-Apr-20	25-Apr-23
2000E-04	20046922	10885403480928	25-Apr-20	25-Apr-23
2000E-04	20055220	10885403480928	04-May-20	04-May-23
2000E-04	20055221	10885403480928	04-May-20	04-May-23
2000E-04	20055222	10885403480928	04-May-20	04-May-23
2000E-04	20055223	10885403480928	04-May-20	04-May-23
2000E-04	20055379	10885403480928	04-May-20	04-May-23
2000E-04	20055380	10885403480928	04-May-20	04-May-23
2000E-04	20055381	10885403480928	05-May-20	05-May-23
2000E-04	20055382	10885403480928	05-May-20	05-May-23
2000E-04	20055503	10885403480928	06-May-20	06-May-23
2000E-04	20055504	10885403480928	08-May-20	08-May-23
2000E-04	20055505	10885403480928	08-May-20	08-May-23
2000E-04	20055506	10885403480928	08-May-20	08-May-23
2000E-04	20055507	10885403480928	08-May-20	08-May-23
2000E-04	20055508	10885403480928	08-May-20	08-May-23
2000E-04	20055880	10885403480928	08-May-20	08-May-23
2000E-04	20055881	10885403480928	09-May-20	09-May-23
2000E-04	20055882	10885403480928	09-May-20	09-May-23
2000E-04	20055883	10885403480928	09-May-20	09-May-23
2000E-04	20055884	10885403480928	12-May-20	12-May-23
2000E-04	20055885	10885403480928	12-May-20	12-May-23
2000E-04	20055886	10885403480928	20-May-20	20-May-23
2000E-04	20055887	10885403480928	20-May-20	20-May-23
2000E-04	20055888	10885403480928	20-May-20	20-May-23
2000E-04	20055889	10885403480928	21-May-20	21-May-23
2000E-04	20055890	10885403480928	21-May-20	21-May-23
2000E-04	20065100	10885403480928	01-Jun-20	01-Jun-23
2000E-04	20065101	10885403480928	01-Jun-20	01-Jun-23
2000E-04	20065102	10885403480928	01-Jun-20	01-Jun-23
2000E-04	20065103	10885403480928	02-Jun-20	02-Jun-23
2000E-04	20065104	10885403480928	01-Jun-20	01-Jun-23
2000E-04	20065105	10885403480928	02-Jun-20	02-Jun-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E-04	20065106	10885403480928	02-Jun-20	02-Jun-23
2000E-04	20065107	10885403480928	02-Jun-20	02-Jun-23
2000E-04	20065108	10885403480928	02-Jun-20	02-Jun-23
2000E-04	20065109	10885403480928	02-Jun-20	02-Jun-23
2000E-04	20065110	10885403480928	12-Jun-20	12-Jun-23
2000E-04	20065111	10885403480928	12-Jun-20	12-Jun-23
2000E-04	20065986	10885403480928	12-Jun-20	12-Jun-23
2000E-04	20065987	10885403480928	12-Jun-20	12-Jun-23
2000E-04	20065988	10885403480928	13-Jun-20	13-Jun-23
2000E-04	20065989	10885403480928	19-Jun-20	19-Jun-23
2000E-04	20066898	10885403480928	22-Jun-20	22-Jun-23
2000E-04	20066936	10885403480928	24-Jun-20	24-Jun-23
2000E-04	20066937	10885403480928	25-Jun-20	25-Jun-23
2000E-04	20066938	10885403480928	25-Jun-20	25-Jun-23
2000E-04	20066939	10885403480928	25-Jun-20	25-Jun-23
2000E-04	20066940	10885403480928	25-Jun-20	25-Jun-23
2000E-04	20067172	10885403480928	26-Jun-20	26-Jun-23
2000E-04	20075176	10885403480928	02-Jul-20	02-Jul-23
2000E-04	20075177	10885403480928	02-Jul-20	02-Jul-23
2000E-04	20075481	10885403480928	07-Jul-20	07-Jul-23
2000E-04	20075482	10885403480928	07-Jul-20	07-Jul-23
2000E-04	20075483	10885403480928	08-Jul-20	08-Jul-23
2000E-04	20075484	10885403480928	08-Jul-20	08-Jul-23
2000E-04	20075623	10885403480928	08-Jul-20	08-Jul-23
2000E-04	20075624	10885403480928	08-Jul-20	08-Jul-23
2000E-04	20075625	10885403480928	08-Jul-20	08-Jul-23
2000E-04	20075793	10885403480928	09-Jul-20	09-Jul-23
2000E-04	20075794	10885403480928	09-Jul-20	09-Jul-23
2000E-04	20076070	10885403480928	10-Jul-20	10-Jul-23
2000E-04	20076071	10885403480928	10-Jul-20	10-Jul-23
2000E-04	20076486	10885403480928	15-Jul-20	15-Jul-23
2000E-04	20076487	10885403480928	15-Jul-20	15-Jul-23
2000E-04	20076488	10885403480928	15-Jul-20	15-Jul-23
2000E-04	20076489	10885403480928	16-Jul-20	16-Jul-23
2000E-04	20076490	10885403480928	16-Jul-20	16-Jul-23
2000E-04	20076491	10885403480928	16-Jul-20	16-Jul-23
2000E-04	20076576	10885403480928	16-Jul-20	16-Jul-23
2000E-04	20076577	10885403480928	16-Jul-20	16-Jul-23
2000E-04	20076578	10885403480928	16-Jul-20	16-Jul-23
2000E-04	20076686	10885403480928	17-Jul-20	17-Jul-23
2000E-04	20076687	10885403480928	17-Jul-20	17-Jul-23
2000E-04	20076769	10885403480928	17-Jul-20	17-Jul-23
2000E-04	20076770	10885403480928	17-Jul-20	17-Jul-23
2000E-04	20076771	10885403480928	17-Jul-20	17-Jul-23
2000E-04	20076772	10885403480928	18-Jul-20	18-Jul-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E-04	20076773	10885403480928	18-Jul-20	18-Jul-23
2000E-04	20076774	10885403480928	18-Jul-20	18-Jul-23
2000E-04	20077355	10885403480928	24-Jul-20	24-Jul-23
2000E-04	20077356	10885403480928	24-Jul-20	24-Jul-23
2000E-04	20077357	10885403480928	27-Jul-20	27-Jul-23
2000E-04	20077358	10885403480928	27-Jul-20	27-Jul-23
2000E-04	20077359	10885403480928	27-Jul-20	27-Jul-23
2000E-04	20077360	10885403480928	28-Jul-20	28-Jul-23
2000E-04	20077361	10885403480928	28-Jul-20	28-Jul-23
2000E-04	20077362	10885403480928	28-Jul-20	28-Jul-23
2000E-04	20077363	10885403480928	29-Jul-20	29-Jul-23
2000E-04	20077364	10885403480928	29-Jul-20	29-Jul-23
2000E-04	20085089	10885403480928	01-Aug-20	01-Aug-23
2000E-04	20085090	10885403480928	01-Aug-20	01-Aug-23
2000E-04	20085091	10885403480928	01-Aug-20	01-Aug-23
2000E-04	20085092	10885403480928	03-Aug-20	03-Aug-23
2000E-04	20085289	10885403480928	03-Aug-20	03-Aug-23
2000E-04	20085290	10885403480928	03-Aug-20	03-Aug-23
2000E-04	20085291	10885403480928	03-Aug-20	03-Aug-23
2000E-04	20085366	10885403480928	03-Aug-20	03-Aug-23
2000E-04	20087446	10885403480928	31-Aug-20	31-Aug-23
2000E-04	20087447	10885403480928	31-Aug-20	31-Aug-23
2000E-04	20087448	10885403480928	31-Aug-20	31-Aug-23
2000E-04	20095317	10885403480928	01-Sep-20	01-Sep-23
2000E-04	20095318	10885403480928	04-Sep-20	04-Sep-23
2000E-04	20105086	10885403480928	01-Oct-20	01-Oct-23
2000E-04	20105087	10885403480928	01-Oct-20	01-Oct-23
2000E-04	20105088	10885403480928	01-Oct-20	01-Oct-23
2000E-04	20105089	10885403480928	01-Oct-20	01-Oct-23
2000E-04	20105571	10885403480928	02-Oct-20	02-Oct-23
2000E-04	20105572	10885403480928	02-Oct-20	02-Oct-23
2000E-04	20105993	10885403480928	17-Oct-20	17-Oct-23
2000E-04	20106356	10885403480928	17-Oct-20	17-Oct-23
2000E-04	20106357	10885403480928	19-Oct-20	19-Oct-23
2000E-04	20106358	10885403480928	19-Oct-20	19-Oct-23
2000E-04	20106359	10885403480928	29-Oct-20	29-Oct-23
2000E-04	20106360	10885403480928	29-Oct-20	29-Oct-23
2000E-04	20106361	10885403480928	29-Oct-20	29-Oct-23
2000E-04	20106362	10885403480928	29-Oct-20	29-Oct-23
2000E-04	20106363	10885403480928	29-Oct-20	29-Oct-23
2000E-04	20106364	10885403480928	29-Oct-20	29-Oct-23
2000E-04	20106402	10885403480928	30-Oct-20	30-Oct-23
2000E-04	20106403	10885403480928	30-Oct-20	30-Oct-23
2000E-04	20106404	10885403480928	30-Oct-20	30-Oct-23
2000E-04	20106405	10885403480928	30-Oct-20	30-Oct-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E-04	20106406	10885403480928	30-Oct-20	30-Oct-23
2000E-04	20106407	10885403480928	30-Oct-20	30-Oct-23
2000E-04	20106727	10885403480928	30-Oct-20	30-Oct-23
2000E-04	20115181	10885403480928	17-Nov-20	17-Nov-23
2000E-04	20115182	10885403480928	02-Nov-20	02-Nov-23
2000E-04	20115183	10885403480928	02-Nov-20	02-Nov-23
2000E-04	20115184	10885403480928	02-Nov-20	02-Nov-23
2000E-04	20115186	10885403480928	02-Nov-20	02-Nov-23
2000E-04	20115203	10885403480928	18-Nov-20	18-Nov-23
2000E-04	20115205	10885403480928	02-Nov-20	02-Nov-23
2000E-04	20115207	10885403480928	03-Nov-20	03-Nov-23
2000E-04	20115208	10885403480928	03-Nov-20	03-Nov-23
2000E-04	20115533	10885403480928	17-Nov-20	17-Nov-23
2000E-04	20115534	10885403480928	22-Nov-20	22-Nov-23
2000E-04	20116576	10885403480928	22-Nov-20	22-Nov-23
2000E-04	20116577	10885403480928	22-Nov-20	22-Nov-23
2000E-04	20116578	10885403480928	23-Nov-20	23-Nov-23
2000E-04	20116579	10885403480928	24-Nov-20	24-Nov-23
2000E-04	20116580	10885403480928	24-Nov-20	24-Nov-23
2000E-04	20116581	10885403480928	24-Nov-20	24-Nov-23
2000E-04	20116582	10885403480928	24-Nov-20	24-Nov-23
2000E-04	20116583	10885403480928	24-Nov-20	24-Nov-23
2000E-04	20116584	10885403480928	24-Nov-20	24-Nov-23
2000E-04	20116585	10885403480928	24-Nov-20	24-Nov-23
2000E-04	20116586	10885403480928	25-Nov-20	25-Nov-23
2000E-04	20116587	10885403480928	25-Nov-20	25-Nov-23
2000E-04	20116588	10885403480928	25-Nov-20	25-Nov-23
2000E-04	20116589	10885403480928	25-Nov-20	25-Nov-23
2000E-04	20116590	10885403480928	26-Nov-20	26-Nov-23
2000E-04	20116591	10885403480928	26-Nov-20	26-Nov-23
2000E-04	20116592	10885403480928	26-Nov-20	26-Nov-23
2000E-04	20116593	10885403480928	26-Nov-20	26-Nov-23
2000E-04	20116594	10885403480928	26-Nov-20	26-Nov-23
2000E-04	20116595	10885403480928	26-Nov-20	26-Nov-23
2000E-04	20126194	10885403480928	15-Dec-20	15-Dec-23
2000E-04	20126195	10885403480928	15-Dec-20	15-Dec-23
2000E-04	20126196	10885403480928	15-Dec-20	15-Dec-23
2000E-04	20126197	10885403480928	16-Dec-20	16-Dec-23
2000E-04	20126198	10885403480928	16-Dec-20	16-Dec-23
2000E-04	20126249	10885403480928	16-Dec-20	16-Dec-23
2000E-04	20126250	10885403480928	16-Dec-20	16-Dec-23
2000E-04	20126253	10885403480928	16-Dec-20	16-Dec-23
2000E-04	20126254	10885403480928	17-Dec-20	17-Dec-23
2000E-04	20126255	10885403480928	17-Dec-20	17-Dec-23
2000E-04	21015660	10885403480928	08-Jan-21	08-Jan-24



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E-04	21015661	10885403480928	08-Jan-21	08-Jan-24
2000E-04	21015724	10885403480928	09-Jan-21	09-Jan-24
2000E-04	21015725	10885403480928	09-Jan-21	09-Jan-24
2000E-04	21015753	10885403480928	11-Jan-21	11-Jan-24
2000E-04	21015754	10885403480928	11-Jan-21	11-Jan-24
2000E-04	21015775	10885403480928	11-Jan-21	11-Jan-24
2000E-04	21016152	10885403480928	25-Jan-21	25-Jan-24
2000E-04	21035775	10885403480928	14-Mar-21	14-Mar-24
2000E-04	21035776	10885403480928	14-Mar-21	14-Mar-24
2000E-04	21035777	10885403480928	14-Mar-21	14-Mar-24
2000E-04	21035778	10885403480928	14-Mar-21	14-Mar-24
2000E-04	21036842	10885403480928	30-Mar-21	30-Mar-24
2000E7D	1015511	07613203013611	11-Mar-20	30-Nov-22
2000E7D	1016651	07613203013611	02-Apr-20	28-Feb-23
2000E7D	1016663	07613203013611	06-Apr-20	31-Mar-23
2000E7D	1016664	07613203013611	08-Apr-20	31-Mar-23
2000E7D	1016761	07613203013611	11-Apr-20	31-Mar-23
2000E7D	1016762	07613203013611	17-Apr-20	31-Mar-23
2000E7D	1017045	07613203013611	06-May-20	31-Mar-23
2000E7D	1017152	07613203013611	08-May-20	30-Apr-23
2000E7D	1017153	07613203013611	12-May-20	30-Apr-23
2000E7D	1017154	07613203013611	15-May-20	30-Apr-23
2000E7D	1017306	07613203013611	02-Jun-20	30-Apr-23
2000E7D	1017679	07613203013611	16-Jun-20	31-May-23
2000E7D	1017680	07613203013611	18-Jun-20	31-May-23
2000E7D	1017681	07613203013611	19-Jun-20	31-May-23
2000E7D	1017964	07613203013611	06-Jul-20	30-Jun-23
2000E7D	1018010	07613203013611	25-Aug-20	30-Jun-23
2000E7D	1018316	07613203013611	01-Sep-20	31-Jul-23
2000E7D	1018317	07613203013611	01-Sep-20	31-Jul-23
2000E7D	1018575	07613203013611	05-Oct-20	31-Aug-23
2000E7D	1018576	07613203013611	01-Oct-20	31-Aug-23
2000E7D	1018577	07613203013611	05-Oct-20	31-Aug-23
2000E7D	1018578	07613203013611	03-Dec-20	31-Aug-23
2000E7D	1018579	07613203013611	00-Jan-00	31-Aug-23
2000E7D	1018583	07613203013611	10-Dec-20	31-Aug-23
2000E7D	1018584	07613203013611	10-Dec-20	31-Aug-23
2000E7D	1019392	07613203013611	15-Dec-20	30-Nov-23
2000E7D	1019393	07613203013611	10-Dec-20	30-Nov-23
2000E7D	1019582	07613203013611	04-Jan-21	30-Nov-23
2000E7D	1019583	07613203013611	07-Jan-21	30-Nov-23
2000E7D	1019584	07613203013611	11-Jan-21	30-Nov-23
2000E7D	1019886	07613203013611	02-Feb-21	31-Dec-23
2000E7D	1019892	07613203013611	08-Feb-21	31-Jan-24
2000E7D	1019925	07613203013611	08-Feb-21	31-Jan-24





Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2000E7D	1020131	07613203013611	22-Feb-21	31-Jan-24
2000E7D	1020967	07613203013611	14-May-21	30-Apr-24
2000E7D	1021008	07613203013611	17-May-21	30-Apr-24
2000E7D	1021015	07613203013611	24-May-21	30-Apr-24
20019E	20016104	07613203012010	15-Jan-20	15-Jan-23
20019E	20016105	07613203012010	16-Jan-20	16-Jan-23
20019E	20016190	07613203012010	16-Jan-20	16-Jan-23
20019E	20016191	07613203012010	17-Jan-20	17-Jan-23
20019E	20016192	07613203012010	20-Jan-20	20-Jan-23
20019E	20016193	07613203012010	20-Jan-20	20-Jan-23
20019E	20016194	07613203012010	20-Jan-20	20-Jan-23
20019E	20016315	07613203012010	21-Jan-20	21-Jan-23
20019E	20016316	07613203012010	23-Jan-20	23-Jan-23
20019E	20016317	07613203012010	23-Jan-20	23-Jan-23
20019E	20016442	07613203012010	23-Jan-20	23-Jan-23
20019E	20016443	07613203012010	24-Jan-20	24-Jan-23
20019E	20016444	07613203012010	24-Jan-20	24-Jan-23
20019E	20016445	07613203012010	24-Jan-20	24-Jan-23
20019E	20016446	07613203012010	24-Jan-20	24-Jan-23
20019E	20025021	07613203012010	06-Feb-20	06-Feb-23
20019E	20025022	07613203012010	07-Feb-20	07-Feb-23
20019E	20025023	07613203012010	07-Feb-20	07-Feb-23
20019E	20025024	07613203012010	08-Feb-20	08-Feb-23
20019E	20035493	07613203012010	05-Mar-20	05-Mar-23
20019E	20035494	07613203012010	06-Mar-20	06-Mar-23
20019E	20035495	07613203012010	10-Mar-20	10-Mar-23
20019E	20035496	07613203012010	13-Mar-20	13-Mar-23
20019E	20035577	07613203012010	14-Mar-20	14-Mar-23
20019E	20035578	07613203012010	17-Mar-20	17-Mar-23
20019E	20045836	07613203012010	11-Apr-20	11-Apr-23
20019E	20045837	07613203012010	13-Apr-20	13-Apr-23
20019E	20045838	07613203012010	13-Apr-20	13-Apr-23
20019E	20045839	07613203012010	14-Apr-20	14-Apr-23
20019E	20045840	07613203012010	14-Apr-20	14-Apr-23
20019E	20045841	07613203012010	16-Apr-20	16-Apr-23
20019E	20045842	07613203012010	16-Apr-20	16-Apr-23
20019E	20045843	07613203012010	17-Apr-20	17-Apr-23
20019E	20055076	07613203012010	04-May-20	04-May-23
20019E	20055077	07613203012010	05-May-20	05-May-23
20019E	20055078	07613203012010	05-May-20	05-May-23
20019E	20055079	07613203012010	06-May-20	06-May-23
20019E	20055080	07613203012010	06-May-20	06-May-23
20019E	20055329	07613203012010	07-May-20	07-May-23
20019E	20055330	07613203012010	07-May-20	07-May-23
20019E	20055331	07613203012010	08-May-20	08-May-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20019E	20055332	07613203012010	09-May-20	09-May-23
20019E	20066104	07613203012010	11-Jun-20	11-Jun-23
20019E	20066105	07613203012010	12-Jun-20	12-Jun-23
20019E	20066106	07613203012010	12-Jun-20	12-Jun-23
20019E	20066107	07613203012010	12-Jun-20	12-Jun-23
20019E	20066108	07613203012010	12-Jun-20	12-Jun-23
20019E	20066109	07613203012010	15-Jun-20	15-Jun-23
20019E	20066228	07613203012010	15-Jun-20	15-Jun-23
20019E	20066229	07613203012010	16-Jun-20	16-Jun-23
20019E	20066230	07613203012010	16-Jun-20	16-Jun-23
20019E	20066231	07613203012010	16-Jun-20	16-Jun-23
20019E	20066446	07613203012010	22-Jun-20	22-Jun-23
20019E	20066447	07613203012010	22-Jun-20	22-Jun-23
20019E	20066448	07613203012010	23-Jun-20	23-Jun-23
20019E	20066717	07613203012010	23-Jun-20	23-Jun-23
20019E	20066718	07613203012010	24-Jun-20	24-Jun-23
20019E	20066719	07613203012010	24-Jun-20	24-Jun-23
20019E	20095132	07613203012010	17-Sep-20	17-Sep-23
20019E	20095133	07613203012010	17-Sep-20	17-Sep-23
20019E	20095134	07613203012010	17-Sep-20	17-Sep-23
20019E	20095135	07613203012010	17-Sep-20	17-Sep-23
20019E	20095136	07613203012010	17-Sep-20	17-Sep-23
20019E	20095137	07613203012010	21-Sep-20	21-Sep-23
20019E	20095138	07613203012010	21-Sep-20	21-Sep-23
20019E	20105084	07613203012010	14-Oct-20	14-Oct-23
20019E	20105085	07613203012010	14-Oct-20	14-Oct-23
20019E	20106229	07613203012010	16-Oct-20	16-Oct-23
20019E	20106230	07613203012010	15-Oct-20	15-Oct-23
20019E	20106231	07613203012010	16-Oct-20	16-Oct-23
20019E	20106570	07613203012010	16-Oct-20	16-Oct-23
20019E	20115014	07613203012010	05-Nov-20	05-Nov-23
20019E	20115015	07613203012010	05-Nov-20	05-Nov-23
20019E	20115016	07613203012010	06-Nov-20	06-Nov-23
20019E	20115017	07613203012010	07-Nov-20	07-Nov-23
20019E	20115018	07613203012010	08-Nov-20	08-Nov-23
20019E	20115019	07613203012010	08-Nov-20	08-Nov-23
20019E	20115020	07613203012010	09-Nov-20	09-Nov-23
20019E	20115021	07613203012010	09-Nov-20	09-Nov-23
20019E	20115022	07613203012010	10-Nov-20	10-Nov-23
20019E	20115023	07613203012010	10-Nov-20	10-Nov-23
20019E	20125948	07613203012010	08-Dec-20	08-Dec-23
20019E	20125949	07613203012010	07-Dec-20	07-Dec-23
20019E	20125950	07613203012010	08-Dec-20	08-Dec-23
20019E	20125972	07613203012010	09-Dec-20	09-Dec-23
20019E	20125973	07613203012010	09-Dec-20	09-Dec-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20019E	20125974	07613203012010	09-Dec-20	09-Dec-23
20019E	20125975	07613203012010	09-Dec-20	09-Dec-23
20019E	20125976	07613203012010	10-Dec-20	10-Dec-23
20019E	20125977	07613203012010	10-Dec-20	10-Dec-23
20019E	20125978	07613203012010	11-Dec-20	11-Dec-23
20019E	20125979	07613203012010	11-Dec-20	11-Dec-23
20019E	20125980	07613203012010	11-Dec-20	11-Dec-23
20019E	20125981	07613203012010	12-Dec-20	12-Dec-23
20019E	20125982	07613203012010	12-Dec-20	12-Dec-23
20019E	20125983	07613203012010	12-Dec-20	12-Dec-23
20019E	20125984	07613203012010	14-Dec-20	14-Dec-23
20019E	20125985	07613203012010	14-Dec-20	14-Dec-23
20019E	20125986	07613203012010	14-Dec-20	14-Dec-23
20019E	20126526	07613203012010	22-Dec-20	22-Dec-23
20019E	20126527	07613203012010	22-Dec-20	22-Dec-23
20019E	20126528	07613203012010	23-Dec-20	23-Dec-23
20019E	20126529	07613203012010	23-Dec-20	23-Dec-23
20019E	20126552	07613203012010	23-Dec-20	23-Dec-23
20019E	20126553	07613203012010	23-Dec-20	23-Dec-23
20019E	20126554	07613203012010	23-Dec-20	23-Dec-23
20019E	21015305	07613203012010	05-Jan-21	05-Jan-24
20019E	21015306	07613203012010	05-Jan-21	05-Jan-24
20019E	21016393	07613203012010	25-Jan-21	25-Jan-24
20019E	21016394	07613203012010	25-Jan-21	25-Jan-24
20019E	21025279	07613203012010	03-Feb-21	03-Feb-24
20019E	21025280	07613203012010	04-Feb-21	04-Feb-24
20019E	21025281	07613203012010	04-Feb-21	04-Feb-24
20019E	21025282	07613203012010	04-Feb-21	04-Feb-24
20019E	21025283	07613203012010	04-Feb-21	04-Feb-24
20019E	21026323	07613203012010	20-Feb-21	20-Feb-24
20019E	21026324	07613203012010	20-Feb-21	20-Feb-24
20019E	21026325	07613203012010	20-Feb-21	20-Feb-24
20019E	21026326	07613203012010	21-Feb-21	21-Feb-24
20019E	21026327	07613203012010	21-Feb-21	21-Feb-24
20019E	21026377	07613203012010	21-Feb-21	21-Feb-24
20019E	21026378	07613203012010	21-Feb-21	21-Feb-24
20019E	21026379	07613203012010	21-Feb-21	21-Feb-24
20019E	21026545	07613203012010	25-Feb-21	25-Feb-24
20019E	21026546	07613203012010	25-Feb-21	25-Feb-24
20019E	21026547	07613203012010	26-Feb-21	26-Feb-24
20019E	21036695	07613203012010	27-Mar-21	27-Mar-24
20019E	21036696	07613203012010	27-Mar-21	27-Mar-24
20019E7D	20015114	07613203014854	02-Jan-20	02-Jan-23
20019E7D	20015115	07613203014854	02-Jan-20	02-Jan-23
20019E7D	20015116	07613203014854	29-Jan-20	29-Jan-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20019E7D	20015117	07613203014854	29-Jan-20	29-Jan-23
20019E7D	20015118	07613203014854	29-Jan-20	29-Jan-23
20019E7D	20025813	07613203014854	14-Feb-20	14-Feb-23
20019E7D	20025814	07613203014854	16-Feb-20	16-Feb-23
20019E7D	20025815	07613203014854	17-Feb-20	17-Feb-23
20019E7D	20025816	07613203014854	17-Feb-20	17-Feb-23
20019E7D	20025817	07613203014854	19-Feb-20	19-Feb-23
20019E7D	20025818	07613203014854	19-Feb-20	19-Feb-23
20019E7D	20025819	07613203014854	20-Feb-20	20-Feb-23
20019E7D	20036187	07613203014854	17-Mar-20	17-Mar-23
20019E7D	20036188	07613203014854	18-Mar-20	18-Mar-23
20019E7D	20036259	07613203014854	18-Mar-20	18-Mar-23
20019E7D	20036260	07613203014854	19-Mar-20	19-Mar-23
20019E7D	20036261	07613203014854	20-Mar-20	20-Mar-23
20019E7D	20036262	07613203014854	21-Mar-20	21-Mar-23
20019E7D	20036327	07613203014854	28-Mar-20	28-Mar-23
20019E7D	20036328	07613203014854	30-Mar-20	30-Mar-23
20019E7D	20036329	07613203014854	31-Mar-20	31-Mar-23
20019E7D	20037133	07613203014854	31-Mar-20	31-Mar-23
20019E7D	20046569	07613203014854	23-Apr-20	23-Apr-23
20019E7D	20046570	07613203014854	23-Apr-20	23-Apr-23
20019E7D	20046615	07613203014854	23-Apr-20	23-Apr-23
20019E7D	20046616	07613203014854	24-Apr-20	24-Apr-23
20019E7D	20046617	07613203014854	24-Apr-20	24-Apr-23
20019E7D	20055967	07613203014854	14-May-20	14-May-23
20019E7D	20055968	07613203014854	14-May-20	14-May-23
20019E7D	20055984	07613203014854	14-May-20	14-May-23
20019E7D	20055985	07613203014854	15-May-20	15-May-23
20019E7D	20055986	07613203014854	16-May-20	16-May-23
20019E7D	20055987	07613203014854	18-May-20	18-May-23
20019E7D	20055988	07613203014854	18-May-20	18-May-23
20019E7D	20055989	07613203014854	20-May-20	20-May-23
20019E7D	20056364	07613203014854	20-May-20	20-May-23
20019E7D	20056365	07613203014854	21-May-20	21-May-23
20019E7D	20056366	07613203014854	21-May-20	21-May-23
20019E7D	20056452	07613203014854	24-May-20	24-May-23
20019E7D	20056453	07613203014854	24-May-20	24-May-23
20019E7D	20056565	07613203014854	24-May-20	24-May-23
20019E7D	20056566	07613203014854	25-May-20	25-May-23
20019E7D	20057064	07613203014854	29-May-20	29-May-23
20019E7D	20075232	07613203014854	03-Jul-20	03-Jul-23
20019E7D	20075235	07613203014854	17-Jul-20	17-Jul-23
20019E7D	20075238	07613203014854	23-Jul-20	23-Jul-23
20019E7D	20075616	07613203014854	21-Jul-20	21-Jul-23
20019E7D	20075617	07613203014854	21-Jul-20	21-Jul-23



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
20019E7D	20085525	07613203014854	06-Aug-20	06-Aug-23
20019E7D	20085526	07613203014854	07-Aug-20	07-Aug-23
20019E7D	20085774	07613203014854	07-Aug-20	07-Aug-23
20019E7D	20085775	07613203014854	10-Aug-20	10-Aug-23
20019E7D	20086097	07613203014854	10-Aug-20	10-Aug-23
20019E7D	20086098	07613203014854	10-Aug-20	10-Aug-23
20019E7D	20086099	07613203014854	10-Aug-20	10-Aug-23
20019E7D	20086100	07613203014854	11-Aug-20	11-Aug-23
20019E7D	20086101	07613203014854	11-Aug-20	11-Aug-23
20019E7D	20086102	07613203014854	12-Aug-20	12-Aug-23
20019E7D	20086427	07613203014854	17-Aug-20	17-Aug-23
20019E7D	20086428	07613203014854	18-Aug-20	18-Aug-23
20019E7D	20086429	07613203014854	21-Aug-20	21-Aug-23
20019E7D	20095126	07613203014854	09-Sep-20	09-Sep-23
20019E7D	20095127	07613203014854	09-Sep-20	09-Sep-23
20019E7D	20095128	07613203014854	10-Sep-20	10-Sep-23
20019E7D	20095129	07613203014854	10-Sep-20	10-Sep-23
20019E7D	20095130	07613203014854	10-Sep-20	10-Sep-23
20019E7D	20095131	07613203014854	10-Sep-20	10-Sep-23
20019E7D	20095146	07613203014854	10-Sep-20	10-Sep-23
20019E7D	20095147	07613203014854	11-Sep-20	11-Sep-23
20019E7D	20095148	07613203014854	11-Sep-20	11-Sep-23
20019E7D	20095150	07613203014854	21-Sep-20	21-Sep-23
20019E7D	20095151	07613203014854	22-Sep-20	22-Sep-23
20019E7D	20096880	07613203014854	22-Sep-20	22-Sep-23
20019E7D	20096881	07613203014854	22-Sep-20	22-Sep-23
20019E7D	20097041	07613203014854	26-Sep-20	26-Sep-23
20019E7D	20097042	07613203014854	10-Sep-20	10-Sep-23
20019E7D	20097043	07613203014854	28-Sep-20	28-Sep-23
20019E7D	20097044	07613203014854	28-Sep-20	28-Sep-23
20019E7D	20097045	07613203014854	29-Sep-20	29-Sep-23
20019E7D	20097046	07613203014854	29-Sep-20	29-Sep-23
20019E7D	20097047	07613203014854	29-Sep-20	29-Sep-23
20019E7D	20105605	07613203014854	07-Oct-20	07-Oct-23
20019E7D	20105606	07613203014854	07-Oct-20	07-Oct-23
20019E7D	20105607	07613203014854	08-Oct-20	08-Oct-23
20019E7D	20106571	07613203014854	20-Oct-20	20-Oct-20
20019E7D	20106572	07613203014854	20-Oct-20	20-Oct-23
20019E7D	20106573	07613203014854	21-Oct-20	21-Oct-23
20019E7D	20106677	07613203014854	21-Oct-20	21-Oct-23
20019E7D	20106678	07613203014854	22-Oct-20	22-Oct-23
20019E7D	20115001	07613203014854	03-Nov-20	03-Nov-23
20019E7D	20115002	07613203014854	04-Nov-20	04-Nov-23
20019E7D	20115003	07613203014854	05-Nov-20	05-Nov-23
20019E7D	20115004	07613203014854	19-Nov-20	19-Nov-23



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
20019E7D	20115005	07613203014854	20-Nov-20	20-Nov-23
20019E7D	20115006	07613203014854	21-Nov-20	21-Nov-23
20019E7D	20115007	07613203014854	22-Nov-20	22-Nov-23
20019E7D	20115008	07613203014854	22-Nov-20	22-Nov-23
20019E7D	20115009	07613203014854	23-Nov-20	23-Nov-23
20019E7D	20115010	07613203014854	24-Nov-20	24-Nov-23
20019E7D	20115011	07613203014854	24-Nov-20	24-Nov-23
20019E7D	20115012	07613203014854	25-Nov-20	25-Nov-23
20019E7D	20116408	07613203014854	25-Nov-20	25-Nov-23
20019E7D	20116409	07613203014854	25-Nov-20	25-Nov-23
20019E7D	20116410	07613203014854	26-Nov-20	26-Nov-23
20019E7D	20116411	07613203014854	27-Nov-20	27-Nov-23
20019E7D	20116815	07613203014854	27-Nov-20	27-Nov-23
20019E7D	20116830	07613203014854	28-Nov-20	28-Nov-23
20019E7D	20116831	07613203014854	28-Nov-20	28-Nov-23
20019E7D	20116832	07613203014854	29-Nov-20	29-Nov-23
20019E7D	20116961	07613203014854	27-Nov-20	27-Nov-23
20019E7D	20117063	07613203014854	30-Nov-20	30-Nov-23
20019E7D	20117064	07613203014854	30-Nov-20	30-Nov-23
20019E7D	20125036	07613203014854	01-Dec-20	01-Dec-23
20019E7D	20125037	07613203014854	01-Dec-20	01-Dec-23
20019E7D	20125245	07613203014854	01-Dec-20	01-Dec-23
20019E7D	20125246	07613203014854	01-Dec-20	01-Dec-23
20019E7D	20125247	07613203014854	02-Dec-20	02-Dec-23
20019E7D	20125248	07613203014854	02-Dec-20	02-Dec-23
20019E7D	20125249	07613203014854	02-Dec-20	02-Dec-23
20019E7D	20125250	07613203014854	03-Dec-20	03-Dec-23
20019E7D	20125454	07613203014854	03-Dec-20	03-Dec-23
20019E7D	20125455	07613203014854	03-Dec-20	03-Dec-23
20019E7D	20125456	07613203014854	03-Dec-20	03-Dec-23
20019E7D	20125457	07613203014854	04-Dec-20	04-Dec-23
20019E7D	20125518	07613203014854	04-Dec-20	04-Dec-23
20019E7D	20125519	07613203014854	04-Dec-20	04-Dec-23
20019E7D	20125520	07613203014854	04-Dec-20	04-Dec-23
20019E7D	20125589	07613203014854	05-Dec-20	05-Dec-23
20019E7D	20125590	07613203014854	05-Dec-20	05-Dec-23
20019E7D	20125591	07613203014854	05-Dec-20	05-Dec-23
20019E7D	20125628	07613203014854	05-Dec-20	05-Dec-23
20019E7D	20125629	07613203014854	06-Dec-20	06-Dec-23
20019E7D	20125630	07613203014854	07-Dec-20	07-Dec-23
20019E7D	20125631	07613203014854	07-Dec-20	07-Dec-23
20019E7D	20125632	07613203014854	08-Dec-20	08-Dec-23
20019E7D	20125633	07613203014854	16-Dec-20	16-Dec-23
20019E7D	20125634	07613203014854	16-Dec-20	16-Dec-23
20019E7D	20125635	07613203014854	16-Dec-20	16-Dec-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20019E7D	20125637	07613203014854	17-Dec-20	17-Dec-23
20019E7D	20125638	07613203014854	17-Dec-20	17-Dec-23
20019E7D	21015733	07613203014854	14-Jan-21	14-Jan-24
20019E7D	21015735	07613203014854	14-Jan-21	14-Jan-24
20019E7D	21015737	07613203014854	16-Jan-21	16-Jan-24
20019E7D	21015877	07613203014854	15-Jan-21	15-Jan-24
20019E7D	21015878	07613203014854	20-Jan-21	20-Jan-24
20019E7D	21015879	07613203014854	20-Jan-21	20-Jan-24
20019E7D	21015989	07613203014854	20-Jan-21	20-Jan-24
20019E7D	21015990	07613203014854	21-Jan-21	21-Jan-24
20019E7D	21015991	07613203014854	21-Jan-21	21-Jan-24
20019E7D	21016007	07613203014854	21-Jan-21	21-Jan-24
20019E7D	21016025	07613203014854	21-Jan-21	21-Jan-24
20019E7D	21016026	07613203014854	21-Jan-21	21-Jan-24
20019E7D	21016056	07613203014854	22-Jan-21	22-Jan-24
20019E7D	21016057	07613203014854	22-Jan-21	22-Jan-24
20019E7D	21016391	07613203014854	23-Jan-21	23-Jan-24
20019E7D	21016392	07613203014854	23-Jan-21	23-Jan-24
20019E7D	21035751	07613203014854	19-Mar-21	19-Mar-24
20019E7D	21035752	07613203014854	19-Mar-21	19-Mar-24
20019E7D	21035753	07613203014854	19-Mar-21	19-Mar-24
20019E7D	21035754	07613203014854	21-Mar-21	21-Mar-24
20019E7D	21035755	07613203014854	21-Mar-21	21-Mar-24
20019E7D	21035756	07613203014854	21-Mar-21	21-Mar-24
20019E7D	21035757	07613203014854	20-Mar-21	20-Mar-24
20019E7D	21035758	07613203014854	20-Mar-21	20-Mar-24
20019E7D	21035759	07613203014854	23-Mar-21	23-Mar-24
20019E7D	21035760	07613203014854	23-Mar-21	23-Mar-24
20019E7D	21045891	07613203014854	17-Apr-21	17-Apr-24
20019E7D	21046139	07613203014854	28-Apr-21	28-Apr-24
20019E7D	21046140	07613203014854	28-Apr-21	28-Apr-24
20019E7D	21046574	07613203014854	29-Apr-21	29-Apr-24
20019E7D	21046581	07613203014854	29-Apr-21	29-Apr-24
20019E7D	21046656	07613203014854	29-Apr-21	29-Apr-24
20019E7D	21046672	07613203014854	30-Apr-21	30-Apr-24
20019E7D	21046673	07613203014854	29-Apr-21	29-Apr-24
20019E7D	21046717	07613203014854	30-Apr-21	30-Apr-24
20021E	20025613	07613203012416	10-Feb-20	10-Feb-23
20021E	20035447	07613203012416	06-Mar-20	06-Mar-23
20021E	20045445	07613203012416	11-Apr-20	11-Apr-23
20021E	20056567	07613203012416	26-May-20	26-May-23
20021E	20075230	07613203012416	02-Jul-20	02-Jul-23
20021E	20075231	07613203012416	02-Jul-20	02-Jul-23
20021E	20096525	07613203012416	15-Sep-20	15-Sep-23
20021E	20096526	07613203012416	15-Sep-20	15-Sep-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20021E	20115443	07613203012416	14-Nov-20	14-Nov-23
20021E	20125734	07613203012416	15-Dec-20	15-Dec-23
20021E	20125735	07613203012416	16-Dec-20	16-Dec-23
20021E	20126420	07613203012416	17-Dec-20	17-Dec-23
20021E	21015632	07613203012416	12-Jan-21	12-Jan-24
20021E	21015745	07613203012416	12-Jan-21	12-Jan-24
20021E	21015746	07613203012416	12-Jan-21	12-Jan-24
20021E7D	20026035	07613203014236	12-Feb-20	12-Feb-23
20021E7D	20026036	07613203014236	13-Feb-20	13-Feb-23
20021E7D	20056669	07613203014236	26-May-20	26-May-23
20021E7D	20056670	07613203014236	27-May-20	27-May-23
20021E7D	20056671	07613203014236	27-May-20	27-May-23
20021E7D	20066232	07613203014236	16-Jun-20	16-Jun-23
20021E7D	20066233	07613203014236	16-Jun-20	16-Jun-23
20021E7D	20066234	07613203014236	17-Jun-20	17-Jun-23
20021E7D	20066235	07613203014236	17-Jun-20	17-Jun-23
20021E7D	20115387	07613203014236	17-Nov-20	17-Nov-23
20021E7D	20115388	07613203014236	17-Nov-20	17-Nov-23
20021E7D	20115389	07613203014236	18-Nov-20	18-Nov-23
20021E7D	20115390	07613203014236	18-Nov-20	18-Nov-23
20021E7D	20115391	07613203014236	18-Nov-20	18-Nov-23
20021E7D	20125048	07613203014236	01-Dec-20	01-Dec-23
20021E7D	20125313	07613203014236	02-Dec-20	02-Dec-23
20021E7D	20125334	07613203014236	01-Dec-20	01-Dec-23
20021E7D	21015747	07613203014236	13-Jan-21	13-Jan-24
20021E7D	21015748	07613203014236	13-Jan-21	13-Jan-24
20021E7D	21015933	07613203014236	13-Jan-21	13-Jan-24
20021E7D	21015934	07613203014236	13-Jan-21	13-Jan-24
20021E7D	21045417	07613203014236	09-Apr-21	09-Apr-24
20035E	20016826	07613203011990	21-Jan-20	21-Jan-23
20035E	20016827	07613203011990	21-Jan-20	21-Jan-23
20035E	20036929	07613203011990	27-Mar-20	27-Mar-23
20035E	20056944	07613203011990	26-May-20	26-May-23
20035E	20076261	07613203011990	15-Jul-20	15-Jul-23
20035E	20095884	07613203011990	11-Sep-20	11-Sep-23
20035E	20115656	07613203011990	04-Nov-20	04-Nov-23
20035E	20117016	07613203011990	30-Nov-20	30-Nov-23
20035E	21029597	07613203011990	08-Feb-21	08-Feb-24
20037E	20055740	07613203012379	18-May-20	18-May-23
20037E	21016395	07613203012379	25-Jan-21	25-Jan-24
20037E7D	20026037	07613203014229	12-Feb-20	12-Feb-23
20037E7D	20026038	07613203014229	13-Feb-20	13-Feb-23
20037E7D	20026039	07613203014229	13-Feb-20	13-Feb-23
20037E7D	20045035	07613203014229	02-Apr-20	02-Apr-23
20037E7D	20097181	07613203014229	30-Sep-20	30-Sep-23





Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20037E7D	20097182	07613203014229	30-Sep-20	30-Sep-23
20037E7D	20115386	07613203014229	17-Nov-20	17-Nov-23
20038E	20016607	07613203011976	17-Jan-20	17-Jan-23
20038E	20025053	07613203011976	01-Feb-20	01-Feb-23
20038E	20025054	07613203011976	04-Feb-20	04-Feb-23
20038E	20025055	07613203011976	05-Feb-20	05-Feb-23
20038E	20035897	07613203011976	12-Mar-20	12-Mar-23
20038E	20045997	07613203011976	11-Apr-20	11-Apr-23
20038E	20045998	07613203011976	13-Apr-20	13-Apr-23
20038E	20045999	07613203011976	14-Apr-20	14-Apr-23
20038E	20046001	07613203011976	14-Apr-20	14-Apr-23
20038E	20056936	07613203011976	28-May-20	28-May-23
20038E	20066144	07613203011976	27-Jun-20	27-Jun-23
20038E	20066145	07613203011976	29-Jun-20	29-Jun-23
20038E	20075640	07613203011976	08-Jul-20	08-Jul-23
20038E	20075946	07613203011976	08-Jul-20	08-Jul-23
20038E	20075947	07613203011976	09-Jul-20	09-Jul-23
20038E	20086316	07613203011976	10-Aug-20	10-Aug-23
20038E	20086446	07613203011976	25-Aug-20	25-Aug-23
20038E	20095356	07613203011976	18-Sep-20	18-Sep-23
20038E	20105634	07613203011976	02-Oct-20	02-Oct-23
20038E	20105635	07613203011976	03-Oct-20	03-Oct-23
20038E	20105636	07613203011976	03-Oct-20	03-Oct-23
20038E	20115305	07613203011976	05-Nov-20	02-Nov-23
20038E	20115437	07613203011976	05-Nov-20	02-Nov-23
20038E	20117017	07613203011976	30-Nov-20	30-Nov-23
20038E	20117018	07613203011976	30-Nov-20	30-Nov-23
20038E	20117019	07613203011976	30-Nov-20	30-Nov-23
20038E	21015446	07613203011976	09-Jan-21	09-Jan-24
20038E	21015447	07613203011976	09-Jan-21	09-Jan-24
20038E	21015534	07613203011976	09-Jan-21	09-Jan-24
20038E	21015535	07613203011976	09-Jan-21	09-Jan-24
20038E	21015536	07613203011976	10-Jan-21	10-Jan-24
20038E	21016268	07613203011976	18-Jan-21	18-Jan-24
20038E	21016269	07613203011976	18-Jan-21	18-Jan-24
20038E	21016416	07613203011976	25-Jan-21	25-Jan-24
20038E	21016417	07613203011976	25-Jan-21	25-Jan-24
20038E	21016418	07613203011976	25-Jan-21	25-Jan-24
20038E	21016419	07613203011976	25-Jan-21	25-Jan-24
20038E	21016420	07613203011976	25-Jan-21	25-Jan-24
20038E	21026521	07613203011976	22-Feb-21	22-Feb-24
20038E	21026522	07613203011976	22-Feb-21	22-Feb-24
20038E	21026523	07613203011976	22-Feb-21	22-Feb-24
20038E7D	20017143	07613203014083	29-Jan-20	29-Jan-23
20038E7D	20017144	07613203014083	31-Jan-20	31-Jan-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20038E7D	20035163	07613203014083	05-Mar-20	05-Mar-23
20038E7D	20035164	07613203014083	05-Mar-20	05-Mar-23
20038E7D	20035255	07613203014083	06-Mar-20	06-Mar-23
20038E7D	20035256	07613203014083	06-Mar-20	06-Mar-23
20038E7D	20035452	07613203014083	06-Mar-20	06-Mar-23
20038E7D	20035453	07613203014083	07-Mar-20	07-Mar-23
20038E7D	20035454	07613203014083	09-Mar-20	09-Mar-23
20038E7D	20077208	07613203014083	24-Jul-20	24-Jul-23
20038E7D	20077209	07613203014083	24-Jul-20	24-Jul-23
20038E7D	20077210	07613203014083	24-Jul-20	24-Jul-23
20038E7D	20077211	07613203014083	24-Jul-20	24-Jul-23
20038E7D	20077248	07613203014083	27-Jul-20	27-Jul-23
20038E7D	20085800	07613203014083	10-Aug-20	10-Aug-23
20038E7D	20085801	07613203014083	12-Aug-20	12-Aug-23
20038E7D	20085802	07613203014083	13-Aug-20	13-Aug-23
20038E7D	20085803	07613203014083	13-Aug-20	13-Aug-23
20038E7D	20085804	07613203014083	13-Aug-20	13-Aug-23
20038E7D	20085805	07613203014083	14-Aug-20	14-Aug-23
20038E7D	20085806	07613203014083	14-Aug-20	14-Aug-23
20038E7D	20086310	07613203014083	15-Aug-20	15-Aug-23
20038E7D	20086311	07613203014083	15-Aug-20	15-Aug-23
20038E7D	20086312	07613203014083	17-Aug-20	17-Aug-23
20038E7D	20086313	07613203014083	18-Aug-20	18-Aug-23
20038E7D	20086314	07613203014083	18-Aug-20	18-Aug-23
20038E7D	20086315	07613203014083	25-Aug-20	25-Aug-23
20038E7D	20095355	07613203014083	17-Sep-20	17-Sep-23
20038E7D	20106470	07613203014083	16-Oct-20	16-Oct-23
20038E7D	20106471	07613203014083	17-Oct-20	17-Oct-23
20038E7D	20106472	07613203014083	17-Oct-20	17-Oct-23
20038E7D	20106473	07613203014083	29-Oct-20	29-Oct-23
20038E7D	20106474	07613203014083	29-Oct-20	29-Oct-23
20038E7D	20115306	07613203014083	05-Nov-20	05-Nov-23
20038E7D	20115307	07613203014083	06-Nov-20	06-Nov-23
20038E7D	20115308	07613203014083	06-Nov-20	06-Nov-23
20038E7D	20115309	07613203014083	06-Nov-20	06-Nov-23
20038E7D	20115310	07613203014083	07-Nov-20	07-Nov-23
20038E7D	20116459	07613203014083	20-Nov-20	20-Nov-23
20038E7D	20116460	07613203014083	20-Nov-20	20-Nov-23
20038E7D	20116506	07613203014083	21-Nov-20	21-Nov-23
20038E7D	20116843	07613203014083	26-Nov-20	26-Nov-23
20038E7D	20116959	07613203014083	27-Nov-20	27-Nov-23
20038E7D	20116960	07613203014083	27-Nov-20	27-Nov-23
20038E7D	20117003	07613203014083	27-Nov-20	27-Nov-23
20038E7D	20117004	07613203014083	28-Nov-20	28-Nov-23
20038E7D	20126275	07613203014083	18-Dec-20	18-Dec-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20038E7D	20126276	07613203014083	19-Dec-20	19-Dec-23
20038E7D	21026315	07613203014083	17-Feb-21	17-Feb-24
20038E7D	21026316	07613203014083	17-Feb-21	17-Feb-24
20038E7D	21036130	07613203014083	19-Mar-21	19-Mar-24
20038E7D	21036131	07613203014083	20-Mar-21	20-Mar-24
20038E7D	21036132	07613203014083	22-Mar-21	22-Mar-24
20039E	20015291	07613203011952	02-Jan-20	02-Jan-23
20039E	20015292	07613203011952	02-Jan-20	02-Jan-23
20039E	20015293	07613203011952	03-Jan-20	03-Jan-23
20039E	20015294	07613203011952	03-Jan-20	03-Jan-23
20039E	20015295	07613203011952	03-Jan-20	03-Jan-23
20039E	20015296	07613203011952	03-Jan-20	03-Jan-23
20039E	20015297	07613203011952	03-Jan-20	03-Jan-23
20039E	20015298	07613203011952	04-Jan-20	04-Jan-23
20039E	20015299	07613203011952	04-Jan-20	04-Jan-23
20039E	20015300	07613203011952	04-Jan-20	04-Jan-23
20039E	20015301	07613203011952	04-Jan-20	04-Jan-23
20039E	20015302	07613203011952	04-Jan-20	04-Jan-23
20039E	20015303	07613203011952	04-Jan-20	04-Jan-23
20039E	20015389	07613203011952	02-Jan-20	02-Jan-23
20039E	20015487	07613203011952	05-Jan-20	05-Jan-23
20039E	20015488	07613203011952	05-Jan-20	05-Jan-23
20039E	20015489	07613203011952	05-Jan-20	05-Jan-23
20039E	20015490	07613203011952	05-Jan-20	05-Jan-23
20039E	20015491	07613203011952	05-Jan-20	05-Jan-23
20039E	20015513	07613203011952	05-Jan-20	05-Jan-23
20039E	20015514	07613203011952	06-Jan-20	06-Jan-23
20039E	20015515	07613203011952	06-Jan-20	06-Jan-23
20039E	20015516	07613203011952	06-Jan-20	06-Jan-23
20039E	20015517	07613203011952	06-Jan-20	06-Jan-23
20039E	20015518	07613203011952	06-Jan-20	06-Jan-23
20039E	20015519	07613203011952	06-Jan-20	06-Jan-23
20039E	20015520	07613203011952	06-Jan-20	06-Jan-23
20039E	20015521	07613203011952	06-Jan-20	06-Jan-23
20039E	20015887	07613203011952	10-Jan-20	10-Jan-23
20039E	20015892	07613203011952	11-Jan-20	11-Jan-23
20039E	20015894	07613203011952	13-Jan-20	13-Jan-23
20039E	20016032	07613203011952	13-Jan-20	13-Jan-23
20039E	20016033	07613203011952	13-Jan-20	13-Jan-23
20039E	20016034	07613203011952	13-Jan-20	13-Jan-23
20039E	20016035	07613203011952	14-Jan-20	14-Jan-23
20039E	20016036	07613203011952	14-Jan-20	14-Jan-23
20039E	20016037	07613203011952	15-Jan-20	15-Jan-23
20039E	20016101	07613203011952	14-Jan-20	14-Jan-23
20039E	20016102	07613203011952	15-Jan-20	15-Jan-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E	20016103	07613203011952	15-Jan-20	15-Jan-23
20039E	20017312	07613203011952	29-Jan-20	29-Jan-23
20039E	20017313	07613203011952	29-Jan-20	29-Jan-23
20039E	20017314	07613203011952	29-Jan-20	29-Jan-23
20039E	20017315	07613203011952	30-Jan-20	30-Jan-23
20039E	20017316	07613203011952	30-Jan-20	30-Jan-23
20039E	20025773	07613203011952	24-Feb-20	24-Feb-23
20039E	20025820	07613203011952	19-Feb-20	19-Feb-23
20039E	20025821	07613203011952	11-Feb-20	11-Feb-23
20039E	20025822	07613203011952	11-Feb-20	11-Feb-23
20039E	20025823	07613203011952	11-Feb-20	11-Feb-23
20039E	20025824	07613203011952	11-Feb-20	11-Feb-23
20039E	20025825	07613203011952	11-Feb-20	11-Feb-23
20039E	20025826	07613203011952	12-Feb-20	12-Feb-23
20039E	20025827	07613203011952	12-Feb-20	12-Feb-23
20039E	20025991	07613203011952	12-Feb-20	12-Feb-23
20039E	20025992	07613203011952	24-Feb-20	24-Feb-23
20039E	20025993	07613203011952	25-Feb-20	25-Feb-23
20039E	20025994	07613203011952	25-Feb-20	25-Feb-23
20039E	20025995	07613203011952	26-Feb-20	26-Feb-23
20039E	20025996	07613203011952	26-Feb-20	26-Feb-23
20039E	20025997	07613203011952	26-Feb-20	26-Feb-23
20039E	20025998	07613203011952	26-Feb-20	26-Feb-23
20039E	20025999	07613203011952	27-Feb-20	27-Feb-23
20039E	20026001	07613203011952	27-Feb-20	27-Feb-23
20039E	20035865	07613203011952	10-Mar-20	10-Mar-23
20039E	20035866	07613203011952	09-Mar-20	09-Mar-23
20039E	20035867	07613203011952	10-Mar-20	10-Mar-23
20039E	20035868	07613203011952	10-Mar-20	10-Mar-23
20039E	20035869	07613203011952	11-Mar-20	11-Mar-23
20039E	20035870	07613203011952	11-Mar-20	11-Mar-23
20039E	20035871	07613203011952	11-Mar-20	11-Mar-23
20039E	20035872	07613203011952	11-Mar-20	11-Mar-23
20039E	20035973	07613203011952	11-Mar-20	11-Mar-23
20039E	20045190	07613203011952	02-Apr-20	02-Apr-23
20039E	20045191	07613203011952	03-Apr-20	03-Apr-23
20039E	20045192	07613203011952	03-Apr-20	03-Apr-23
20039E	20045193	07613203011952	03-Apr-20	03-Apr-23
20039E	20045308	07613203011952	06-Apr-20	06-Apr-23
20039E	20045309	07613203011952	06-Apr-20	06-Apr-23
20039E	20045310	07613203011952	06-Apr-20	06-Apr-23
20039E	20045311	07613203011952	06-Apr-20	06-Apr-23
20039E	20045439	07613203011952	07-Apr-20	07-Apr-23
20039E	20045440	07613203011952	07-Apr-20	07-Apr-23
20039E	20045441	07613203011952	07-Apr-20	07-Apr-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E	20045442	07613203011952	07-Apr-20	07-Apr-23
20039E	20045443	07613203011952	07-Apr-20	07-Apr-23
20039E	20045444	07613203011952	08-Apr-20	08-Apr-23
20039E	20045702	07613203011952	08-Apr-20	08-Apr-23
20039E	20045703	07613203011952	08-Apr-20	08-Apr-23
20039E	20045704	07613203011952	08-Apr-20	08-Apr-23
20039E	20045705	07613203011952	09-Apr-20	09-Apr-23
20039E	20045706	07613203011952	09-Apr-20	09-Apr-23
20039E	20045707	07613203011952	09-Apr-20	09-Apr-23
20039E	20045708	07613203011952	09-Apr-20	09-Apr-23
20039E	20045709	07613203011952	17-Apr-20	17-Apr-23
20039E	20045710	07613203011952	17-Apr-20	17-Apr-23
20039E	20045711	07613203011952	18-Apr-20	18-Apr-23
20039E	20046120	07613203011952	20-Apr-20	20-Apr-23
20039E	20046121	07613203011952	20-Apr-20	20-Apr-23
20039E	20046240	07613203011952	20-Apr-20	20-Apr-23
20039E	20046241	07613203011952	20-Apr-20	20-Apr-23
20039E	20046242	07613203011952	20-Apr-20	20-Apr-23
20039E	20046243	07613203011952	20-Apr-20	20-Apr-23
20039E	20046382	07613203011952	21-Apr-20	21-Apr-23
20039E	20046383	07613203011952	21-Apr-20	21-Apr-23
20039E	20046384	07613203011952	21-Apr-20	21-Apr-23
20039E	20046456	07613203011952	21-Apr-20	21-Apr-23
20039E	20046457	07613203011952	21-Apr-20	21-Apr-23
20039E	20046458	07613203011952	21-Apr-20	21-Apr-23
20039E	20046459	07613203011952	21-Apr-20	21-Apr-23
20039E	20046460	07613203011952	22-Apr-20	22-Apr-23
20039E	20046662	07613203011952	22-Apr-20	22-Apr-23
20039E	20046663	07613203011952	22-Apr-20	22-Apr-23
20039E	20046672	07613203011952	22-Apr-20	22-Apr-23
20039E	20046682	07613203011952	22-Apr-20	22-Apr-23
20039E	20046683	07613203011952	24-Apr-20	24-Apr-23
20039E	20046684	07613203011952	25-Apr-20	25-Apr-23
20039E	20046685	07613203011952	25-Apr-20	25-Apr-23
20039E	20046779	07613203011952	25-Apr-20	25-Apr-23
20039E	20046780	07613203011952	25-Apr-20	25-Apr-23
20039E	20046781	07613203011952	25-Apr-20	25-Apr-23
20039E	20046875	07613203011952	27-Apr-20	27-Apr-23
20039E	20046876	07613203011952	27-Apr-20	27-Apr-23
20039E	20055136	07613203011952	01-May-20	01-May-23
20039E	20055137	07613203011952	01-May-20	01-May-23
20039E	20055138	07613203011952	01-May-20	01-May-23
20039E	20055139	07613203011952	01-May-20	01-May-23
20039E	20055140	07613203011952	02-May-20	02-May-23
20039E	20055141	07613203011952	04-May-20	04-May-23



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
20039E	20055142	07613203011952	16-May-20	16-May-23
20039E	20055143	07613203011952	04-May-20	04-May-23
20039E	20055144	07613203011952	16-May-20	16-May-23
20039E	20055145	07613203011952	17-May-20	17-May-23
20039E	20055146	07613203011952	17-May-20	17-May-23
20039E	20065646	07613203011952	09-Jun-20	09-Jun-23
20039E	20065647	07613203011952	10-Jun-20	10-Jun-23
20039E	20065648	07613203011952	10-Jun-20	10-Jun-23
20039E	20065851	07613203011952	10-Jun-20	10-Jun-23
20039E	20065852	07613203011952	10-Jun-20	10-Jun-23
20039E	20065853	07613203011952	10-Jun-20	10-Jun-23
20039E	20065854	07613203011952	10-Jun-20	10-Jun-23
20039E	20065855	07613203011952	10-Jun-20	10-Jun-23
20039E	20065856	07613203011952	10-Jun-20	10-Jun-23
20039E	20065857	07613203011952	11-Jun-20	11-Jun-23
20039E	20065858	07613203011952	11-Jun-20	11-Jun-23
20039E	20065859	07613203011952	11-Jun-20	11-Jun-23
20039E	20065860	07613203011952	11-Jun-20	11-Jun-23
20039E	20065861	07613203011952	11-Jun-20	11-Jun-23
20039E	20065862	07613203011952	11-Jun-20	11-Jun-23
20039E	20075405	07613203011952	04-Jul-20	04-Jul-23
20039E	20075406	07613203011952	04-Jul-20	04-Jul-23
20039E	20075407	07613203011952	04-Jul-20	04-Jul-23
20039E	20075408	07613203011952	04-Jul-20	04-Jul-23
20039E	20075472	07613203011952	06-Jul-20	06-Jul-23
20039E	20075473	07613203011952	06-Jul-20	06-Jul-23
20039E	20075474	07613203011952	07-Jul-20	07-Jul-23
20039E	20075475	07613203011952	07-Jul-20	07-Jul-23
20039E	20075476	07613203011952	07-Jul-20	07-Jul-23
20039E	20075477	07613203011952	07-Jul-20	07-Jul-23
20039E	20075478	07613203011952	07-Jul-20	07-Jul-23
20039E	20075479	07613203011952	07-Jul-20	07-Jul-23
20039E	20075480	07613203011952	07-Jul-20	07-Jul-23
20039E	20075784	07613203011952	07-Jul-20	07-Jul-23
20039E	20075785	07613203011952	08-Jul-20	08-Jul-23
20039E	20075786	07613203011952	08-Jul-20	08-Jul-23
20039E	20075787	07613203011952	08-Jul-20	08-Jul-23
20039E	20075788	07613203011952	08-Jul-20	08-Jul-23
20039E	20075789	07613203011952	08-Jul-20	08-Jul-23
20039E	20075850	07613203011952	08-Jul-20	08-Jul-23
20039E	20075851	07613203011952	14-Jul-20	14-Jul-23
20039E	20075852	07613203011952	14-Jul-20	14-Jul-23
20039E	20075919	07613203011952	14-Jul-20	14-Jul-23
20039E	20075920	07613203011952	14-Jul-20	14-Jul-23
20039E	20075921	07613203011952	16-Jul-20	16-Jul-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E	20075922	07613203011952	16-Jul-20	16-Jul-23
20039E	20075923	07613203011952	15-Jul-20	15-Jul-23
20039E	20076062	07613203011952	16-Jul-20	16-Jul-23
20039E	20076063	07613203011952	16-Jul-20	16-Jul-23
20039E	20076064	07613203011952	16-Jul-20	16-Jul-23
20039E	20076065	07613203011952	16-Jul-20	16-Jul-23
20039E	20076066	07613203011952	15-Jul-20	15-Jul-23
20039E	20076067	07613203011952	16-Jul-20	16-Jul-23
20039E	20076068	07613203011952	16-Jul-20	16-Jul-23
20039E	20076069	07613203011952	16-Jul-20	16-Jul-23
20039E	20076242	07613203011952	17-Jul-20	17-Jul-23
20039E	20076243	07613203011952	17-Jul-20	17-Jul-23
20039E	20076244	07613203011952	17-Jul-20	17-Jul-23
20039E	20076245	07613203011952	17-Jul-20	17-Jul-23
20039E	20076572	07613203011952	18-Jul-20	18-Jul-23
20039E	20076573	07613203011952	18-Jul-20	18-Jul-23
20039E	20076796	07613203011952	17-Jul-20	17-Jul-23
20039E	20076797	07613203011952	18-Jul-20	18-Jul-23
20039E	20076798	07613203011952	19-Jul-20	19-Jul-23
20039E	20076799	07613203011952	20-Jul-20	20-Jul-23
20039E	20076800	07613203011952	20-Jul-20	20-Jul-23
20039E	20076801	07613203011952	20-Jul-20	20-Jul-23
20039E	20076802	07613203011952	20-Jul-20	20-Jul-23
20039E	20076803	07613203011952	20-Jul-20	20-Jul-23
20039E	20076903	07613203011952	20-Jul-20	20-Jul-23
20039E	20095296	07613203011952	11-Sep-20	11-Sep-23
20039E	20095297	07613203011952	11-Sep-20	11-Sep-23
20039E	20095298	07613203011952	12-Sep-20	12-Sep-23
20039E	20095299	07613203011952	12-Sep-20	12-Sep-23
20039E	20095996	07613203011952	11-Sep-20	11-Sep-23
20039E	20095997	07613203011952	12-Sep-20	12-Sep-23
20039E	20095998	07613203011952	12-Sep-20	12-Sep-23
20039E	20095999	07613203011952	12-Sep-20	12-Sep-23
20039E	20096000	07613203011952	12-Sep-20	12-Sep-23
20039E	20096001	07613203011952	14-Sep-20	14-Sep-23
20039E	20096002	07613203011952	14-Sep-20	14-Sep-23
20039E	20096003	07613203011952	14-Sep-20	14-Sep-23
20039E	20096004	07613203011952	14-Sep-20	14-Sep-23
20039E	20096321	07613203011952	14-Sep-20	14-Sep-23
20039E	20096322	07613203011952	14-Sep-20	14-Sep-23
20039E	20096323	07613203011952	14-Sep-20	14-Sep-23
20039E	20096324	07613203011952	15-Sep-20	15-Sep-23
20039E	20096325	07613203011952	15-Sep-20	15-Sep-23
20039E	20096326	07613203011952	15-Sep-20	15-Sep-23
20039E	20096327	07613203011952	15-Sep-20	15-Sep-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E	20096328	07613203011952	15-Sep-20	15-Sep-23
20039E	20096329	07613203011952	15-Sep-20	15-Sep-23
20039E	20096330	07613203011952	16-Sep-20	16-Sep-23
20039E	20096331	07613203011952	16-Sep-20	16-Sep-23
20039E	20096332	07613203011952	16-Sep-20	16-Sep-23
20039E	20096333	07613203011952	16-Sep-20	16-Sep-23
20039E	20096334	07613203011952	16-Sep-20	16-Sep-23
20039E	20096335	07613203011952	17-Sep-20	17-Sep-23
20039E	20096336	07613203011952	17-Sep-20	17-Sep-23
20039E	20096337	07613203011952	17-Sep-20	17-Sep-23
20039E	20096338	07613203011952	17-Sep-20	17-Sep-23
20039E	20096339	07613203011952	17-Sep-20	17-Sep-23
20039E	20096340	07613203011952	17-Sep-20	17-Sep-23
20039E	20096341	07613203011952	18-Sep-20	18-Sep-23
20039E	20096342	07613203011952	18-Sep-20	18-Sep-23
20039E	20096343	07613203011952	18-Sep-20	18-Sep-23
20039E	20096344	07613203011952	18-Sep-20	18-Sep-23
20039E	20096345	07613203011952	18-Sep-20	18-Sep-23
20039E	20096346	07613203011952	18-Sep-20	18-Sep-23
20039E	20096347	07613203011952	19-Sep-20	19-Sep-23
20039E	20096348	07613203011952	19-Sep-20	19-Sep-23
20039E	20096349	07613203011952	19-Sep-20	19-Sep-23
20039E	20096350	07613203011952	19-Sep-20	19-Sep-23
20039E	20096351	07613203011952	19-Sep-20	19-Sep-23
20039E	20096352	07613203011952	19-Sep-20	19-Sep-23
20039E	20096353	07613203011952	19-Sep-20	19-Sep-23
20039E	20096354	07613203011952	19-Sep-20	19-Sep-23
20039E	20096355	07613203011952	20-Sep-20	20-Sep-23
20039E	20096356	07613203011952	20-Sep-20	20-Sep-23
20039E	20096357	07613203011952	20-Sep-20	20-Sep-23
20039E	20096358	07613203011952	20-Sep-20	20-Sep-23
20039E	20096359	07613203011952	20-Sep-20	20-Sep-23
20039E	20096932	07613203011952	22-Sep-20	22-Sep-23
20039E	20096933	07613203011952	22-Sep-20	22-Sep-23
20039E	20096934	07613203011952	23-Sep-20	23-Sep-23
20039E	20096935	07613203011952	23-Sep-20	23-Sep-23
20039E	20096936	07613203011952	23-Sep-20	23-Sep-23
20039E	20096937	07613203011952	23-Sep-20	23-Sep-23
20039E	20096938	07613203011952	23-Sep-20	23-Sep-23
20039E	20096939	07613203011952	23-Sep-20	23-Sep-23
20039E	20096961	07613203011952	24-Sep-20	24-Sep-23
20039E	20096994	07613203011952	25-Sep-20	25-Sep-23
20039E	20096995	07613203011952	25-Sep-20	25-Sep-23
20039E	20097259	07613203011952	30-Sep-20	30-Sep-23
20039E	20105081	07613203011952	01-Oct-20	01-Oct-23





Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E	20105082	07613203011952	01-Oct-20	01-Oct-23
20039E	20105083	07613203011952	02-Oct-20	02-Oct-23
20039E	20105298	07613203011952	02-Oct-20	02-Oct-23
20039E	20105299	07613203011952	02-Oct-20	02-Oct-23
20039E	20105300	07613203011952	02-Oct-20	02-Oct-23
20039E	20105305	07613203011952	02-Oct-20	02-Oct-23
20039E	20105306	07613203011952	02-Oct-20	02-Oct-23
20039E	20105307	07613203011952	03-Oct-20	03-Oct-23
20039E	20105595	07613203011952	05-Oct-20	05-Oct-23
20039E	20105596	07613203011952	05-Oct-20	05-Oct-23
20039E	20105597	07613203011952	05-Oct-20	05-Oct-23
20039E	20105766	07613203011952	06-Oct-20	06-Oct-23
20039E	20105767	07613203011952	06-Oct-20	06-Oct-23
20039E	20105768	07613203011952	06-Oct-20	06-Oct-23
20039E	20105769	07613203011952	06-Oct-20	06-Oct-23
20039E	20105770	07613203011952	06-Oct-20	06-Oct-23
20039E	20105771	07613203011952	08-Oct-20	08-Oct-23
20039E	20105772	07613203011952	09-Oct-20	09-Oct-23
20039E	20105773	07613203011952	09-Oct-20	09-Oct-23
20039E	20105774	07613203011952	09-Oct-20	09-Oct-23
20039E	20105775	07613203011952	09-Oct-20	09-Oct-23
20039E	20105882	07613203011952	10-Oct-20	10-Oct-23
20039E	20105883	07613203011952	10-Oct-20	10-Oct-23
20039E	20105884	07613203011952	10-Oct-20	10-Oct-23
20039E	20105885	07613203011952	12-Oct-20	12-Oct-23
20039E	20105886	07613203011952	12-Oct-20	12-Oct-23
20039E	20105887	07613203011952	12-Oct-20	12-Oct-23
20039E	20105888	07613203011952	12-Oct-20	12-Oct-23
20039E	20105889	07613203011952	13-Oct-20	13-Oct-23
20039E	20105890	07613203011952	13-Oct-20	13-Oct-23
20039E	20105891	07613203011952	13-Oct-20	13-Oct-23
20039E	20105892	07613203011952	13-Oct-20	13-Oct-23
20039E	20105893	07613203011952	13-Oct-20	13-Oct-23
20039E	20106333	07613203011952	13-Oct-20	13-Oct-23
20039E	20106457	07613203011952	16-Oct-20	16-Oct-23
20039E	20106458	07613203011952	17-Oct-20	17-Oct-23
20039E	20106459	07613203011952	17-Oct-20	17-Oct-23
20039E	20106460	07613203011952	19-Oct-20	19-Oct-23
20039E	20106461	07613203011952	19-Oct-20	19-Oct-23
20039E	20106462	07613203011952	19-Oct-20	19-Oct-23
20039E	20106463	07613203011952	20-Oct-20	20-Oct-23
20039E	20106464	07613203011952	20-Oct-20	20-Oct-23
20039E	20106465	07613203011952	20-Oct-20	20-Oct-23
20039E	20115402	07613203011952	04-Nov-20	04-Nov-23
20039E	20115403	07613203011952	04-Nov-20	04-Nov-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E	20115404	07613203011952	05-Nov-20	05-Nov-23
20039E	20115405	07613203011952	05-Nov-20	05-Nov-23
20039E	20115406	07613203011952	05-Nov-20	05-Nov-23
20039E	20115407	07613203011952	09-Nov-20	09-Nov-23
20039E	20115408	07613203011952	09-Nov-20	09-Nov-23
20039E	20115409	07613203011952	09-Nov-20	09-Nov-23
20039E	20115410	07613203011952	09-Nov-20	09-Nov-23
20039E	20115411	07613203011952	09-Nov-20	09-Nov-23
20039E	20115412	07613203011952	09-Nov-20	09-Nov-23
20039E	20115413	07613203011952	10-Nov-20	10-Nov-23
20039E	20115414	07613203011952	10-Nov-20	10-Nov-23
20039E	20115415	07613203011952	10-Nov-20	10-Nov-23
20039E	20115416	07613203011952	10-Nov-20	10-Nov-23
20039E	20115417	07613203011952	05-Nov-20	05-Nov-23
20039E	20115418	07613203011952	05-Nov-20	05-Nov-23
20039E	20115419	07613203011952	05-Nov-20	05-Nov-23
20039E	20115420	07613203011952	05-Nov-20	05-Nov-23
20039E	20115421	07613203011952	05-Nov-20	05-Nov-23
20039E	20115422	07613203011952	06-Nov-20	06-Nov-23
20039E	20115423	07613203011952	06-Nov-20	06-Nov-23
20039E	20115424	07613203011952	06-Nov-20	06-Nov-23
20039E	20115425	07613203011952	06-Nov-20	06-Nov-23
20039E	20115426	07613203011952	06-Nov-20	06-Nov-23
20039E	20115427	07613203011952	06-Nov-20	06-Nov-23
20039E	20115428	07613203011952	07-Nov-20	07-Nov-23
20039E	20115429	07613203011952	07-Nov-20	07-Nov-23
20039E	20115444	07613203011952	07-Nov-20	07-Nov-23
20039E	20115445	07613203011952	07-Nov-20	07-Nov-23
20039E	20115446	07613203011952	07-Nov-20	07-Nov-23
20039E	20115447	07613203011952	07-Nov-20	07-Nov-23
20039E	20115448	07613203011952	08-Nov-20	08-Nov-23
20039E	20115449	07613203011952	08-Nov-20	08-Nov-23
20039E	20115450	07613203011952	09-Nov-20	09-Nov-23
20039E	20115795	07613203011952	11-Nov-20	11-Nov-23
20039E	20115796	07613203011952	11-Nov-20	11-Nov-23
20039E	20115797	07613203011952	11-Nov-20	11-Nov-23
20039E	20115798	07613203011952	11-Nov-20	11-Nov-23
20039E	20115799	07613203011952	11-Nov-20	11-Nov-23
20039E	20115800	07613203011952	10-Nov-20	10-Nov-23
20039E	20115801	07613203011952	11-Nov-20	11-Nov-23
20039E	20116168	07613203011952	11-Nov-20	11-Nov-23
20039E	20116169	07613203011952	11-Nov-20	11-Nov-23
20039E	20116170	07613203011952	19-Nov-20	19-Nov-23
20039E	20116171	07613203011952	19-Nov-20	19-Nov-23
20039E	20116172	07613203011952	19-Nov-20	19-Nov-23



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
20039E	20116173	07613203011952	20-Nov-20	20-Nov-23
20039E	20116174	07613203011952	20-Nov-20	20-Nov-23
20039E	20116175	07613203011952	20-Nov-20	20-Nov-23
20039E	20116176	07613203011952	21-Nov-20	21-Nov-23
20039E	20116523	07613203011952	20-Nov-20	20-Nov-23
20039E	20116524	07613203011952	21-Nov-20	21-Nov-23
20039E	20125584	07613203011952	04-Dec-20	04-Dec-23
20039E	20125585	07613203011952	04-Dec-20	04-Dec-23
20039E	20125662	07613203011952	04-Dec-20	04-Dec-23
20039E	20125663	07613203011952	05-Dec-20	05-Dec-23
20039E	20125664	07613203011952	05-Dec-20	05-Dec-23
20039E	20125665	07613203011952	05-Dec-20	05-Dec-23
20039E	20125765	07613203011952	07-Dec-20	07-Dec-23
20039E	20125766	07613203011952	07-Dec-20	07-Dec-23
20039E	20125767	07613203011952	07-Dec-20	07-Dec-23
20039E	20125768	07613203011952	07-Dec-20	07-Dec-23
20039E	20125769	07613203011952	07-Dec-20	07-Dec-23
20039E	20125795	07613203011952	07-Dec-20	07-Dec-23
20039E	20125796	07613203011952	08-Dec-20	08-Dec-23
20039E	20125797	07613203011952	08-Dec-20	08-Dec-23
20039E	20125798	07613203011952	08-Dec-20	08-Dec-23
20039E	20125799	07613203011952	08-Dec-20	08-Dec-23
20039E	20125881	07613203011952	08-Dec-20	08-Dec-23
20039E	20125882	07613203011952	08-Dec-20	08-Dec-23
20039E	20125883	07613203011952	09-Dec-20	09-Dec-23
20039E	20125884	07613203011952	09-Dec-20	09-Dec-23
20039E	20125885	07613203011952	09-Dec-20	09-Dec-23
20039E	20125922	07613203011952	09-Dec-20	09-Dec-23
20039E	20125923	07613203011952	10-Dec-20	10-Dec-23
20039E	20125924	07613203011952	10-Dec-20	10-Dec-23
20039E	20125925	07613203011952	10-Dec-20	10-Dec-23
20039E	20125926	07613203011952	10-Dec-20	10-Dec-23
20039E	20125927	07613203011952	10-Dec-20	10-Dec-23
20039E	20125928	07613203011952	10-Dec-20	10-Dec-23
20039E	20125930	07613203011952	10-Dec-20	10-Dec-23
20039E	20125931	07613203011952	11-Dec-20	11-Dec-23
20039E	20125932	07613203011952	11-Dec-20	11-Dec-23
20039E	20125933	07613203011952	11-Dec-20	11-Dec-23
20039E	20125934	07613203011952	15-Dec-20	15-Dec-23
20039E	20126083	07613203011952	15-Dec-20	15-Dec-23
20039E	20126084	07613203011952	18-Dec-20	18-Dec-23
20039E	20126085	07613203011952	18-Dec-20	18-Dec-23
20039E	20126086	07613203011952	18-Dec-20	18-Dec-23
20039E	20126087	07613203011952	18-Dec-20	18-Dec-23
20039E	20126088	07613203011952	18-Dec-20	18-Dec-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E	20126089	07613203011952	18-Dec-20	18-Dec-23
20039E	20126090	07613203011952	18-Dec-20	18-Dec-23
20039E	21015161	07613203011952	06-Jan-21	06-Jan-24
20039E	21015162	07613203011952	06-Jan-21	06-Jan-24
20039E	21015163	07613203011952	07-Jan-21	07-Jan-24
20039E	21015164	07613203011952	07-Jan-21	07-Jan-24
20039E	21015165	07613203011952	05-Jan-21	05-Jan-24
20039E	21015166	07613203011952	05-Jan-21	05-Jan-24
20039E	21015167	07613203011952	06-Jan-21	06-Jan-24
20039E	21015168	07613203011952	06-Jan-21	06-Jan-24
20039E	21015169	07613203011952	06-Jan-21	06-Jan-24
20039E	21015265	07613203011952	07-Jan-21	07-Jan-24
20039E	21015266	07613203011952	07-Jan-21	07-Jan-24
20039E	21015267	07613203011952	08-Jan-21	08-Jan-24
20039E	21015268	07613203011952	08-Jan-21	08-Jan-24
20039E	21015269	07613203011952	08-Jan-21	08-Jan-24
20039E	21015270	07613203011952	08-Jan-21	08-Jan-24
20039E	21015457	07613203011952	09-Jan-21	09-Jan-24
20039E	21015458	07613203011952	09-Jan-21	09-Jan-24
20039E	21015459	07613203011952	09-Jan-21	09-Jan-24
20039E	21015460	07613203011952	11-Jan-21	11-Jan-24
20039E	21015461	07613203011952	11-Jan-21	11-Jan-24
20039E	21015462	07613203011952	11-Jan-21	11-Jan-24
20039E	21015463	07613203011952	11-Jan-21	11-Jan-24
20039E	21015464	07613203011952	12-Jan-21	12-Jan-24
20039E	21025789	07613203011952	09-Feb-21	09-Feb-24
20039E	21025790	07613203011952	09-Feb-21	09-Feb-24
20039E	21025791	07613203011952	09-Feb-21	09-Feb-24
20039E	21025792	07613203011952	09-Feb-21	09-Feb-24
20039E	21025793	07613203011952	09-Feb-21	09-Feb-24
20039E	21025794	07613203011952	09-Feb-21	09-Feb-24
20039E	21025795	07613203011952	09-Feb-21	09-Feb-24
20039E	21025797	07613203011952	09-Feb-21	09-Feb-24
20039E	21025798	07613203011952	09-Feb-21	09-Feb-24
20039E	21025799	07613203011952	10-Feb-21	10-Feb-24
20039E	21025800	07613203011952	10-Feb-21	10-Feb-24
20039E	21025801	07613203011952	10-Feb-21	10-Feb-24
20039E	21025802	07613203011952	10-Feb-21	10-Feb-24
20039E	21025803	07613203011952	10-Feb-21	10-Feb-24
20039E	21025804	07613203011952	10-Feb-21	10-Feb-24
20039E	21025805	07613203011952	18-Feb-21	18-Feb-24
20039E	21025806	07613203011952	18-Feb-21	18-Feb-24
20039E	21025807	07613203011952	18-Feb-21	18-Feb-24
20039E	21025808	07613203011952	13-Feb-21	13-Feb-24
20039E	21025809	07613203011952	18-Feb-21	18-Feb-24



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E	21025810	07613203011952	18-Feb-21	18-Feb-24
20039E	21025881	07613203011952	18-Feb-21	18-Feb-24
20039E	21025882	07613203011952	18-Feb-21	18-Feb-24
20039E	21025886	07613203011952	19-Feb-21	19-Feb-24
20039E	21025887	07613203011952	19-Feb-21	19-Feb-24
20039E	21025892	07613203011952	19-Feb-21	19-Feb-24
20039E	21035761	07613203011952	11-Mar-21	11-Mar-24
20039E	21035762	07613203011952	11-Mar-21	11-Mar-24
20039E	21035763	07613203011952	11-Mar-21	11-Mar-24
20039E	21035764	07613203011952	11-Mar-21	11-Mar-24
20039E	21035765	07613203011952	11-Mar-21	11-Mar-24
20039E	21035766	07613203011952	12-Mar-21	12-Mar-24
20039E	21035767	07613203011952	12-Mar-21	12-Mar-24
20039E	21035924	07613203011952	12-Mar-21	12-Mar-24
20039E	21035930	07613203011952	13-Mar-21	13-Mar-24
20039E	21035931	07613203011952	13-Mar-21	13-Mar-24
20039E	21035932	07613203011952	14-Mar-21	14-Mar-24
20039E	21035933	07613203011952	14-Mar-21	14-Mar-24
20039E	21035934	07613203011952	17-Mar-21	17-Mar-24
20039E	21035935	07613203011952	17-Mar-21	17-Mar-24
20039E	21035936	07613203011952	17-Mar-21	17-Mar-24
20039E	21035937	07613203011952	17-Mar-21	17-Mar-24
20039E	21035938	07613203011952	17-Mar-21	17-Mar-24
20039E	21035939	07613203011952	17-Mar-21	17-Mar-24
20039E	21035940	07613203011952	17-Mar-21	17-Mar-24
20039E	21035941	07613203011952	18-Mar-21	18-Mar-24
20039E	21035942	07613203011952	18-Mar-21	18-Mar-24
20039E	21035943	07613203011952	18-Mar-21	18-Mar-24
20039E	21045429	07613203011952	07-Apr-21	07-Apr-24
20039E	21045430	07613203011952	08-Apr-21	08-Apr-24
20039E	21045431	07613203011952	08-Apr-21	08-Apr-24
20039E	21045432	07613203011952	08-Apr-21	08-Apr-24
20039E	21045433	07613203011952	08-Apr-21	08-Apr-24
20039E	21045434	07613203011952	08-Apr-21	08-Apr-24
20039E	21045435	07613203011952	08-Apr-21	08-Apr-24
20039E	21045436	07613203011952	08-Apr-21	08-Apr-24
20039E7D	20016252	07613203014243	16-Jan-20	16-Jan-23
20039E7D	20016253	07613203014243	16-Jan-20	16-Jan-23
20039E7D	20016254	07613203014243	16-Jan-20	16-Jan-23
20039E7D	20016255	07613203014243	17-Jan-20	17-Jan-23
20039E7D	20016256	07613203014243	17-Jan-20	17-Jan-23
20039E7D	20016257	07613203014243	17-Jan-20	17-Jan-23
20039E7D	20016258	07613203014243	17-Jan-20	17-Jan-23
20039E7D	20016259	07613203014243	25-Jan-20	25-Jan-23
20039E7D	20016260	07613203014243	25-Jan-20	25-Jan-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E7D	20016261	07613203014243	26-Jan-20	26-Jan-23
20039E7D	20016540	07613203014243	26-Jan-20	26-Jan-23
20039E7D	20016541	07613203014243	26-Jan-20	26-Jan-23
20039E7D	20016542	07613203014243	27-Jan-20	27-Jan-23
20039E7D	20016543	07613203014243	27-Jan-20	27-Jan-23
20039E7D	20016544	07613203014243	27-Jan-20	27-Jan-23
20039E7D	20016689	07613203014243	27-Jan-20	27-Jan-23
20039E7D	20016690	07613203014243	27-Jan-20	27-Jan-23
20039E7D	20016691	07613203014243	27-Jan-20	27-Jan-23
20039E7D	20016692	07613203014243	27-Jan-20	27-Jan-23
20039E7D	20016693	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017011	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017012	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017013	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017014	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017015	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017073	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017074	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017075	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20017076	07613203014243	28-Jan-20	28-Jan-23
20039E7D	20027157	07613203014243	28-Feb-20	28-Feb-23
20039E7D	20027158	07613203014243	28-Feb-20	28-Feb-23
20039E7D	20027159	07613203014243	28-Feb-20	28-Feb-23
20039E7D	20027160	07613203014243	28-Feb-20	28-Feb-23
20039E7D	20035250	07613203014243	05-Mar-20	05-Mar-23
20039E7D	20035251	07613203014243	05-Mar-20	05-Mar-23
20039E7D	20035685	07613203014243	12-Mar-20	12-Mar-23
20039E7D	20035686	07613203014243	12-Mar-20	12-Mar-23
20039E7D	20035687	07613203014243	12-Mar-20	12-Mar-23
20039E7D	20036089	07613203014243	13-Mar-20	13-Mar-23
20039E7D	20036090	07613203014243	13-Mar-20	13-Mar-23
20039E7D	20036091	07613203014243	14-Mar-20	14-Mar-23
20039E7D	20036092	07613203014243	17-Mar-20	17-Mar-23
20039E7D	20036184	07613203014243	17-Mar-20	17-Mar-23
20039E7D	20036185	07613203014243	18-Mar-20	18-Mar-23
20039E7D	20036186	07613203014243	18-Mar-20	18-Mar-23
20039E7D	20036434	07613203014243	18-Mar-20	18-Mar-23
20039E7D	20036458	07613203014243	18-Mar-20	18-Mar-23
20039E7D	20036459	07613203014243	18-Mar-20	18-Mar-23
20039E7D	20036460	07613203014243	18-Mar-20	18-Mar-23
20039E7D	20036461	07613203014243	18-Mar-20	18-Mar-23
20039E7D	20036462	07613203014243	19-Mar-20	19-Mar-23
20039E7D	20036463	07613203014243	19-Mar-20	19-Mar-23
20039E7D	20036464	07613203014243	21-Mar-20	21-Mar-23
20039E7D	20056672	07613203014243	27-May-20	27-May-23



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
20039E7D	20056698	07613203014243	27-May-20	27-May-23
20039E7D	20056699	07613203014243	28-May-20	28-May-23
20039E7D	20056700	07613203014243	28-May-20	28-May-23
20039E7D	20056701	07613203014243	28-May-20	28-May-23
20039E7D	20056702	07613203014243	29-May-20	29-May-23
20039E7D	20066601	07613203014243	18-Jun-20	18-Jun-23
20039E7D	20066602	07613203014243	18-Jun-20	18-Jun-23
20039E7D	20066603	07613203014243	19-Jun-20	19-Jun-23
20039E7D	20066604	07613203014243	19-Jun-20	19-Jun-23
20039E7D	20066605	07613203014243	19-Jun-20	19-Jun-23
20039E7D	20067001	07613203014243	24-Jun-20	24-Jun-23
20039E7D	20067002	07613203014243	24-Jun-20	24-Jun-23
20039E7D	20067003	07613203014243	24-Jun-20	24-Jun-23
20039E7D	20067004	07613203014243	24-Jun-20	24-Jun-23
20039E7D	20067005	07613203014243	25-Jun-20	25-Jun-23
20039E7D	20067006	07613203014243	25-Jun-20	25-Jun-23
20039E7D	20067007	07613203014243	25-Jun-20	25-Jun-23
20039E7D	20067008	07613203014243	25-Jun-20	25-Jun-23
20039E7D	20067218	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067219	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067220	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067221	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067222	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067223	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067224	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067225	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067226	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20067227	07613203014243	30-Jun-20	30-Jun-23
20039E7D	20075225	07613203014243	01-Jul-20	01-Jul-23
20039E7D	20075226	07613203014243	01-Jul-20	01-Jul-23
20039E7D	20075228	07613203014243	01-Jul-20	01-Jul-23
20039E7D	20075229	07613203014243	02-Jul-20	02-Jul-23
20039E7D	20085419	07613203014243	04-Aug-20	04-Aug-23
20039E7D	20085420	07613203014243	04-Aug-20	04-Aug-23
20039E7D	20085510	07613203014243	04-Aug-20	04-Aug-23
20039E7D	20085511	07613203014243	04-Aug-20	04-Aug-23
20039E7D	20085512	07613203014243	04-Aug-20	04-Aug-23
20039E7D	20085768	07613203014243	05-Aug-20	05-Aug-23
20039E7D	20085769	07613203014243	05-Aug-20	05-Aug-23
20039E7D	20085829	07613203014243	05-Aug-20	05-Aug-23
20039E7D	20085830	07613203014243	05-Aug-20	05-Aug-23
20039E7D	20085831	07613203014243	05-Aug-20	05-Aug-23
20039E7D	20085832	07613203014243	05-Aug-20	05-Aug-23
20039E7D	20085833	07613203014243	05-Aug-20	05-Aug-23
20039E7D	20085834	07613203014243	06-Aug-20	06-Aug-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E7D	20085835	07613203014243	06-Aug-20	06-Aug-23
20039E7D	20106574	07613203014243	22-Oct-20	22-Oct-23
20039E7D	20106575	07613203014243	22-Oct-20	22-Oct-23
20039E7D	20106576	07613203014243	22-Oct-20	22-Oct-23
20039E7D	20106577	07613203014243	22-Oct-20	22-Oct-23
20039E7D	20106578	07613203014243	23-Oct-20	23-Oct-23
20039E7D	20106579	07613203014243	23-Oct-20	23-Oct-23
20039E7D	20115320	07613203014243	03-Nov-20	03-Nov-23
20039E7D	20115321	07613203014243	03-Nov-20	03-Nov-23
20039E7D	20115322	07613203014243	03-Nov-20	03-Nov-23
20039E7D	20115323	07613203014243	12-Nov-20	12-Nov-23
20039E7D	20115324	07613203014243	12-Nov-20	12-Nov-23
20039E7D	20115325	07613203014243	18-Nov-20	18-Nov-23
20039E7D	20115326	07613203014243	18-Nov-20	18-Nov-23
20039E7D	20115327	07613203014243	18-Nov-20	18-Nov-23
20039E7D	20115329	07613203014243	12-Nov-20	12-Nov-23
20039E7D	20115330	07613203014243	13-Nov-20	13-Nov-23
20039E7D	20115331	07613203014243	13-Nov-20	13-Nov-23
20039E7D	20115332	07613203014243	13-Nov-20	13-Nov-23
20039E7D	20115333	07613203014243	13-Nov-20	13-Nov-23
20039E7D	20115334	07613203014243	18-Nov-20	18-Nov-23
20039E7D	20115335	07613203014243	12-Nov-20	12-Nov-23
20039E7D	20115336	07613203014243	12-Nov-20	12-Nov-23
20039E7D	20115337	07613203014243	13-Nov-20	13-Nov-23
20039E7D	20115338	07613203014243	14-Nov-20	14-Nov-23
20039E7D	20115339	07613203014243	14-Nov-20	14-Nov-23
20039E7D	20115340	07613203014243	18-Nov-20	18-Nov-23
20039E7D	20115345	07613203014243	12-Nov-20	12-Nov-23
20039E7D	20115346	07613203014243	12-Nov-20	12-Nov-23
20039E7D	20115351	07613203014243	10-Nov-20	10-Nov-23
20039E7D	20116509	07613203014243	18-Nov-20	18-Nov-23
20039E7D	20116510	07613203014243	18-Nov-20	18-Nov-23
20039E7D	20116511	07613203014243	19-Nov-20	19-Nov-23
20039E7D	20125315	07613203014243	02-Dec-20	02-Dec-23
20039E7D	20125316	07613203014243	02-Dec-20	02-Dec-23
20039E7D	20125458	07613203014243	02-Dec-20	02-Dec-23
20039E7D	20125459	07613203014243	02-Dec-20	02-Dec-23
20039E7D	20125460	07613203014243	02-Dec-20	02-Dec-23
20039E7D	20125461	07613203014243	03-Dec-20	03-Dec-23
20039E7D	20125523	07613203014243	03-Dec-20	03-Dec-23
20039E7D	20125524	07613203014243	03-Dec-20	03-Dec-23
20039E7D	20125525	07613203014243	03-Dec-20	03-Dec-23
20039E7D	20125582	07613203014243	03-Dec-20	03-Dec-23
20039E7D	20125583	07613203014243	04-Dec-20	04-Dec-23
20039E7D	21016029	07613203014243	18-Jan-21	18-Jan-24





Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039E7D	21016030	07613203014243	18-Jan-21	18-Jan-24
20039E7D	21016031	07613203014243	18-Jan-21	18-Jan-24
20039E7D	21016032	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016033	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016034	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016035	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016036	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016037	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016038	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016233	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016234	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016235	07613203014243	19-Jan-21	19-Jan-24
20039E7D	21016236	07613203014243	20-Jan-21	20-Jan-24
20039E7D	21016310	07613203014243	20-Jan-21	20-Jan-24
20039E7D	21016311	07613203014243	20-Jan-21	20-Jan-24
20039E7D	21016312	07613203014243	20-Jan-21	20-Jan-24
20039E7D	21016313	07613203014243	20-Jan-21	20-Jan-24
20039E7D	21016314	07613203014243	20-Jan-21	20-Jan-24
20039E7D	21016315	07613203014243	21-Jan-21	21-Jan-24
20039E7D	21016316	07613203014243	21-Jan-21	21-Jan-24
20039E7D	21016317	07613203014243	21-Jan-21	21-Jan-24
20039E7D	21016428	07613203014243	26-Jan-21	26-Jan-24
20039E7D	21016429	07613203014243	26-Jan-21	26-Jan-24
20039E7D	21016430	07613203014243	27-Jan-21	27-Jan-24
20039E7D	21016431	07613203014243	27-Jan-21	27-Jan-24
20039E7D	21016432	07613203014243	27-Jan-21	27-Jan-24
20039E7D	21016598	07613203014243	27-Jan-21	27-Jan-24
20039E7D	21016599	07613203014243	27-Jan-21	27-Jan-24
20039E7D	21016600	07613203014243	27-Jan-21	27-Jan-24
20039E7D	21016601	07613203014243	27-Jan-21	27-Jan-24
20039E7D	21016656	07613203014243	28-Jan-21	28-Jan-24
20039E7D	21016657	07613203014243	28-Jan-21	28-Jan-24
20039E7D	21016658	07613203014243	28-Jan-21	28-Jan-24
20039E7D	21025039	07613203014243	04-Feb-21	04-Feb-24
20039E7D	21025040	07613203014243	05-Feb-21	05-Feb-24
20039E7D	21025041	07613203014243	05-Feb-21	05-Feb-24
20039E7D	21025042	07613203014243	05-Feb-21	05-Feb-24
20039E7D	21025043	07613203014243	05-Feb-21	05-Feb-24
20039E7D	21025512	07613203014243	05-Feb-21	05-Feb-24
20039E7D	21025513	07613203014243	05-Feb-21	05-Feb-24
20039E7D	21025514	07613203014243	05-Feb-21	05-Feb-24
20039E7D	21025515	07613203014243	05-Feb-21	05-Feb-24
20039E7D	21025519	07613203014243	06-Feb-21	06-Feb-24
20039E7D	21025520	07613203014243	06-Feb-21	06-Feb-24
20039E7D	21025521	07613203014243	08-Feb-21	08-Feb-24



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
20039E7D	21025522	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025523	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025524	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025525	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025747	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025748	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025749	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025750	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025751	07613203014243	08-Feb-21	08-Feb-24
20039E7D	21025752	07613203014243	21-Feb-21	21-Feb-24
20039E7D	21025753	07613203014243	21-Feb-21	21-Feb-24
20039E7D	21025754	07613203014243	21-Feb-21	21-Feb-24
20039E7D	21026189	07613203014243	23-Feb-21	23-Feb-24
20039E7D	21026236	07613203014243	23-Feb-21	23-Feb-24
20039E7D	21026239	07613203014243	23-Feb-21	23-Feb-24
20039E7D	21026240	07613203014243	23-Feb-21	23-Feb-24
20039E7D	21026242	07613203014243	24-Feb-21	24-Feb-24
20039E7D	21026243	07613203014243	24-Feb-21	24-Feb-24
20039E7D	21026244	07613203014243	24-Feb-21	24-Feb-24
20039E7D	21035339	07613203014243	03-Mar-21	03-Mar-24
20039E7D	21035351	07613203014243	03-Mar-21	03-Mar-24
20039E7D	21035352	07613203014243	03-Mar-21	03-Mar-24
20039E7D	21035355	07613203014243	03-Mar-21	03-Mar-24
20039E7D	21035356	07613203014243	03-Mar-21	03-Mar-24
20039E7D	21035438	07613203014243	05-Mar-21	05-Mar-24
20039E7D	21035443	07613203014243	05-Mar-21	05-Mar-24
20039E7D	21035445	07613203014243	06-Mar-21	06-Mar-24
20039E7D	21035446	07613203014243	06-Mar-21	06-Mar-24
20039E7D	21035649	07613203014243	06-Mar-21	06-Mar-24
20039E7D	21035650	07613203014243	06-Mar-21	06-Mar-24
20039E7D	21035651	07613203014243	06-Mar-21	06-Mar-24
20039E7D	21035652	07613203014243	08-Mar-21	08-Mar-24
20039E7D	21035653	07613203014243	13-Mar-21	13-Mar-24
20039E7D	21035654	07613203014243	13-Mar-21	13-Mar-24
20039E7D	21035655	07613203014243	13-Mar-21	13-Mar-24
20039E7D	21035656	07613203014243	13-Mar-21	13-Mar-24
20039E7D	21045121	07613203014243	08-Apr-21	08-Apr-24
20039E7D	21045122	07613203014243	08-Apr-21	08-Apr-24
20039E7D	21045123	07613203014243	08-Apr-21	08-Apr-24
20039E7D	21045124	07613203014243	08-Apr-21	08-Apr-24
20039E7D	21045583	07613203014243	09-Apr-21	09-Apr-24
20039E7D	21045584	07613203014243	09-Apr-21	09-Apr-24
20039E7D	21045586	07613203014243	09-Apr-21	09-Apr-24
20039E7D	21045587	07613203014243	09-Apr-21	09-Apr-24
20039R	20026981	10885403275050	27-Feb-20	-



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20039R	20026983	10885403275050	27-Feb-20	-
20039R	20026984	10885403275050	28-Feb-20	-
20039R	20026985	10885403275050	28-Feb-20	-
20039R	20026987	10885403275050	28-Feb-20	-
20039R	20026988	10885403275050	28-Feb-20	-
20039R	20026989	10885403275050	28-Feb-20	-
20039R	20026990	10885403275050	28-Feb-20	-
20039R	20026992	10885403275050	28-Feb-20	-
20049E7D	20026042	07613203014151	14-Feb-20	14-Feb-23
20057E	20096040	10885403234897	30-Sep-20	30-Sep-23
20057E	20106553	10885403234897	23-Oct-20	23-Oct-23
20059E	20015143	07613203012317	02-Jan-20	02-Jan-23
20059E	20015144	07613203012317	02-Jan-20	02-Jan-23
20059E	20015145	07613203012317	02-Jan-20	02-Jan-23
20059E	20015146	07613203012317	03-Jan-20	03-Jan-23
20059E	20015147	07613203012317	03-Jan-20	03-Jan-23
20059E	20015148	07613203012317	03-Jan-20	03-Jan-23
20059E	20015149	07613203012317	03-Jan-20	03-Jan-23
20059E	20015228	07613203012317	04-Jan-20	04-Jan-23
20059E	20015229	07613203012317	04-Jan-20	04-Jan-23
20059E	20015230	07613203012317	05-Jan-20	05-Jan-23
20059E	20015231	07613203012317	08-Jan-20	08-Jan-23
20059E	20026814	07613203012317	26-Feb-20	26-Feb-23
20059E	20026815	07613203012317	26-Feb-20	26-Feb-23
20059E	20026816	07613203012317	27-Feb-20	27-Feb-23
20059E	20035236	07613203012317	03-Mar-20	03-Mar-23
20059E	20035237	07613203012317	03-Mar-20	03-Mar-23
20059E	20035238	07613203012317	03-Mar-20	03-Mar-23
20059E	20035912	07613203012317	11-Mar-20	11-Mar-23
20059E	20035913	07613203012317	13-Mar-20	13-Mar-23
20059E	20035914	07613203012317	12-Mar-20	12-Mar-23
20059E	20035915	07613203012317	13-Mar-20	13-Mar-23
20059E	20035916	07613203012317	14-Mar-20	14-Mar-23
20059E	20035980	07613203012317	15-Mar-20	15-Mar-23
20059E	20035981	07613203012317	14-Mar-20	14-Mar-23
20059E	20035982	07613203012317	14-Mar-20	14-Mar-23
20059E	20046571	07613203012317	25-Apr-20	25-Apr-23
20059E	20046572	07613203012317	25-Apr-20	25-Apr-23
20059E	20046573	07613203012317	27-Apr-20	27-Apr-23
20059E	20046574	07613203012317	27-Apr-20	27-Apr-23
20059E	20046622	07613203012317	27-Apr-20	27-Apr-23
20059E	20057093	07613203012317	31-May-20	31-May-23
20059E	20065337	07613203012317	03-Jun-20	03-Jun-23
20059E	20065338	07613203012317	08-Jun-20	08-Jun-23
20059E	20065339	07613203012317	06-Jun-20	06-Jun-23



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
20059E	20065615	07613203012317	07-Jun-20	07-Jun-23
20059E	20065616	07613203012317	08-Jun-20	08-Jun-23
20059E	20065617	07613203012317	08-Jun-20	08-Jun-23
20059E	20065729	07613203012317	09-Jun-20	09-Jun-23
20059E	20065730	07613203012317	09-Jun-20	09-Jun-23
20059E	20067103	07613203012317	29-Jun-20	29-Jun-23
20059E	20067259	07613203012317	30-Jun-20	30-Jun-23
20059E	20075461	07613203012317	06-Jul-20	06-Jul-23
20059E	20075462	07613203012317	08-Jul-20	08-Jul-23
20059E	20076621	07613203012317	15-Jul-20	15-Jul-23
20059E	20076812	07613203012317	16-Jul-20	16-Jul-23
20059E	20076814	07613203012317	23-Jul-20	23-Jul-23
20059E	20077145	07613203012317	30-Jul-20	30-Jul-23
20059E	20077146	07613203012317	30-Jul-20	30-Jul-23
20059E	20077147	07613203012317	30-Jul-20	30-Jul-23
20059E	20077269	07613203012317	30-Jul-20	30-Jul-23
20059E	20077270	07613203012317	30-Jul-20	30-Jul-23
20059E	20086539	07613203012317	19-Aug-20	19-Aug-23
20059E	20086540	07613203012317	24-Aug-20	24-Aug-23
20059E	20086541	07613203012317	24-Aug-20	24-Aug-23
20059E	20086542	07613203012317	24-Aug-20	24-Aug-23
20059E	20086543	07613203012317	24-Aug-20	24-Aug-23
20059E	20086544	07613203012317	25-Aug-20	25-Aug-23
20059E	20086545	07613203012317	25-Aug-20	25-Aug-23
20059E	20086546	07613203012317	25-Aug-20	25-Aug-23
20059E	20086931	07613203012317	24-Aug-20	24-Aug-23
20059E	20095443	07613203012317	01-Sep-20	01-Sep-23
20059E	20095444	07613203012317	01-Sep-20	01-Sep-23
20059E	20095445	07613203012317	01-Sep-20	01-Sep-23
20059E	20095446	07613203012317	02-Sep-20	02-Sep-23
20059E	20095468	07613203012317	02-Sep-20	02-Sep-23
20059E	20095974	07613203012317	08-Sep-20	08-Sep-23
20059E	20095975	07613203012317	08-Sep-20	08-Sep-23
20059E	20095976	07613203012317	08-Sep-20	08-Sep-23
20059E	20096014	07613203012317	07-Sep-20	07-Sep-23
20059E	20096015	07613203012317	07-Sep-20	07-Sep-23
20059E	20099546	07613203012317	15-Sep-20	15-Sep-23
20059E	20099547	07613203012317	15-Sep-20	15-Sep-23
20059E	20099548	07613203012317	15-Sep-20	15-Sep-23
20059E	20099554	07613203012317	21-Sep-20	21-Sep-23
20059E	20099563	07613203012317	23-Sep-20	23-Sep-23
20059E	20099567	07613203012317	23-Sep-20	23-Sep-23
20059E	20099568	07613203012317	24-Sep-20	24-Sep-23
20059E	20099570	07613203012317	28-Sep-20	28-Sep-23
20059E	20099571	07613203012317	29-Sep-20	29-Sep-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20059E	20109501	07613203012317	08-Oct-20	08-Oct-23
20059E	20109502	07613203012317	08-Oct-20	08-Oct-23
20059E	20109504	07613203012317	08-Oct-20	08-Oct-23
20059E	20109505	07613203012317	12-Oct-20	12-Oct-23
20059E	20109580	07613203012317	23-Oct-20	23-Oct-23
20059E	20109581	07613203012317	26-Oct-20	26-Oct-23
20059E	20109582	07613203012317	26-Oct-20	26-Oct-23
20059E	20109583	07613203012317	27-Oct-20	27-Oct-23
20059E	20109584	07613203012317	27-Oct-20	27-Oct-23
20059E	20109585	07613203012317	28-Oct-20	28-Oct-23
20059E	20109586	07613203012317	28-Oct-20	28-Oct-23
20059E	20109587	07613203012317	29-Oct-20	29-Oct-23
20059E	20109588	07613203012317	29-Oct-20	29-Oct-23
20059E	20109589	07613203012317	29-Oct-20	29-Oct-23
20059E	20109590	07613203012317	30-Oct-20	30-Oct-23
20059E	20119624	07613203012317	18-Nov-20	18-Nov-23
20059E	20119625	07613203012317	19-Nov-20	19-Nov-23
20059E	20119626	07613203012317	19-Nov-20	19-Nov-23
20059E	20119627	07613203012317	20-Nov-20	20-Nov-23
20059E	20119628	07613203012317	23-Nov-20	23-Nov-23
20059E	20119629	07613203012317	23-Nov-20	23-Nov-23
20059E	20119630	07613203012317	24-Nov-20	24-Nov-23
20059E	20119631	07613203012317	24-Nov-20	24-Nov-23
20059E	20119632	07613203012317	25-Nov-20	25-Nov-23
20059E	20119633	07613203012317	18-Nov-20	18-Nov-23
20059E	20119634	07613203012317	18-Nov-20	18-Nov-23
20059E	20119635	07613203012317	18-Nov-20	18-Nov-23
20059E	21019733	07613203012317	05-Jan-21	05-Jan-24
20059E	21019734	07613203012317	05-Jan-21	05-Jan-24
20059E	21019735	07613203012317	05-Jan-21	05-Jan-24
20059E	21019736	07613203012317	05-Jan-21	05-Jan-24
20059E	21029625	07613203012317	11-Feb-21	11-Feb-24
20059E	21029626	07613203012317	12-Feb-21	12-Feb-24
20059E	21029627	07613203012317	15-Feb-21	15-Feb-24
20059E	21029628	07613203012317	16-Feb-21	16-Feb-24
20059E	21029629	07613203012317	16-Feb-21	16-Feb-24
20059E	21029630	07613203012317	16-Feb-21	16-Feb-24
20059E	21029631	07613203012317	17-Feb-21	17-Feb-24
20059E	21029632	07613203012317	18-Feb-21	18-Feb-24
20059E	21029633	07613203012317	18-Feb-21	18-Feb-24
20059E	21029634	07613203012317	18-Feb-21	18-Feb-24
20059E	21029635	07613203012317	19-Feb-21	19-Feb-24
20059E	21029636	07613203012317	19-Feb-21	19-Feb-24
20059E	21029688	07613203012317	12-Feb-21	12-Feb-24
20059E	21029689	07613203012317	15-Feb-21	15-Feb-24



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20059E	21029690	07613203012317	15-Feb-21	15-Feb-24
20059E	21029691	07613203012317	15-Feb-21	15-Feb-24
20059E	21029692	07613203012317	15-Feb-21	15-Feb-24
20059E	21029693	07613203012317	15-Feb-21	15-Feb-24
20059E-0006	20016653	07613203011891	24-Jan-20	24-Jan-23
20059E-0006	20016654	07613203011891	24-Jan-20	24-Jan-23
20059E-0006	20056799	07613203011891	27-May-20	27-May-23
20059E-0006	20086132	07613203011891	10-Aug-20	10-Aug-23
20059E-0006	20106389	07613203011891	14-Oct-20	14-Oct-23
20059E-0006	20106390	07613203011891	14-Oct-20	14-Oct-23
20059E-0006	20106409	07613203011891	15-Oct-20	15-Oct-23
20059E-0006	20116233	07613203011891	12-Nov-20	12-Nov-23
20059E-0006	20116234	07613203011891	13-Nov-20	13-Nov-23
20059E-0006	21036467	07613203011891	24-Mar-21	24-Mar-24
20059E-0006	21036468	07613203011891	24-Mar-21	24-Mar-24
20061E	20096296	07613203011877	22-Sep-20	22-Sep-23
20061E	20106551	07613203011877	20-Oct-20	20-Oct-23
20061E7D	20125787	07613203014106	05-Dec-20	05-Dec-23
20061E7D	21016498	07613203014106	25-Jan-21	25-Jan-24
20062E	20017023	07613203011853	25-Jan-20	25-Jan-23
20062E	20096038	07613203011853	18-Sep-20	18-Sep-23
20062E	20115375	07613203011853	07-Nov-20	07-Nov-23
20062E	20115376	07613203011853	07-Nov-20	07-Nov-23
20062E7D	20116554	07613203014120	18-Nov-20	18-Nov-23
20062E7D	20116555	07613203014120	18-Nov-20	18-Nov-23
20062E7D	20125592	07613203014120	09-Dec-20	09-Dec-23
20062E7D	20125593	07613203014120	10-Dec-20	10-Dec-23
20062E7D	21036533	07613203014120	24-Mar-21	24-Mar-24
20159E	20016205	07613203012270	15-Jan-20	15-Jan-23
20159E	20027096	07613203012270	28-Feb-20	28-Feb-23
20159E	20056457	07613203012270	29-May-20	29-May-23
20159E	20075676	07613203012270	08-Jul-20	08-Jul-23
20159E	20085789	07613203012270	08-Aug-20	08-Aug-23
20159E	20096455	07613203012270	15-Sep-20	15-Sep-23
20159E	20096547	07613203012270	15-Sep-20	15-Sep-23
20159E	20105632	07613203012270	04-Oct-20	04-Oct-23
20159E	20105633	07613203012270	05-Oct-20	05-Oct-23
20159E	20106516	07613203012270	15-Oct-20	15-Oct-23
20159E	20116316	07613203012270	18-Nov-20	18-Nov-23
20159E	20116317	07613203012270	19-Nov-20	19-Nov-23
20159E	21016579	07613203012270	25-Jan-21	25-Jan-24
20159E	21039822	07613203012270	30-Mar-21	30-Mar-24
20159E7D	20105459	07613203014090	02-Oct-20	02-Oct-23
20159E7D	20107206	07613203014090	31-Oct-20	31-Oct-23
20338E	20095143	10885403234927	14-Sep-20	14-Sep-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
20338E	20095144	10885403234927	15-Sep-20	15-Sep-23
21059E	20097132	10885403234972	29-Sep-20	29-Sep-23
22003E-07	20016655	10885403242045	24-Jan-20	24-Jan-21
22003E-07	20016656	10885403242045	25-Jan-20	25-Jan-21
22003E-07	20016657	10885403242045	25-Jan-20	25-Jan-21
22003E-07	20017033	10885403242045	27-Jan-20	27-Jan-21
22003E-07	20017034	10885403242045	27-Jan-20	27-Jan-21
22003E-07	20017035	10885403242045	27-Jan-20	27-Jan-21
22003E-07	20017036	10885403242045	28-Jan-20	28-Jan-21
22003E-07	20017093	10885403242045	28-Jan-20	28-Jan-21
22003E-07	20017094	10885403242045	29-Jan-20	29-Jan-21
22003E-07	20017095	10885403242045	29-Jan-20	29-Jan-21
22003E-07	20027009	10885403242045	27-Feb-20	27-Feb-21
22003E-07	20027010	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20027011	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20027012	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20027013	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20027014	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20027136	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20027137	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20027138	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20027139	10885403242045	28-Feb-20	28-Feb-21
22003E-07	20035235	10885403242045	02-Mar-20	02-Mar-21
22003E-07	20035272	10885403242045	03-Mar-20	03-Mar-21
22003E-07	20035273	10885403242045	03-Mar-20	03-Mar-21
22003E-07	20036333	10885403242045	19-Mar-20	19-Mar-21
22003E-07	20036334	10885403242045	20-Mar-20	20-Mar-21
22003E-07	20036335	10885403242045	20-Mar-20	20-Mar-21
22003E-07	20036374	10885403242045	20-Mar-20	20-Mar-21
22003E-07	20036375	10885403242045	20-Mar-20	20-Mar-21
22003E-07	20036376	10885403242045	20-Mar-20	20-Mar-21
22003E-07	20036638	10885403242045	20-Mar-20	20-Mar-21
22003E-07	20036639	10885403242045	21-Mar-20	21-Mar-21
22003E-07	20036640	10885403242045	21-Mar-20	21-Mar-21
22003E-07	20036641	10885403242045	21-Mar-20	21-Mar-21
22003E-07	20037265	10885403242045	27-Mar-20	27-Mar-21
22003E-07	20037378	10885403242045	31-Mar-20	31-Mar-21
22003E-07	20037379	10885403242045	31-Mar-20	31-Mar-21
22003E-07	20037380	10885403242045	31-Mar-20	31-Mar-21
22003E-07	20069048	10885403242045	24-Jun-20	24-Jun-21
22003E-07	20069049	10885403242045	24-Jun-20	24-Jun-21
22003E-07	20069050	10885403242045	25-Jun-20	25-Jun-21
22003E-07	20079754	10885403242045	22-Jul-20	22-Jul-21
22003E-07	20079755	10885403242045	22-Jul-20	22-Jul-21
22003E-07	20089527	10885403242045	11-Aug-20	11-Aug-21



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
22003E-07	20089538	10885403242045	13-Aug-20	13-Aug-21
22003E-07	20089548	10885403242045	14-Aug-20	14-Aug-21
22003E-07	20089549	10885403242045	14-Aug-20	14-Aug-21
22003E-07	20089553	10885403242045	17-Aug-20	17-Aug-21
22003E-07	20099512	10885403242045	02-Sep-20	02-Sep-21
22003E-07	20109525	10885403242045	16-Oct-20	16-Oct-21
22003E-07	20109526	10885403242045	16-Oct-20	16-Oct-21
22003E-07	20109527	10885403242045	19-Oct-20	19-Oct-21
22003E-07	20109528	10885403242045	19-Oct-20	19-Oct-21
22003E-07	20109557	10885403242045	19-Oct-20	19-Oct-21
22003E-07	20109558	10885403242045	20-Oct-20	20-Oct-21
22003E-07	20109559	10885403242045	20-Oct-20	20-Oct-21
22003E-07	20109560	10885403242045	21-Oct-20	21-Oct-21
22003E-07	20109561	10885403242045	21-Oct-20	21-Oct-21
22003E-07	20109562	10885403242045	21-Oct-20	21-Oct-21
22003E-07	20109563	10885403242045	22-Oct-20	22-Oct-21
22003E-07	20109564	10885403242045	22-Oct-20	22-Oct-21
22003E-07	20109565	10885403242045	22-Oct-20	22-Oct-21
22003E-07	20109579	10885403242045	22-Oct-20	22-Oct-21
22003E-07	20119557	10885403242045	03-Nov-20	03-Nov-21
22003E-07	20119558	10885403242045	03-Nov-20	03-Nov-21
22003E-07	20119559	10885403242045	05-Nov-20	05-Nov-21
22003E-07	20119560	10885403242045	06-Nov-20	06-Nov-21
22003E-07	20119561	10885403242045	11-Nov-20	11-Nov-21
22003E-07	20119562	10885403242045	12-Nov-20	12-Nov-21
22003E-07	20119614	10885403242045	13-Nov-20	13-Nov-21
22003E-07	20119615	10885403242045	13-Nov-20	13-Nov-21
22003E-07	20119616	10885403242045	17-Nov-20	17-Nov-21
22003E-07	20119617	10885403242045	18-Nov-20	18-Nov-21
22003E-07	20119618	10885403242045	18-Nov-20	18-Nov-21
22003E-07	20119619	10885403242045	19-Nov-20	19-Nov-21
22003E-07	20119636	10885403242045	20-Nov-20	20-Nov-21
22003E-07	20119637	10885403242045	23-Nov-20	23-Nov-21
22003E-07	20119638	10885403242045	24-Nov-20	24-Nov-21
22003E-07	20119639	10885403242045	25-Nov-20	25-Nov-21
22003E-07	21039511	10885403242045	03-Jan-21	03-Jan-22
22005E-07T	20036006	10885403250781	21-Mar-20	21-Mar-23
22005E-07T	20065679	10885403250781	30-Jun-20	30-Jun-23
22005E-07T	20115292	10885403250781	11-Nov-20	11-Nov-23
22059E	20076719	10885403235016	23-Jul-20	23-Jul-23
22059E	20076720	10885403235016	23-Jul-20	23-Jul-23
22059E	20085945	10885403235016	10-Aug-20	10-Aug-23
22059E	20097234	10885403235016	30-Sep-20	30-Sep-23
22059E	21016366	10885403235016	19-Jan-21	19-Jan-24
22059E	21016367	10885403235016	19-Jan-21	19-Jan-24





Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2300E	20066502	07613203011679	25-Jun-20	25-Jun-23
2300E	20066503	07613203011679	26-Jun-20	26-Jun-23
2300E	20066504	07613203011679	26-Jun-20	26-Jun-23
2300E	20076251	07613203011679	30-Jul-20	30-Jul-23
2300E	20076352	07613203011679	31-Jul-20	31-Jul-23
2300E	20086483	07613203011679	15-Aug-20	15-Aug-23
2300E	20086484	07613203011679	18-Aug-20	18-Aug-23
2300E	20086485	07613203011679	19-Aug-20	19-Aug-23
2300E	20086486	07613203011679	20-Aug-20	20-Aug-23
2300E	20087129	07613203011679	24-Aug-20	24-Aug-23
2300E	20087130	07613203011679	29-Aug-20	29-Aug-23
2300E	20096280	07613203011679	11-Sep-20	11-Sep-23
2300E	20096281	07613203011679	11-Sep-20	11-Sep-23
2300E	20096282	07613203011679	11-Sep-20	11-Sep-23
2300E	20096283	07613203011679	12-Sep-20	12-Sep-23
2300E	20097163	07613203011679	30-Sep-20	30-Sep-23
2300E	20105074	07613203011679	05-Oct-20	05-Oct-23
2300E	20105075	07613203011679	06-Oct-20	06-Oct-23
2300E-0500	20025450	10885403235023	14-Feb-20	14-Feb-23
2300E-0500	20025451	10885403235023	14-Feb-20	14-Feb-23
2300E-0500	20025452	10885403235023	14-Feb-20	14-Feb-23
2300E-0500	20025531	10885403235023	14-Feb-20	14-Feb-23
2300E-0500	20025532	10885403235023	18-Feb-20	18-Feb-23
2300E-0500	20025533	10885403235023	19-Feb-20	19-Feb-23
2300E-0500	20025534	10885403235023	20-Feb-20	20-Feb-23
2300E-0500	20026527	10885403235023	20-Feb-20	20-Feb-23
2300E-0500	20026529	10885403235023	20-Feb-20	20-Feb-23
2300E-0500	20026532	10885403235023	22-Feb-20	22-Feb-23
2300E-0500	20026533	10885403235023	22-Feb-20	22-Feb-23
2300E-0500	20026534	10885403235023	24-Feb-20	24-Feb-23
2300E-0500	20065580	10885403235023	09-Jun-20	09-Jun-23
2300E-0500	20065581	10885403235023	10-Jun-20	10-Jun-23
2300E-0500	20065582	10885403235023	09-Jun-20	09-Jun-23
2300E-0500	20075239	10885403235023	03-Jul-20	03-Jul-23
2300E-0500	20075240	10885403235023	07-Jul-20	07-Jul-23
2300E-0500	20075241	10885403235023	08-Jul-20	08-Jul-23
2300E-0500	20075242	10885403235023	09-Jul-20	09-Jul-23
2300E-0500	20075447	10885403235023	19-Jul-20	19-Jul-23
2300E-0500	20075448	10885403235023	19-Jul-20	19-Jul-23
2300E-0500	20075449	10885403235023	19-Jul-20	19-Jul-23
2300E-0500	20075450	10885403235023	20-Jul-20	20-Jul-23
2300E-0500	20075451	10885403235023	21-Jul-20	21-Jul-23
2300E-0500	20075452	10885403235023	22-Jul-20	22-Jul-23
2300E-0500	20075453	10885403235023	22-Jul-20	22-Jul-23
2300E-0500	20075618	10885403235023	11-Jul-20	11-Jul-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2300E-0500	20075619	10885403235023	12-Jul-20	12-Jul-23
2300E-0500	20075790	10885403235023	28-Jul-20	28-Jul-23
2300E-0500	20075791	10885403235023	28-Jul-20	28-Jul-23
2300E-0500	20075792	10885403235023	30-Jul-20	30-Jul-23
2300E-0500	20075924	10885403235023	28-Jul-20	28-Jul-23
2300E-0500	20085362	10885403235023	04-Aug-20	04-Aug-23
2300E-0500	20085363	10885403235023	04-Aug-20	04-Aug-23
2300E-0500	20085364	10885403235023	06-Aug-20	06-Aug-23
2300E-0500	20085365	10885403235023	06-Aug-20	06-Aug-23
2300E-0500	20085421	10885403235023	10-Aug-20	10-Aug-23
2300E-0500	20085422	10885403235023	10-Aug-20	10-Aug-23
2300E-0500	20085522	10885403235023	10-Aug-20	10-Aug-23
2300E-0500	20085523	10885403235023	10-Aug-20	10-Aug-23
2300E-0500	20095302	10885403235023	07-Sep-20	07-Sep-23
2300E-0500	20095303	10885403235023	07-Sep-20	07-Sep-23
2300E-0500	20095310	10885403235023	08-Sep-20	08-Sep-23
2300E-0500	20095311	10885403235023	09-Sep-20	09-Sep-23
2300E-0500	20095312	10885403235023	11-Sep-20	11-Sep-23
2300E-0500	20095533	10885403235023	08-Sep-20	08-Sep-23
2300E-0500	20095534	10885403235023	08-Sep-20	08-Sep-23
2300E-0500	20116500	10885403235023	24-Nov-20	24-Nov-23
2300E-0500	20116501	10885403235023	24-Nov-20	24-Nov-23
2300E-0500	21016470	10885403235023	26-Jan-21	26-Jan-24
2300E-0500	21016471	10885403235023	26-Jan-21	26-Jan-24
2300E-0500	21016593	10885403235023	30-Jan-21	30-Jan-24
2309E	20035697	07613203011655	07-Mar-20	07-Mar-23
2309E	20066252	07613203011655	17-Jun-20	17-Jun-23
2309E	20066399	07613203011655	24-Jun-20	24-Jun-23
2309E	20076260	07613203011655	25-Jul-20	25-Jul-23
2309E	20095373	07613203011655	07-Sep-20	07-Sep-23
2309E	20095887	07613203011655	10-Sep-20	10-Sep-23
2309E	20105859	07613203011655	29-Sep-20	29-Sep-23
2309E	20105860	07613203011655	29-Oct-20	29-Oct-23
2309E	20106710	07613203011655	29-Oct-20	29-Oct-23
2309E-0006	20095374	07613203011655	19-Sep-20	19-Sep-23
2309E-0006	20095375	07613203011655	21-Sep-20	21-Sep-23
2309E-0006	20096302	07613203011655	22-Sep-20	22-Sep-23
2309E-0006	20105754	07613203011655	05-Oct-20	05-Oct-23
2309E-0006	20105755	07613203011655	06-Oct-20	06-Oct-23
2309E-0006	20105756	07613203011655	06-Oct-20	06-Oct-23
2309E-0006	20105757	07613203011655	06-Oct-20	06-Oct-23
2309E-0006	20105758	07613203011655	07-Oct-20	07-Oct-23
2309E-0006	20115470	07613203011655	02-Nov-20	02-Nov-23
2440-0600	20015347	07613203011372	23-Jan-20	23-Jan-23
2440-0600	20016918	07613203011372	23-Jan-20	23-Jan-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2440-0600	20016919	07613203011372	23-Jan-20	23-Jan-23
2440-0600	20016960	07613203011372	24-Jan-20	24-Jan-23
2440-0600	20016961	07613203011372	24-Jan-20	24-Jan-23
2441-0006	20077170	07613203020176	23-Jul-20	23-Jul-23
2441-0007	20035102	10885403235047	02-Mar-20	02-Mar-23
2441-0007	20035133	10885403235047	02-Mar-20	02-Mar-23
2441-0007	20035134	10885403235047	03-Mar-20	03-Mar-23
2441-0007	20035281	10885403235047	03-Mar-20	03-Mar-23
2441-0007	20035295	10885403235047	04-Mar-20	04-Mar-23
2441-0007	20035296	10885403235047	04-Mar-20	04-Mar-23
2441-0007	20035368	10885403235047	05-Mar-20	05-Mar-23
2441-0007	20035526	10885403235047	05-Mar-20	05-Mar-23
2441-0007	20035527	10885403235047	05-Mar-20	05-Mar-23
2441-0007	20035603	10885403235047	06-Mar-20	06-Mar-23
2441-0007	20035604	10885403235047	06-Mar-20	06-Mar-23
2441-0007	20056757	10885403235047	22-May-20	22-May-23
2441-0007	20065191	10885403235047	01-Jun-20	01-Jun-23
2441-0007	20065192	10885403235047	02-Jun-20	02-Jun-23
2441-0007	20065281	10885403235047	02-Jun-20	02-Jun-23
2441-0007	20065282	10885403235047	02-Jun-20	02-Jun-23
2441-0007	20065298	10885403235047	02-Jun-20	02-Jun-23
2441-0007	20065299	10885403235047	03-Jun-20	03-Jun-23
2441-0007	20065346	10885403235047	03-Jun-20	03-Jun-23
2441-0007	20065347	10885403235047	03-Jun-20	03-Jun-23
2441-0007	20065362	10885403235047	03-Jun-20	03-Jun-23
2441-0007	20065759	10885403235047	08-Jun-20	08-Jun-23
2441-0007	20065944	10885403235047	09-Jun-20	09-Jun-23
2441-0007	20065945	10885403235047	10-Jun-20	10-Jun-23
2441-0007	20065946	10885403235047	10-Jun-20	10-Jun-23
2441-0007	20075258	10885403235047	02-Jul-20	02-Jul-23
2441-0007	20075510	10885403235047	03-Jul-20	03-Jul-23
2441-0007	20075511	10885403235047	04-Jul-20	04-Jul-23
2441-0007	20075512	10885403235047	04-Jul-20	04-Jul-23
2441-0007	20075577	10885403235047	06-Jul-20	06-Jul-23
2441-0007	20075643	10885403235047	06-Jul-20	06-Jul-23
2441-0007	20075644	10885403235047	06-Jul-20	06-Jul-23
2441-0007	20075660	10885403235047	07-Jul-20	07-Jul-23
2441-0007	20075719	10885403235047	07-Jul-20	07-Jul-23
2441-0007	20075720	10885403235047	07-Jul-20	07-Jul-23
2441-0007	20075721	10885403235047	07-Jul-20	07-Jul-23
2441-0007	20075887	10885403235047	08-Jul-20	08-Jul-23
2441-0007	20075888	10885403235047	08-Jul-20	08-Jul-23
2441-0007	20075889	10885403235047	08-Jul-20	08-Jul-23
2441-0007	20075890	10885403235047	08-Jul-20	08-Jul-23
2441-0007	20076847	10885403235047	17-Jul-20	17-Jul-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2441-0007	20076848	10885403235047	17-Jul-20	17-Jul-23
2441-0007	20076872	10885403235047	18-Jul-20	18-Jul-23
2441-0007	20076873	10885403235047	18-Jul-20	18-Jul-23
2441-0007	20076874	10885403235047	18-Jul-20	18-Jul-23
2441-0007	20076965	10885403235047	20-Jul-20	20-Jul-23
2441-0007	20076966	10885403235047	20-Jul-20	20-Jul-23
2441-0007	20076967	10885403235047	21-Jul-20	21-Jul-23
2441-0007	20076968	10885403235047	21-Jul-20	21-Jul-23
2441-0007	20076969	10885403235047	21-Jul-20	21-Jul-23
2441-0007	20077074	10885403235047	21-Jul-20	21-Jul-23
2441-0007	20077075	10885403235047	22-Jul-20	22-Jul-23
2441-0007	20085132	10885403235047	01-Aug-20	01-Aug-23
2441-0007	20085133	10885403235047	03-Aug-20	03-Aug-23
2441-0007	20085134	10885403235047	03-Aug-20	03-Aug-23
2441-0007	20085135	10885403235047	03-Aug-20	03-Aug-23
2441-0007	20085697	10885403235047	04-Aug-20	04-Aug-23
2441-0007	20085698	10885403235047	04-Aug-20	04-Aug-23
2441-0007	20085699	10885403235047	04-Aug-20	04-Aug-23
2441-0007	20085729	10885403235047	04-Aug-20	04-Aug-23
2441-0007	20085730	10885403235047	05-Aug-20	05-Aug-23
2441-0007	20085731	10885403235047	05-Aug-20	05-Aug-23
2441-0007	20085858	10885403235047	05-Aug-20	05-Aug-23
2441-0007	20085859	10885403235047	05-Aug-20	05-Aug-23
2441-0007	20085875	10885403235047	06-Aug-20	06-Aug-23
2441-0007	20085876	10885403235047	06-Aug-20	06-Aug-23
2441-0007	20086019	10885403235047	06-Aug-20	06-Aug-23
2441-0007	20086020	10885403235047	06-Aug-20	06-Aug-23
2441-0007	20086036	10885403235047	07-Aug-20	07-Aug-23
2441-0007	20086037	10885403235047	07-Aug-20	07-Aug-23
2441-0007	20086213	10885403235047	07-Aug-20	07-Aug-23
2441-0007	20086214	10885403235047	07-Aug-20	07-Aug-23
2441-0007	20086227	10885403235047	08-Aug-20	08-Aug-23
2441-0007	20086228	10885403235047	08-Aug-20	08-Aug-23
2441-0007	20086229	10885403235047	10-Aug-20	10-Aug-23
2441-0007	20086338	10885403235047	10-Aug-20	10-Aug-23
2441-0007	20086339	10885403235047	18-Aug-20	18-Aug-23
2441-0007	20086353	10885403235047	18-Aug-20	18-Aug-23
2441-0007	20086354	10885403235047	18-Aug-20	18-Aug-23
2441-0007	20086491	10885403235047	19-Aug-20	19-Aug-23
2441-0007	20086924	10885403235047	19-Aug-20	19-Aug-23
2441-0007	20086960	10885403235047	19-Aug-20	19-Aug-23
2441-0007	20086965	10885403235047	26-Aug-20	26-Aug-23
2441-0007	20086966	10885403235047	24-Aug-20	24-Aug-23
2441-0007	20087333	10885403235047	28-Aug-20	28-Aug-23
2441-0007	20087334	10885403235047	28-Aug-20	28-Aug-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
2441-0007	20087341	10885403235047	28-Aug-20	28-Aug-23
2441-0007	20087412	10885403235047	29-Aug-20	29-Aug-23
2441-0007	20096174	10885403235047	09-Sep-20	09-Sep-23
2441-0007	20096183	10885403235047	09-Sep-20	09-Sep-23
2441-0007	20096260	10885403235047	10-Sep-20	10-Sep-23
2441-0007	20096366	10885403235047	24-Sep-20	24-Sep-23
2441-0007	20096400	10885403235047	25-Sep-20	25-Sep-23
2441-0007	20096401	10885403235047	25-Sep-20	25-Sep-23
2441-0007	20097100	10885403235047	28-Sep-20	28-Sep-23
2441-0007	20097101	10885403235047	28-Sep-20	28-Sep-23
2441-0007	20097102	10885403235047	29-Sep-20	29-Sep-23
2441-0007	20097157	10885403235047	29-Sep-20	29-Sep-23
2444-0007	20015346	10885403274121	30-Jan-20	30-Jan-23
2447-0007	20056758	10885403235078	25-May-20	25-May-23
2447-0007	20056830	10885403235078	26-May-20	26-May-23
2447-0007	20056850	10885403235078	26-May-20	26-May-23
2447-0007	20056851	10885403235078	27-May-20	27-May-23
2447-0007	20066701	10885403235078	23-Jun-20	23-Jun-23
2447-0007	20066724	10885403235078	23-Jun-20	23-Jun-23
2447-0007	20066775	10885403235078	23-Jun-20	23-Jun-23
2447-0007	20066776	10885403235078	23-Jun-20	23-Jun-23
2447-0007	20066777	10885403235078	24-Jun-20	24-Jun-23
2447-0007	20105671	10885403235078	06-Oct-20	06-Oct-23
2447-0007	20105706	10885403235078	06-Oct-20	06-Oct-23
2447-0007	20125218	10885403235078	08-Dec-20	08-Dec-23
2447-0007	20125219	10885403235078	09-Dec-20	09-Dec-23
42163E	21015553	10885403235313	13-Jan-21	13-Jan-24
42407E-05	20085151	10885403245480	10-Aug-20	10-Aug-23
60103E	20026728	07613203000574	21-Feb-20	21-Feb-23
60103E	20026729	07613203000574	21-Feb-20	21-Feb-23
60103E	20026730	07613203000574	22-Feb-20	22-Feb-23
60103E	20026852	07613203000574	24-Feb-20	24-Feb-23
60103E	20026853	07613203000574	25-Feb-20	25-Feb-23
60103E	20026865	07613203000574	26-Feb-20	26-Feb-23
60103E	20026866	07613203000574	26-Feb-20	26-Feb-23
60103E	20026882	07613203000574	27-Feb-20	27-Feb-23
60103E	20026883	07613203000574	27-Feb-20	27-Feb-23
60103E	20037109	07613203000574	30-Mar-20	30-Mar-23
60103E	20037110	07613203000574	30-Mar-20	30-Mar-23
60103E	20037111	07613203000574	31-Mar-20	31-Mar-21
60103E	20066210	07613203000574	16-Jun-20	16-Jun-23
60103E	20066351	07613203000574	16-Jun-20	16-Jun-23
60103E	20066482	07613203000574	17-Jun-20	17-Jun-23
60103E	20066549	07613203000574	18-Jun-20	18-Jun-23
60103E	20066582	07613203000574	18-Jun-20	18-Jun-23



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
60103E	20076167	07613203000574	09-Jul-20	09-Jul-23
60103E	20076187	07613203000574	09-Jul-20	09-Jul-23
60103E	20076306	07613203000574	10-Jul-20	10-Jul-23
60103E	20076307	07613203000574	10-Jul-20	10-Jul-23
60103E	20076319	07613203000574	11-Jul-20	11-Jul-23
60103E	20076320	07613203000574	13-Jul-20	13-Jul-23
60103E	20076321	07613203000574	13-Jul-20	13-Jul-23
60103E	20076322	07613203000574	14-Jul-20	14-Jul-23
60103E	20076400	07613203000574	14-Jul-20	14-Jul-23
60103E	20076401	07613203000574	15-Jul-20	15-Jul-23
60103E	20076507	07613203000574	15-Jul-20	15-Jul-23
60103E	20076508	07613203000574	16-Jul-20	16-Jul-23
60103E	20105324	07613203000574	07-Oct-20	07-Oct-23
60103E	20105670	07613203000574	07-Oct-20	07-Oct-23
60103E	20116236	07613203000574	13-Nov-20	13-Nov-23
60103E	20116237	07613203000574	13-Nov-20	13-Nov-23
60103E	20116238	07613203000574	14-Nov-20	14-Nov-23
60103E	20116270	07613203000574	26-Nov-20	26-Nov-23
60103E	20116271	07613203000574	27-Nov-20	27-Nov-23
60103E	20116369	07613203000574	27-Nov-20	27-Nov-23
60103E	20117006	07613203000574	28-Nov-20	28-Nov-23
60103E	20117007	07613203000574	30-Nov-20	30-Nov-23
60103E	20117008	07613203000574	30-Nov-20	30-Nov-23
60103E	20125276	07613203000574	01-Dec-20	01-Dec-23
60103E	20125277	07613203000574	01-Dec-20	01-Dec-23
60103E	20125465	07613203000574	02-Dec-20	02-Dec-23
60103E	20125466	07613203000574	02-Dec-20	02-Dec-23
60103E	20125490	07613203000574	03-Dec-20	03-Dec-23
60103E	20125547	07613203000574	04-Dec-20	04-Dec-23
60103E	20125598	07613203000574	04-Dec-20	04-Dec-23
60103E	21015686	07613203000574	18-Jan-21	18-Jan-24
60103E	21035102	07613203000574	02-Mar-21	02-Mar-24
60103E	21035103	07613203000574	03-Mar-21	03-Mar-24
60103E	21035268	07613203000574	03-Mar-21	03-Mar-24
60103E	21035269	07613203000574	04-Mar-21	04-Mar-24
60103E	21035387	07613203000574	05-Mar-21	05-Mar-24
60103E	21035388	07613203000574	05-Mar-21	05-Mar-24
60103E	21035522	07613203000574	06-Mar-21	06-Mar-24
63441E	20076662	07613203010610	30-Jul-20	30-Jul-23
63441E	20096261	07613203010610	10-Sep-20	10-Sep-23
63441E	21046587	07613203010610	29-Apr-21	29-Apr-24
72103E-0006	20025345	07613203010191	06-Feb-20	06-Feb-23
72103E-0006	20025483	07613203010191	27-Feb-20	27-Feb-23
72103E-0006	20037040	07613203010191	25-Mar-20	25-Mar-23
72103E-0006	20037041	07613203010191	26-Mar-20	26-Mar-23



Catalog Number	Lot Number	UDI (GTIN, DI + PI) if applicable	Date of Manufacture	Expiry Date
72103E-0006	20037042	07613203010191	26-Mar-20	26-Mar-23
72103E-0006	20086605	07613203010191	19-Aug-20	19-Aug-23
72103E-0006	20086606	07613203010191	20-Aug-20	20-Aug-23
72103E-0006	20086696	07613203010191	24-Aug-20	24-Aug-23
72103E-0006	20095927	07613203010191	08-Sep-20	08-Sep-23
72103E-0006	20096020	07613203010191	09-Sep-20	09-Sep-23
72103E-0006	20105323	07613203010191	02-Oct-20	02-Oct-23
72103E-0006	20125599	07613203010191	05-Dec-20	05-Dec-23
72103E-0006	20125695	07613203010191	06-Dec-20	06-Dec-23
72103E-0006	20125709	07613203010191	07-Dec-20	07-Dec-23
72103E-0006	21015395	07613203010191	07-Jan-21	07-Jan-24
82104E	20096121	07613203009652	18-Sep-20	18-Sep-23
82104E	20115539	07613203009652	03-Nov-20	03-Nov-23
82113E	21025938	07613203009591	11-Feb-21	11-Feb-24
82113E	21025939	07613203009591	12-Feb-21	12-Feb-24
82113E	21025940	07613203009591	12-Feb-21	12-Feb-24
82113E	21025941	07613203009591	13-Feb-21	13-Feb-24
82113E	21035175	07613203009591	04-Mar-21	04-Mar-24
82113E	21035176	07613203009591	05-Mar-21	05-Mar-24
82113E-0006	20066405	07613203020596	19-Jun-20	19-Jun-23
82113E-0006	20075901	07613203020596	11-Jul-20	11-Jul-23
82113E-0006	20075902	07613203020596	11-Jul-20	11-Jul-23
82113E-0006	20076376	07613203020596	11-Jul-20	11-Jul-23
82113E-0006	20076377	07613203020596	11-Jul-20	11-Jul-23
82113E-0006	20116469	07613203020596	19-Nov-20	19-Nov-23
82113E-0006	21016457	07613203020596	26-Jan-21	26-Jan-24
82113E-0006	21016458	07613203020596	26-Jan-21	26-Jan-24
82113E-0006	21016459	07613203020596	27-Jan-21	27-Jan-24
82113E-0006	21016460	07613203020596	27-Jan-21	27-Jan-24
82113E-0006	21016690	07613203020596	27-Jan-21	27-Jan-24
82113E-0006	21016691	07613203020596	27-Jan-21	27-Jan-24
82113E-0006	21035306	07613203020596	05-Mar-21	05-Mar-24
82113E-0006	21035418	07613203020596	06-Mar-21	06-Mar-24
82113E-0006	21035419	07613203020596	06-Mar-21	06-Mar-24
82113E-0006	21035420	07613203020596	08-Mar-21	08-Mar-24
82113E-0006	21035421	07613203020596	08-Mar-21	08-Mar-24
82113E-0006	21035422	07613203020596	08-Mar-21	08-Mar-24
82113E-0006	21035423	07613203020596	09-Mar-21	09-Mar-24
82113E-0006	21035424	07613203020596	09-Mar-21	09-Mar-24
82114E	20066404	10885403235474	29-Jun-20	29-Jun-23
82114E	21035174	10885403235474	04-Mar-21	04-Mar-24
82115E	20096120	07613203009577	19-Sep-20	19-Sep-23
82115E	20106171	07613203009577	29-Oct-20	29-Oct-23
82115E	20106172	07613203009577	29-Oct-20	29-Oct-23
82115E	20106173	07613203009577	30-Oct-20	30-Oct-23



<b>Catalog Number</b>	<b>Lot Number</b>	<b>UDI (GTIN, DI + PI) if applicable</b>	<b>Date of Manufacture</b>	<b>Expiry Date</b>
82115E	20106174	07613203009577	30-Oct-20	30-Oct-23
82115E	20106175	07613203009577	30-Oct-20	30-Oct-23
C24104E	20017400	10885403235535	31-Jan-20	31-Jan-23
C24104E	20017401	10885403235535	31-Jan-20	31-Jan-23
C24104E	20036962	10885403235535	24-Mar-20	24-Mar-23
C25015E	20085617	10885403235634	07-Aug-20	07-Aug-23
CM42500E-07	20025968	10885403473746	22-Feb-20	22-Feb-21
CM42500E-07	20026013	10885403473746	24-Feb-20	24-Feb-23
CM42500E-07	20035834	10885403473746	10-Mar-20	10-Mar-21
CM42500E-07	20035835	10885403473746	10-Mar-20	10-Mar-21
CS42522E-07	20026132	10885403473753	24-Feb-20	24-Feb-23
CS42522E-07	20026133	10885403473753	25-Feb-20	25-Feb-23
MV0450	20035873	10885403250859	13-Mar-20	13-Mar-23
MV0450	20035874	10885403250859	18-Mar-20	18-Mar-23
MV0450	20045033	10885403250859	02-Apr-20	02-Apr-23
MV0450	20095305	10885403250859	03-Sep-20	03-Sep-23
MV0450	20095306	10885403250859	03-Sep-20	03-Sep-23
MV0450	20095840	10885403250859	04-Sep-20	04-Sep-23
MV0450	20106343	10885403250859	15-Oct-20	15-Oct-23
MV0450	20115582	10885403250859	24-Nov-20	24-Nov-23