

May 31, 2017

To: Surgeons

Subject: **URGENT MEDICAL DEVICE FIELD SAFETY NOTICE (REMOVAL) - LOT SPECIFIC**

Affected Product: Trauma Implants, Reference Attachment 2

Zimmer is initiating a field action for sterile-packaged implants packaged in a single package configuration due to packaging design verification test failures. Specifically, multiple test samples from the configuration failed simulated distribution and shipping testing. The devices impacted are generally the heaviest outlier sizes within the respective product family. Below are photos representative of failures seen during this testing.

Our records indicate that you may have received one or more of the affected products.

It is very likely that any damage to the cavity during surgery. The package insert (instructions) contains a section on sterility. It instructs the user to inspect the cavity is damaged.

Zimmer is initiating a recall for sterile-packaged implants packaged in a single package

the outer tray would be detected prior to surgery with the device or device system contains any damage and not use the device if any seal or




Fig 1. Outer Tray Flange Broken



Fig 2. Inner Tray cracks near the plug

Risks

Immediate health consequences (injuries or illness) that may result from use of or exposure to the product issue.	Most Probable	Highest Severity
	Extension of surgery 10-30 min	Infection
Long range health consequences (injuries or illness) that may result from use of or exposure to the product issue.	Most Probable	Highest Severity
	None	Revision surgery due to Infection

- If compromised packaging is detected during surgery, it may result in a slight delay in surgery to obtain another implant.
- In the unlikely event that compromised packaging is not detected prior to or during surgery by the user, there is a risk of infection, which could result in revision surgery.

Hospital Responsibilities:

1. Review this notification and ensure affected team members are aware of the contents.
2. Complete the Certification of Acknowledgement portion of **Attachment 1**
 - a. Return a digital copy to fielddaction.emea@zimmerbiomet.com within three (3) days.
3. Assist your Zimmer Biomet sales representative quarantine all affected product.
4. Please notify surgeon's that have implanted the affected product within the past 12 months.
5. If after reviewing this notice you have further questions or concerns, please contact your Zimmer Biomet representative. Alternatively, your questions may be sent by email to fielddaction.emea@zimmerbiomet.com.

Surgeon Responsibilities:

1. Review this notification for awareness of the contents.
2. There are no specific patient monitoring instructions related to this recall that are recommended beyond your existing surgical follow up protocol.
3. Complete Attachment 1 – Certificate of Acknowledgement.
 - a. Return a digital copy to fielddaction.emea@zimmerbiomet.com.

Other Information

This field action was reported to the U.S. Food and Drug Administration and will be reported to other Competent Authorities, Notified Bodies, and Regulatory Authorities as required.

Please keep Zimmer Biomet informed of any adverse events associated with this product or any other Zimmer Biomet product by emailing winterthur.per@zimmerbiomet.com.



The undersigned confirms that this notice has been delivered to the appropriate Competent Authorities, as per MEDDEV 2.12-1 revision 8.

We thank you for your cooperation and regret any inconvenience caused by this field action.





ATTACHMENT 1: Certificate of Acknowledgement

ZFA2017-139

By signing below, I acknowledge that the required actions have been taken in accordance with the Recall Notice.

Hospital Facility **Surgeon** (Please select)

Printed Name: _____ Signature: _____

Title: _____ Telephone: () _____ - _____ Date: ____/____/____

Facility Name: _____

Facility Address: _____

City: _____ State: _____ Zip: _____

Country: _____

Note: This form must be returned to Zimmer Biomet before this action can be considered closed for your account. It is important that you complete this form and email a copy to:
fieldaction.emea@zimmerbiomet.com.

Attachment 2

All lot numbers expiring before May 2027 are affected

Item Number	Description
00119309016	VERSA-FXII SC TUBE PLT 90DX16H
00119309516	VERSA-FXII SC TUBE PLT 95DX16H
00119309518	VERSA-FXII SC TUBE PLT 95DX18H
00119309520	VERSA-FXII SC TUBE PLT 95DX20H
00119309522	VERSA-FXII SC TUBE PLT 95DX22H
00119313016	VERSA-FX II STD TUBE 130DX16H
00119313018	VERSA-FX II STD TUBE 130DX18H
00119313020	VERSA-FX II STD TUBE 130DX20H
00119313516	VERSA-FX II STD TUBE 135DX16H
00119313518	VERSA-FX II STD TUBE 135DX18H
00119313520	VERSA-FX II STD TUBE 135DX20H
00119314016	VERSA-FX II STD TUBE 140DX16H
00119314018	VERSA-FX II STD TUBE 140DX18H
00119314020	VERSA-FX II STD TUBE 140DX20H
00119314516	VERSA-FX II STD TUBE 145DX16H
00119314518	VERSA-FX II STD TUBE 145DX18H
00119314520	VERSA-FX II STD TUBE 145DX20H
00119315016	VERSA-FX II STD TUBE 150DX16H
00119315018	VERSA-FX II STD TUBE 150DX18H
00119315020	VERSA-FX II STD TUBE 150DX20H
00224001048	R-F IM NAIL 10MMDX48CM LG
00224001142	R-F IM NAIL 11MMDX42CM
00224001144	R-F IM NAIL 11MMDX44CM
00224001146	R-F IM NAIL 11MMDX46CM LG
00224001148	R-F IM NAIL 11MMDX48CM LG
00224001246	R-F IM NAIL 12MMDX46CM LG
00224001248	R-F IM NAIL 12MMDX48CM LG
00224001340	R-F IM NAIL 13MMDX40CM
00224001342	R-F IM NAIL 13MMDX42CM
00224001344	R-F IM NAIL 13MMDX44CM
00224001346	R-F IM NAIL 13MMDX46CM LG
00224001348	R-F IM NAIL 13MMDX48CM LG
00224001434	R-F IM NAIL 14MMDX34CM
00224001436	R-F IM NAIL 14MMDX36CM
00224001438	R-F IM NAIL 14MMDX38CM
00224001440	R-F IM NAIL 14MMDX40CM

Item Number	Description
00224001442	R-F IM NAIL 14MMDX42CM
00224001444	R-F IM NAIL 14MMDX44CM
00224001446	R-F IM NAIL 14MMDX46CM LG
00224001448	R-F IM NAIL 14MMDX48CM LG
00225230015	FEM IM NAIL 15MMDX30CM
00225230016	FEM IM NAIL 16MMDX30CM
00225232014	FEM IM NAIL 14MMDX32CM
00225232015	FEM IM NAIL 15MMDX32CM
00225232016	FEM IM NAIL 16MMDX32CM
00225234014	FEM IM NAIL 14MMDX34CM
00225234015	FEM IM NAIL 15MMDX34CM
00225234016	FEM IM NAIL 16MMDX34CM
00225236014	FEM IM NAIL 14MMDX36CM
00225236015	FEM IM NAIL 15MMDX36CM
00225238014	FEM IM NAIL 14MMDX38CM
00225238015	FEM IM NAIL 15MMDX38CM
00225240011	FEM IM NAIL 11MMDX40CM
00225240013	FEM IM NAIL 13MMDX40CM
00225240014	FEM IM NAIL 14MMDX40CM
00225240015	FEM IM NAIL 15MMDX40CM
00225242011	FEM IM NAIL 11MMDX42CM
00225242012	FEM IM NAIL 12MMDX42CM
00225242013	FEM IM NAIL 13MMDX42CM
00225242014	FEM IM NAIL 14MMDX42CM
00225242015	FEM IM NAIL 15MMDX42CM
00225244011	FEM IM NAIL 11MMDX44CM
00225244012	FEM IM NAIL 12MMDX44CM
00225244013	FEM IM NAIL 13MMDX44CM
00225244014	FEM IM NAIL 14MMDX44CM
00225246010	FEM IM NAIL 10MMDX46CM
00225246011	FEM IM NAIL 11MMDX46CM
00225246012	FEM IM NAIL 12MMDX46CM
00225246013	FEM IM NAIL 13MMDX46CM
00225246014	FEM IM NAIL 14MMDX46CM
00225248010	FEM IM NAIL 10MMDX48CM
00225248011	FEM IM NAIL 11MMDX48CM

Item Number	Description
00225248012	FEM IM NAIL 12MMDX48CM
00225248013	FEM IM NAIL 13MMDX48CM
00225248014	FEM IM NAIL 14MMDX48CM
00225250013	FEM IM NAIL 13MMDX50CM
00225250014	FEM IM NAIL 14MMDX50CM
00225326015	TIBIAL I/M NAIL 15MMDX26CM
00225328015	TIBIAL I/M NAIL 15MMDX28CM
00225330015	TIBIAL I/M NAIL 15MMDX30CM
00225332014	TIBIAL I/M NAIL 14MMDX32CM
00225332015	TIBIAL I/M NAIL 15MMDX32CM
00225334014	TIBIAL I/M NAIL 14MMDX34CM
00225334015	TIBIAL I/M NAIL 15MMDX34CM
00225336014	TIBIAL I/M NAIL 14MMDX36CM
00225336015	TIBIAL I/M NAIL 15MMDX36CM
00225338014	TIBIAL I/M NAIL 14MMDX38CM
00225338015	TIBIAL I/M NAIL 15MMDX38CM
00225340013	TIBIAL I/M NAIL 13MMDX40CM
00225340014	TIBIAL I/M NAIL 14MMDX40CM
00225340015	TIBIAL I/M NAIL 15MMDX40CM
00225342013	TIBIAL I/M NAIL 13MMDX42CM
00225342014	TIBIAL I/M NAIL 14MMDX42CM
00225342015	TIBIAL I/M NAIL 15MMDX42CM
00225344011	TIBIAL I/M NAIL 11MMDX44CM
00225344013	TIBIAL I/M NAIL 13MMDX44CM
00225344014	TIBIAL I/M NAIL 14MMDX44CM
00225346011	TIBIAL I/M NAIL 11MMDX46CM
00225346012	TIBIAL I/M NAIL 12MMDX46CM
00225346013	TIBIAL I/M NAIL 13MMDX46CM
00225730004	FEM IM NAIL 14MMDX30CM LEFT
00225730014	FEM IM NAIL 14MMDX30CM RIGHT
00225732003	FEM IM NAIL 13MMDX32CM LEFT
00225732004	FEM IM NAIL 14MMDX32CM LEFT
00225732013	FEM IM NAIL 13MMDX32CM RIGHT
00225732014	FEM IM NAIL 14MMDX32CM RIGHT
00225734001	FEM IM NAIL 11MMDX34CM LEFT
00225734003	FEM IM NAIL 13MMDX34CM LEFT
00225734004	FEM IM NAIL 14MMDX34CM LEFT
00225734011	FEM IM NAIL 11MMDX34CM RIGHT
00225734013	FEM IM NAIL 13MMDX34CM RIGHT
00225734014	FEM IM NAIL 14MMDX34CM RIGHT

Item Number	Description
00225736001	FEM IM NAIL 11MMDX36CM LEFT
00225736002	FEM IM NAIL 12MMDX36CM LEFT
00225736003	FEM IM NAIL 13MMDX36CM LEFT
00225736004	FEM IM NAIL 14MMDX36CM LEFT
00225736011	FEM IM NAIL 11MMDX36CM RIGHT
00225736012	FEM IM NAIL 12MMDX36CM RIGHT
00225736013	FEM IM NAIL 13MMDX36CM RIGHT
00225736014	FEM IM NAIL 14MMDX36CM RIGHT
00225738000	FEM IM NAIL 10MMDX38CM LEFT
00225738001	FEM IM NAIL 11MMDX38CM LEFT
00225738002	FEM IM NAIL 12MMDX38CM LEFT
00225738003	FEM IM NAIL 13MMDX38CM LEFT
00225738004	FEM IM NAIL 14MMDX38CM LEFT
00225738010	FEM IM NAIL 10MMDX38CM RIGHT
00225738011	FEM IM NAIL 11MMDX38CM RIGHT
00225738012	FEM IM NAIL 12MMDX38CM RIGHT
00225738013	FEM IM NAIL 13MMDX38CM RIGHT
00225738014	FEM IM NAIL 14MMDX38CM RIGHT
00225740000	FEM IM NAIL 10MMDX40CM LEFT
00225740001	FEM IM NAIL 11MMDX40CM LEFT
00225740002	FEM IM NAIL 12MMDX40CM LEFT
00225740003	FEM IM NAIL 13MMDX40CM LEFT
00225740004	FEM IM NAIL 14MMDX40CM LEFT
00225740010	FEM IM NAIL 10MMDX40CM RIGHT
00225740011	FEM IM NAIL 11MMDX40CM RIGHT
00225740012	FEM IM NAIL 12MMDX40CM RIGHT
00225740013	FEM IM NAIL 13MMDX40CM RIGHT
00225740014	FEM IM NAIL 14MMDX40CM RIGHT
00225742000	FEM IM NAIL 10MMDX42CM LEFT
00225742001	FEM IM NAIL 11MMDX42CM LEFT
00225742002	FEM IM NAIL 12MMDX42CM LEFT
00225742003	FEM IM NAIL 13MMDX42CM LEFT
00225742004	FEM IM NAIL 14MMDX42CM LEFT
00225742010	FEM IM NAIL 10MMDX42CM RIGHT
00225742011	FEM IM NAIL 11MMDX42CM RIGHT
00225742012	FEM IM NAIL 12MMDX42CM RIGHT
00225742013	FEM IM NAIL 13MMDX42CM RIGHT
00225742014	FEM IM NAIL 14MMDX42CM RIGHT
00225744000	FEM IM NAIL 10MMDX44CM LEFT
00225744001	FEM IM NAIL 11MMDX44CM LEFT

Item Number	Description
00225744002	FEM IM NAIL 12MMDX44CM LEFT
00225744003	FEM IM NAIL 13MMDX44CM LEFT
00225744004	FEM IM NAIL 14MMDX44CM LEFT
00225744010	FEM IM NAIL 10MMDX44CM RIGHT
00225744011	FEM IM NAIL 11MMDX44CM RIGHT
00225744012	FEM IM NAIL 12MMDX44CM RIGHT
00225744013	FEM IM NAIL 13MMDX44CM RIGHT
00225744014	FEM IM NAIL 14MMDX44CM RIGHT
00225746000	FEM IM NAIL 10MMDX46CM LEFT
00225746001	FEM IM NAIL 11MMDX46CM LEFT
00225746002	FEM IM NAIL 12MMDX46CM LEFT
00225746003	FEM IM NAIL 13MMDX46CM LEFT
00225746004	FEM IM NAIL 14MMDX46CM LEFT
00225746010	FEM IM NAIL 10MMDX46CM RIGHT
00225746011	FEM IM NAIL 11MMDX46CM RIGHT
00225746012	FEM IM NAIL 12MMDX46CM RIGHT
00225746013	FEM IM NAIL 13MMDX46CM RIGHT
00225746014	FEM IM NAIL 14MMDX46CM RIGHT
00225748000	FEM IM NAIL 10MMDX48CM LEFT
00225748001	FEM IM NAIL 11MMDX48CM LEFT
00225748002	FEM IM NAIL 12MMDX48CM LEFT

Item Number	Description
00225748003	FEM IM NAIL 13MMDX48CM LEFT
00225748010	FEM IM NAIL 10MMDX48CM RIGHT
00225748011	FEM IM NAIL 11MMDX48CM RIGHT
00225748012	FEM IM NAIL 12MMDX48CM RIGHT
00225748013	FEM IM NAIL 13MMDX48CM RIGHT
00225750000	FEM IM NAIL 10MMDX50CM LEFT
00225750001	FEM IM NAIL 11MMDX50CM LEFT
00225750002	FEM IM NAIL 12MMDX50CM LEFT
00225750003	FEM IM NAIL 13MMDX50CM LEFT
00225750010	FEM IM NAIL 10MMDX50CM RIGHT
00225750011	FEM IM NAIL 11MMDX50CM RIGHT
00225750012	FEM IM NAIL 12MMDX50CM RIGHT
00225750013	FEM IM NAIL 13MMDX50CM RIGHT
47249146014	Z NAIL PF 14MM X 46CM UNIV
47249246014	Z NAIL GT 14MM X 46CM R
47249246114	Z NAIL GT 14MM X 46CM L
47249538014	Z NAIL TIBIA 14MM X 38CM UNIV
47249540014	Z NAIL TIBIA 14MM X 40CM UNIV
47249542014	Z NAIL TIBIA 14MM X 42CM UNIV
47249544014	Z NAIL TIBIA 14MM X 44CM UNIV
47494502601	4.5 BROAD SCP PLT 26-H STER