



August 17th, 2018

To: Surgeons/ Hospitals

Subject: **URGENT MEDICAL DEVICE FIELD SAFETY NOTICE – REMOVAL**

Reference: **ZFA2018-00210**

Affected Product: HipLOC™ System Plate Barrel

Attachment 2 – Affected Product batch numbers

Item Number	Item Description
2110-130-003	HIPLOC 130DG PLATE 3 HOLE 38MM BARREL
2110-130-004	HIPLOC 130DG PLATE 4 HOLE 38MM BARREL
2110-130-005	HIPLOC 130DG PLATE 5 HOLE 38MM BARREL
2110-135-002	HIPLOC 135DG PLATE 2 HOLE 38MM BARREL
2110-135-003	HIPLOC 135DG PLATE 3 HOLE 38MM BARREL
2110-135-004	HIPLOC 135DG PLATE 4 HOLE 38MM BARREL
2110-135-005	HIPLOC 135DG PLATE 5 HOLE 38MM BARREL
2110-135-006	HIPLOC 135DG PLATE 6 HOLE 38MM BARREL
2110-135-008	HIPLOC 135DG PLATE 8 HOLE 38MM BARREL
2110-140-002	HIPLOC 140DG PLATE 2 HOLE 38MM BARREL
2110-140-003	HIPLOC 140DG PLATE 3 HOLE 38MM BARREL
2110-140-004	HIPLOC 140DG PLATE 4 HOLE 38MM BARREL
2110-140-005	HIPLOC 140DG PLATE 5 HOLE 38MM BARREL
2110-140-006	HIPLOC 140DG PLATE 6 HOLE 38MM BARREL
2110-150-003	HIPLOC 150DG PLATE 3 HOLE 38MM BARREL

2110-150-004	HIPLOC 150DG PLATE 4 HOLE 38MM BARREL
2110-150-012	HIPLOC 150DG PLATE 12 HOLE 38MM BARREL
2111-130-004	HIPLOC 130DG PLATE 4 HOLE 25MM BARREL
2111-135-003	HIPLOC 135DG PLATE 3 HOLE 25MM BARREL
2111-135-005	HIPLOC 135DG PLATE 5 HOLE 25MM BARREL
2111-135-007	HIPLOC 135DG PLATE 7 HOLE 25MM BARREL
2111-150-002	HIPLOC 150DG PLATE 2 HOLE 25MM BARREL

Biomet UK Ltd is conducting a medical device Field Safety Corrective Action (removal) for certain HipLOC™ barrels following certain complaints regarding an undersized bore on HipLOC plate barrel resulting in the lag screw not fitting through the barrel. The issue is detectable by the user immediately prior its usage when attempting to fit the lag screw into the HipLOC barrel. For the affected item/batch number combinations please see the above in combination with the listed batch numbers in attachment 2.



Picture 1: View of the affected device

Risks		
<i>Describe immediate health consequences (injuries or illness) that may result from use of or exposure to the product issue.</i>	Most Probable	Highest Severity
	<i>Delay in surgery <30 minutes to obtain replacement part.</i>	<i>Abandonment of surgery if no replacement part available.</i>
<i>Describe long range health consequences (injuries or illness) that may result from use of or exposure to the product issue.</i>	Most Probable	Highest Severity
	<i>None</i>	<i>Repeated anesthesia if initial surgery abandoned</i>

Our records indicate that you may have received one or more of the affected products. The affected units were distributed from January 2014.

Hospital Responsibilities:

1. Review this notification and ensure that affected personnel are aware of the contents.
2. If you have affected product at your facility, assist your Zimmer Biomet sales representative and quarantine all affected product. Your Zimmer Biomet sales representative will remove the affected product from your facility.
3. Complete **Attachment 1 – Certificate of Acknowledgement** and send to fielddaction.emea@zimmerbiomet.com. This form must be returned even if you do not have affected products at your facility.
4. Retain a copy of the acknowledgement form with your field action records in the event of a compliance audit of your facility's documentation.
5. If you have further questions or concerns after reviewing this notice, please contact your Zimmer Biomet representative.

Other Information

This medical device Field Safety Notice was reported to all relevant Competent Authorities and the related Notified Body as required under the applicable regulations for Medical Devices as per MEDDEV 2.12-1 in Europe.

Please keep Zimmer Biomet informed of any adverse events associated with this product or any other Zimmer Biomet product by emailing per.uk@zimmerbiomet.com or to your local Zimmer Biomet contact.

Please be aware that the names of user facilities notified are routinely provided to the Competent Authorities for audit purposes.

The undersigned confirms that this notice has been delivered to the appropriate Regulatory Agencies.

We would like to thank you for your co-operation in advance and regret any inconveniences caused by this field action.



ATTACHMENT 1
Certificate of Acknowledgement

IMMEDIATE RESPONSE REQUIRED – TIME SENSITIVE ACTION NEEDED

Affected Product: HipLOC™ Compression Hip Screw System- Plate Barrel

Field Action Reference: ZFA2018-00210

Please return the completed form to your Zimmer Biomet contact person:
fieldaction.emea@zimmerbiomet.com

I received and understood the Field Safety Notice.

Regarding the products:

All inventories for the affected products have been checked and following products are to be returned:

Product Reference	Lot Reference	Number of products returned

OR

The affected products which are unavailable for return have been implanted

By signing below, I acknowledge that the required actions have been taken in accordance with the Field Safety Notice.

Hospital Facility **Surgeon** *(Please check one as applicable)*

Printed Name: _____ **Signature:** _____ **Date:** ____/____/____

Title: _____ **Telephone:** () ____ - _____

Facility Name: _____ **Facility Address:** _____

City: _____ **ZIP:** _____ **Country:** _____

ATTACHMENT 2

Affected Product batch numbers

Material	Material Description	Batch
2110-130-003	STD BRL 3-HOLE PLATE 130 DEG	3311642
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3216454
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3222007
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3222008
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3222009
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3222010
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3222011
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3232282
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3323232
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3330209
2110-130-004	STD BRL 4-HOLE PLATE 130 DEG	3330212
2110-130-005	HIPLOC 130D PL SS 5HL 38BARREL	3231398
2110-130-005	HIPLOC 130D PL SS 5HL 38BARREL	3323238
2110-130-005	HIPLOC 130D PL SS 5HL 38BARREL	3382591
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3222014
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3222015
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3222016
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3222019
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3222020
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3222021
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3271382
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3271383
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3271384
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3271385
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3271386
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3311652
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3311653
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3323239
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680879
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680880
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680881
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680882

2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680883
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680884
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680885
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680886
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3680887
2110-135-002	STD BRL 2-HOLE PLATE 135 DEG	3700268
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3222033
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3222034
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3222035
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3222036
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3231401
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3231402
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271307
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271311
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271317
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271318
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271319
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271320
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271321
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271322
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3271323
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3314404
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3314405
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3314406
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680819
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680820
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680821
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680822
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680823
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680824
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680825
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680826
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680827
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680829
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3680830
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3700272
2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	3700273

2110-135-003	STD BRL 3-HOLE PLATE 135 DEG	R3314406A
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216516
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216517
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216518
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216519
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216520
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216526
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216527
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216528
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216529
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216530
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216531
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216532
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216533
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216534
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216535
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216536
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3216537
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3222064
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3222065
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3222066
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3222067
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270125
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270126
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270127
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270128
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270129
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270130
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270131
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270132
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270133
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270134
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270135
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270136
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270137
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270138
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3270139

2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271353
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271354
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271355
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271357
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271358
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271359
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271360
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271361
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271362
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271364
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3271365
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3307037
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3307038
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3307039
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3307040
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3307047
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3323245
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3323246
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3344888
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355791
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355792
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355794
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355795
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355796
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355797
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355798
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355799
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3355806
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3410850
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3411236
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3411237
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3411238
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3411239
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3411398
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3411399
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3411400
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3411401

2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3412005
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3412006
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3412007
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3412008
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3412009
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3412010
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3474628
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3474629
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3474630
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3474631
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3474632
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3475377
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3475378
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3476133
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3476134
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3476135
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3476136
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3476137
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3482911
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3482912
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3482913
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609491
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609492
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609493
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609494
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609496
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609497
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609498
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609499
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609500
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609501
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609502
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609504
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609508
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3609509
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619241
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619243

2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619245
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619246
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619248
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619249
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619250
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619251
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619252
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619253
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619254
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3619256
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623352
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623353
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623354
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623355
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623356
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623357
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623358
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623359
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623360
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623361
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623362
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623363
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623364
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623365
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3623366
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3649859
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3649860
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700274
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700281
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700283
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700284
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700286
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700287
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700289
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700291
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700293
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700297



2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700299
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700301
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700303
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700307
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700608
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700609
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700610
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700611
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700612
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700613
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3700614
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3732476
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3732477
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3749659
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3749660
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3749661
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3749662
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3749663
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3749664
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3749665
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3823936
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3823937
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3823938
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3823939
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3829535
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3829536
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3843019
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3843020
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3843021
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3843022
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3843023
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3843024
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3847671
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3847672
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3852121
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3871317
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3875764

2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3879198
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3879199
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3879200
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3879201
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3879202
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3889525
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3889526
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3889527
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3889528
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3889529
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3889530
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3892930
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3892931
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3892932
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897093
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897094
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897095
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897096
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897097
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897098
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897099
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897100
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897101
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897102
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897103
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897104
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897105
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897106
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897107
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897108
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897109
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897111
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897112
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897114
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897115
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897116
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3897117

2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3902246
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3902247
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3902248
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3902249
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3902250
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3902251
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3914227
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3914228
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3914229
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918335
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918336
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918337
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918338
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918339
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918340
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918341
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918342
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918343
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918344
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918345
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918346
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918347
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918348
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918349
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918350
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918351
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918352
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918353
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918354
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918355
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918356
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918357
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918358
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918359
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918360
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918361
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918362



2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918363
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918364
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918365
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918366
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918367
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918368
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918369
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918370
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918371
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918372
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	3918373
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6025532
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6062927
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097518
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097519
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097520
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097521
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097522
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097523
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097524
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097525
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097526
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097527
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097528
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097529
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097530
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097531
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097532
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097533
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097534
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097535
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097536
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097537
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097538
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097539
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097540
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097541

2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097542
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097543
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097544
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097545
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097546
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097547
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097548
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097549
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097550
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097551
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097552
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6097553
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6102404
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6102405
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6102406
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6102407
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6102776
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6179081
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6179120
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6182758
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6182759
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6182761
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6187959
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6188952
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6192487
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6192488
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6192489
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6192490
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6192491
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6204649
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6218703
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6224834
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6224835
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6224864
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6224865
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6225930
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6225931

2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6225932
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6229925
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6229926
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6229927
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6229928
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6229929
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6229930
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6239135
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264804
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264805
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264806
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264807
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264808
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264809
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264810
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264811
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264812
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264813
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6264814
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334438
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334440
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334441
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334442
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334443
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334444
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334445
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334446
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6334447
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6345442
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6345443
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6345444
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	6345448
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	R3852119A
2110-135-004	STD BRL 4-HOLE PLATE 135 DEG	R3852120A
2110-135-005	HIPLOC 135DG PLATE SS 5 HOLE 38MM BARREL	3271328
2110-135-005	HIPLOC 135DG PLATE SS 5 HOLE 38MM BARREL	3271329

2110-135-005	HIPLOC 135DG PLATE SS 5 HOLE 38MM BARREL	3271330
2110-135-005	HIPLOC 135DG PLATE SS 5 HOLE 38MM BARREL	3271331
2110-135-005	HIPLOC 135DG PLATE SS 5 HOLE 38MM BARREL	3271332
2110-135-005	HIPLOC 135DG PLATE SS 5 HOLE 38MM BARREL	3271333
2110-135-006	STD BRL 6-HOLE PLATE 135 DEG	3216622
2110-135-006	STD BRL 6-HOLE PLATE 135 DEG	3216623
2110-135-006	STD BRL 6-HOLE PLATE 135 DEG	3271388
2110-135-006	STD BRL 6-HOLE PLATE 135 DEG	3271389
2110-135-008	STD BRL 8-HOLE PLATE 135 DEG	3368204
2110-140-002	STD BRL 2-HOLE PLATE 140 DEG	3135691
2110-140-003	STD BRL 3-HOLE PLATE 140 DEG	3135693
2110-140-003	STD BRL 3-HOLE PLATE 140 DEG	3170092
2110-140-003	STD BRL 3-HOLE PLATE 140 DEG	3466423
2110-140-003	STD BRL 3-HOLE PLATE 140 DEG	3466424
2110-140-003	STD BRL 3-HOLE PLATE 140 DEG	3525134
2110-140-003	STD BRL 3-HOLE PLATE 140 DEG	3535477
2110-140-003	STD BRL 3-HOLE PLATE 140 DEG	3639952
2110-140-003	STD BRL 3-HOLE PLATE 140 DEG	3897110
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3095980
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3135695
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3170094
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3201370
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3296682
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3311665
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3369699
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3466425
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3482917
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3622113
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3664778
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3680647
2110-140-004	STD BRL 4-HOLE PLATE 140 DEG	3916008
2110-140-005	HIPLOC 140D PL SS 5HL 38BARREL	3664779
2110-140-005	HIPLOC 140D PL SS 5HL 38BARREL	3741838
2110-140-005	HIPLOC 140D PL SS 5HL 38BARREL	3916009



2110-140-006	STD BRL 6-HOLE PLATE 140 DEG	3639953
2110-150-003	HIPLOC 150D PL SS 3HL 38BARREL	3201372
2110-150-003	HIPLOC 150D PL SS 3HL 38BARREL	3222131
2110-150-004	HIPLOC 150D PL SS 4HL 38BARREL	3710377
2110-150-012	HIPLC 150D PL SS 12HL 38BARREL	2993408
2111-130-004	SHT BRL 4-HOLE PLATE 130 DEG	3847673
2111-130-004	SHT BRL 4-HOLE PLATE 130 DEG	3976183
2111-130-004	SHT BRL 4-HOLE PLATE 130 DEG	6192486
2111-135-003	SHT BRL 3-HOLE PLATE 135 DEG	3113626
2111-135-003	SHT BRL 3-HOLE PLATE 135 DEG	3355815
2111-135-003	SHT BRL 3-HOLE PLATE 135 DEG	3502221
2111-135-003	SHT BRL 3-HOLE PLATE 135 DEG	R3113626A
2111-135-005	HIPLOC 135D PL SS 5HL 25BARREL	3152144
2111-135-005	HIPLOC 135D PL SS 5HL 25BARREL	3222132
2111-135-005	HIPLOC 135D PL SS 5HL 25BARREL	R3152144A
2111-135-005	HIPLOC 135D PL SS 5HL 25BARREL	R3222132A
2111-135-007	HIPLOC 135D PL SS 7HL 25BARREL	3502233