

# **Recall Action Notification**

Philips HeartStart FRx and HS1 OnSite Automated External Defibrillators (AEDs)

### Important information on the System for Australian Recall Actions

The TGA publishes information about therapeutic goods supplied in the Australian market that have been subject to a recall action in a publicly searchable database.

Recall action means action taken by the responsible entity (being the person who is responsible for taking the recall action) to resolve a problem with therapeutic goods supplied in the Australian market that have, or may potentially have, deficiencies relating to safety, quality, efficacy (performance) or presentation.

- Recall actions include: the permanent removal of therapeutic goods from supply in the market, the
  taking of corrective action in relation to therapeutic goods (such as repair, modification, adjustment or
  relabelling) and, in the case of medical devices that have been implanted into patients, the issuing of a
  hazard alert containing information for health practitioners on how to manage patients.
- More information about Australian recall actions is available at <a href="http://tga.gov.au/safety/recalls-about.htm">http://tga.gov.au/safety/recalls-about.htm</a>
- If you are taking a medicine, using a medical device or have had a medical device implanted into you, that is the subject to a recall action, and you have any concerns you should seek advice from a health professional. <a href="http://www.healthdirect.org.au/">http://www.healthdirect.org.au/</a>

#### About the release of this information

While reasonable care is taken to ensure that the information is an accurate record of recall actions that responsible entities have reported to the TGA or of which the TGA has become aware, the TGA does not guarantee or warrant the accuracy, reliability, completeness or currency of the information or its usefulness in achieving any purpose.

To the fullest extent permitted by law, including but not limited to section 61A of the Therapeutic Goods Act 1989, the TGA will not be liable for any loss, damage, cost or expense incurred in or arising by reason of any person relying on this information.

The information contained in the SARA database is released under s 61(5C) of the Therapeutic Goods Act 1989. Copyright restrictions apply to the System of Australian Recall actions (SARA) <a href="http://tga.gov.au/about/website-copyright.htm">http://tga.gov.au/about/website-copyright.htm</a>.

## Recall detail

Type of Product <sup>i</sup>	Medical Device
TGA Recall Referenceii	RC-2019-RN-01388-1
Product Name/Description <sup>iii</sup>	Philips HeartStart FRx and HS1 OnSite Automated External Defibrillators (AEDs)
	FRx Model: 861304
	Serial Number: B17G-05063
	HS1 Model: M5066A Serial Numbers: A15F-00203, A15F-00836, A16A-05193, A18B-01898, A18B-04580
	ARTG 92346 (Philips Electronics Australia Ltd Defibrillator - Non-rechargeable semi-automatic external defibrillator)
Recall Action Leveliv	Hospital
Recall Action Classification <sup>v</sup>	Class I
Recall Action Commencement Datevi	20/09/2019
Responsible Entity <sup>vii</sup>	Philips Electronics Australia Ltd
Reason / Issue <sup>viii</sup>	During a quality review of the supplier's documentation, Philips determined that the supplier had not followed Philips-specified procedures when repairing Printed Circuit Assembly (PCAs) that initially failed required testing. After being repaired, these PCAs passed all required testing at the supplier's facility prior to being shipped to Philips and assembled into AEDs. However, Philips quality standards require that the PCAs both be repaired according to specified procedures and pass testing.
	A higher failure rate on the impacted PCAs may be possible, but there is no expectation that these will lead to an increased rate of harm due to the design mitigations to detect component failures.
Recall Actionix	Recall
Recall Action Instructions <sup>x</sup>	Philips are advising they will replace affected AEDs at no charge. A Philips representative will be in contact with customers to arrange the replacement of affected devices.
	Customers are requested to keep all AED accessories and transfer them to the replacement AED (pads, battery, etc).
	While awaiting for a replacement AED, customers may continue to use affected devices. In the event the device fails routine self tests (daily, weekly and monthly), built in alarms will be activated indicating the device is not ready for use.
Contact Information <sup>xi</sup>	1800 251 400 - Philips Customer Care Centre

#### **Footnotes**

- <sup>i</sup> Type of Product: Medicine, Medical Device, or Biological
- ii TGA Recall Reference: Unique number given by the TGA
- iii Product Name/Description: Brand name (including active ingredient for medicines) and may include generic reference for the kind of medical devices. Includes all necessary information such as affected: catalogue / model and / or batch / serial numbers.
- iv Recall Action Level: The level to which the recall action is to be undertaken. This is based on the significance of the risk and the channels through which the goods have been distributed. The recall action levels are / Wholesale / Hospital / Retail / Consumer.
  - Wholesale includes wholesalers and state purchasing authorities.
  - Hospital includes nursing homes and institutions, hospital pharmacists, ambulance services, blood and tissue banks and laboratories as well as wholesale as appropriate.
  - Retail includes retail pharmacists, medical, dental and other health care professionals as well as wholesale and hospital as appropriate.
  - Consumer includes patients and consumers, as well as wholesale, hospital and retail levels as appropriate.
- <sup>v</sup> Recall Action Classification: Recall actions of therapeutic goods are classified based on the potential risk the deficiency poses to patients / consumers. They are classified as Class I, Class II or Class III.
  - Class I recall action occurs when the product deficiency is potentially life-threatening or could cause a serious
    risk to health.
  - Class II recall action occurs when the product deficiency could cause illness, injury or result in mistreatment, but are not Class I.
  - Class III recall action occurs when the product deficiency may not pose a significant hazard to health, but action may be initiated for other reasons eg. quality related issues.
- vi Recall Action Commencement Date: The date the recall strategy and communication was agreed by the TGA.
- vii Responsible Entity: Sponsor / Supplier / Importer responsible for the recall actions.
- viii Reason / Issue: Reason for the recall action.
- ix Recall Action: Recall action is an action taken to resolve a problem with a therapeutic good already supplied in the market for which there are issues or deficiencies in relation to safety, quality, efficacy (performance) or presentation. There are three distinct recall actions recall, recall for product correction and hazard alert.
  - Recall The permanent removal of an affected therapeutic good from supply or use in the market.
  - Recall for product correction Repair, modification, adjustment or re-labelling of a therapeutic good. The corrective action may take place at the user's premises or any other agreed location.
  - Hazard alert Information issued to healthcare professionals about issues or deficiencies relating to an implanted medical device or biological product and advice about the ongoing management of patients.
- <sup>x</sup> Recall Action Instructions: What the customer should do.
- xi Contact Information: Who the customer should contact for additional information and clarification regarding the recall action.