

Cobalt Chromium Coronary Stent System

Yukon[®] CC

Ordering information*

Small vessel design (SV)

product matrix

Balloon Ø [mm]	Stent length [mm] & Article number							
	8	12	16	18	21	24	28	32
Ø 2,00	T-CMG2008	T-CMG2012	T-CMG2016	T-CMG2018	T-CMG2021	T-CMG2024	T-CMG2028	T-CMG2032
Ø 2,25	T-CMG22508	T-CMG22512	T-CMG22516	T-CMG22518	T-CMG22521	T-CMG22524	T-CMG22528	T-CMG22532
Ø 2,50	T-CMG2508	T-CMG2512	T-CMG2516	T-CMG2518	T-CMG2521	T-CMG2524	T-CMG2528	T-CMG2532

Medium vessel design (MV)

Balloon Ø [mm]	Stent length [mm] & Article number								
	8	12	16	18	21	24	28	32	40
Ø 2,75	T-CMG27508	T-CMG27512	T-CMG27516	T-CMG27518	T-CMG27521	T-CMG27524	T-CMG27528	T-CMG27532	T-CMG27540
Ø 3,00	T-CMG3008	T-CMG3012	T-CMG3016	T-CMG3018	T-CMG3021	T-CMG3024	T-CMG3028	T-CMG3032	T-CMG3040
Ø 3,50	T-CMG3508	T-CMG3512	T-CMG3516	T-CMG3518	T-CMG3521	T-CMG3524	T-CMG3528	T-CMG3532	T-CMG3540
Ø 4,00	T-CMG4008	T-CMG4012	T-CMG4016	T-CMG4018	T-CMG4021	T-CMG4024	T-CMG4028	T-CMG4032	T-CMG4040

* Please contact our customer service for available sizes

Compliance chart

Balloon Ø [mm]	Inflation pressure [bar or 10 ⁻⁵ Pascal]														
	6	7	8	9	10	NP					RBP				
Ø 2,00	1,83	1,87	1,90	1,93	1,96	2,00	2,03	2,06	2,10	2,13	2,16	2,20	2,23	2,26	2,29
Ø 2,25	2,08	2,11	2,14	2,18	2,21	2,25	2,28	2,31	2,35	2,38	2,42	2,45	2,48	2,52	2,55
Ø 2,50	2,33	2,36	2,40	2,43	2,47	2,50	2,53	2,57	2,60	2,64	2,67	2,70	2,74	2,77	2,81
Ø 2,75	2,58	2,61	2,65	2,68	2,71	2,75	2,78	2,81	2,85	2,88	2,91	2,94	2,98	3,01	3,04
Ø 3,00	2,81	2,85	2,89	2,92	2,96	3,00	3,04	3,07	3,11	3,15	3,18	3,22	3,26	3,29	3,33
Ø 3,50	3,29	3,34	3,38	3,42	3,46	3,50	3,55	3,59	3,63	3,67	3,71	3,76	3,80	3,84	3,88
Ø 4,00	3,75	3,80	3,85	3,90	3,95	4,00	4,06	4,11	4,16	4,21	4,26	4,31	4,36	4,41	4,46

CE 0124

© 2010 Translumina GmbH

All Rights Reserved

® YUKON, PEARL Surface and translumina are registered trademarks of translumina GmbH

Please refer to the Instruction for Use supplied with these devices for indications, contraindications, adverse effects, suggested procedures, warnings and precautions.

Translumina GmbH | Neue Rottenburger Strasse 50
72379 Hechingen | Germany | info@translumina.de |
t + 49 7471 98 94 - 0 | f + 49 7471 98 94 - 360



www.translumina.de

translumina[®]