

Attachment to Cyberonics Declaration of Conformity Cyberonics Product List

The following list identifies the Cyberonics part numbers for CE Mark products that Cyberonics intends to distribute in conformity with the Active Implantable Medical Device Directive (90/385/EEC).

These products are further identified with either a lot number or serial number that can be traced back to the applicable top-level assembly number. In all cases, the final sales assembly document number for each product is listed. These assembly documents contain multiple tabulations that indicate single language or multi-lingual product assemblies that are then indicated on a shipment matrix, which is used at the time of order fulfillment. Please note that the package labels with an assembly level printed on them may indicate a sub-assembly level of the final sales level for CE Mark product.

| Product/Model | Part Number (Ref) |
|---|-------------------------|
| <i>Sterile Product</i> | |
| Model 102 | 10-0006-4200 |
| Model 102R | 10-0007-5800 |
| Model 103 | 10-0008-6503 |
| Model 104 | 10-0008-6604 |
| Model 105 | 10-0008-9205 |
| Model 106 | 10-0009-1300 |
| Model 302 (2mm) | 10-0007-0704 |
| Model 302 (3mm) | 10-0007-0706 |
| Model 303 (2mm) | 10-0006-7104 |
| Model 303 (3mm) | 10-0006-7106 |
| Model 304 (2mm) | 10-0007-9504 |
| Model 304 (3mm) | 10-0007-9506 |
| Model 402 | 10-0006-7800 |
| Model 502 | 10-0007-0800 |
| <i>Non-Sterile Product</i> | |
| Model 201 | 10-0008-0600 |
| Model 220 | 10-0008-3200 |
| Model 250, version 8.1 | 10-0009-3300 (Envelope) |
| | 10-0009-3400 (HH X50) |
| | 10-0009-3500 (HH X5) |
| Model 250 Version 11.0 Tablet Computer | 10-0009-7900 |

← Accepted