

## Attachment to Cyberonics Declaration of Conformity Cyberonics Product List

The following list identifies the Cyberonics part numbers for CE Mark products that Cyberonics intends to distribute in conformity with the Active Implantable Medical Device Directive (90/385/EEC).

These products are further identified with either a lot number or serial number that can be traced back to the applicable top-level assembly number. In all cases, the final sales assembly document number for each product is listed. These assembly documents contain multiple tabulations that indicate single language or multi-lingual product assemblies that are then indicated on a shipment matrix, which is used at the time of order fulfillment. Please note that the package labels with an assembly level printed on them may indicate a sub-assembly level of the final sales level for CE Mark product.

Product/Model	Part Number (Ref)
<i>Sterile Product</i>	
Model 102	10-0006-4200
Model 102R	10-0007-5800
Model 103	10-0008-6503
Model 104	10-0008-6604
Model 105	10-0008-9205
Model 106	10-0009-1300
Model 302 (2mm)	10-0007-0704
Model 302 (3mm)	10-0007-0706
Model 303 (2mm)	10-0006-7104
Model 303 (3mm)	10-0006-7106
Model 304 (2mm)	10-0007-9504
Model 304 (3mm)	10-0007-9506
Model 402	10-0006-7800
Model 502	10-0007-0800
<i>Non-Sterile Product</i>	
Model 201	10-0008-0600
Model 220	10-0008-3200
Model 250, version 8.1	10-0009-3300 (Envelope)
	10-0009-3400 (HH X50)
	10-0009-3500 (HH X5)
Model 250 Version 11.0 Tablet Computer	10-0009-7900

← Accepted